AirLife™ Disposable Manometer

Instructions for Use

REF	Catalogue number	R_X ONLY	U.S. Federal law restricts this device to sales by or on the order of a physician
LOT	Lot Number		
QTY	Quantity		Manufacturer
	Not made with natural rubber latex		Single Use
	Consult Instructions for use		Date of Manufacture
	Caution		

EN English

Intended Use: The disposable pressure manometer is used to provide visual indication of patient's airway pressure during ventilation. It may be attached to the pressure manometer port of the resuscitation bag or hyperinflation bag.

Caution: If manometer is not used, the pressure port must be capped.

Caution: Check function and accuracy of the manometer prior to use, including following repeated use of on the same patient. A minimal amount of leakage of airway gas is normal. The effects of this leakage must be evaluated for each patient.

Caution: This device must be used only by persons who have received adequate training in the use of airway pressure monitoring and ventilation devices.

Caution: During use, the patient's condition must be monitored. If for any reason the performance of this device is erratic, it must immediately be replaced or removed and the pressure port be capped.

Instructions for Use

1. Locate the pressure sensing port.
2. Occlude the patient port of ventilation device, pressurize the manometer to 60 cmH₂O (maximum) and release. Note that the indicator returns smoothly to zero. If indicator does not return to zero, replace manometer.
3. Connect resuscitator to artificial airway and ventilate patient.
4. The "O" ring on the body of the manometer may be set as a point of reference to facilitate the desired pressure reading.
5. When repeated use is required on the same patient, check the function and accuracy of the manometer prior to each use.
6. Discard after use.
7. When not in use, the pressure sensing port must be capped.

Manometer accuracy

 AirLife
2710 Northridge Dr. NW, Suite A
Grand Rapids, MI 49544 USA
www.myAirLife.com

(+/-) 1 cmH₂O from 0-10 cmH₂O (hPa)
(+/-) 2 cmH₂O from 10-40 cmH₂O (hPa)
(+/-) 3 cmH₂O from 40 cmH₂O (hPa)

This product is not designed or validated to be reused in different patients. Reuse may cause a risk of cross-contamination, affect the measurement accuracy, system performance, or cause a malfunction as a result of the product being physically damaged due to cleaning, disinfection, re-sterilization, or reuse.

The reuse of the original device on the same patient is acceptable. The user must ensure that the product is not damaged or contaminated between usages on the same patient.

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