AirLife™ Ventilator Circuits

Instruction For Use

**REF** Catalogue number

**Rx ONLY** U.S. Federal law restricts this device to sales by or on the order of a physician

**LOT** Lot Number

**Manufacturer**

**QTY** Quantity

**Date of manufacture**

**LATEX** Not made with natural rubber latex

**Single Use**

**Consult Instructions for Use**

**Caution**

**PHT** This device contains Phthalates

**DEHP**

EN English

Intended Use: A ventilator circuit is used as a conduit for gases between a ventilator and a patient during ventilation of the patient.

Test circuit prior to use:

- Prior to any patient use, ensure that all connections (and, if applicable, other accessories, such as water traps) are tight and have not loosened during transit. All circuits should be pressure tested for obstructions, occlusions or leaks in accordance with the ventilator manufacturer's specifications.
- If any port openings on the patient wye are not utilized, ensure that the attached cover/port caps are secured to prevent leaks.
- Ventilator monitoring and warning systems should be used as recommended by the ventilator manufacturer.
- Condensation accumulation within the circuit should be continually monitored during use and evacuated routinely to reduce the hazard of accidental aspiration by the patient.

Directions for use:

If applicable, attach ball hanger to support arm of ventilator.

- Inspiratory Supply Line:** Attach the connector of the patient inspiratory limb to the outlet of the humidifier or ventilator. Connect the fresh gas supply line to the inlet port of the humidifier. Connect the opposite end of the fresh gas supply line to the ventilator.

**AU REP**
AirLife Australia Holdings PTY LTD
PO Box 97
North Ryde BC, NSW, 1670
Australia

**AirLife**
2710 Northridge Dr. NW, Suite A
Grand Rapids, MI 49544 USA
www.myAirLife.com
Made in Mexico

AirLife®

- Expiratory Supply Line:** Connect exhalation valve drive line to the exhalation valve stem on the valve. Attach the other end to the "expiratory valve" fitting on the ventilator. Make sure that connections are secure. If the circuit does not have the exhalation valve drive line, then connect the 22mm I.D. expiratory supply line to the expiratory port on the ventilator. Neonatal circuits use 15mm O.D. fittings. *If the proximal airway sensing line is utilized, securely attach one end of the sensing line to the barb on the right arm of the "Y". Securely attach the opposite end of the proximal airway sensing line to the appropriate outlet. Plug the barb if the proximal airway line is not utilized.
- Temperature Probe Adapter "T":** Proximal airway temperatures can be monitored by placing temperature probe in the port of the cuff on the left or inspiratory arm of the "Y". Plug the port on the cuff if temperature monitoring is not performed. If the humidifier utilizes a dual temperature probe set, replace the non ported 30°vent connector with the ported 30°vent connector provided.
- Water Trap:** To remove water accumulation, twist jar clockwise and remove. Water trap will seal with no loss of ventilation to patient. Empty water from jar and reassemble jar to water trap top by twisting in a counterclockwise direction. Ensure that all water traps are tight prior to patient use.

WARNING:

- This circuit may be delivered with supply lines disconnected and "Y" ports unplugged. The connections must be completed prior to patient use. Improper assembly may result in ventilator malfunction or inadequate ventilation.
- Ventilator/machine monitoring and warning systems should be used as recommended by the ventilator manufacturer. Before attaching this circuit to the patient, ensure that all appropriate alarm systems are operational and turned on. Before instituting mechanical ventilation, ensure that the circuit and ventilator have been tested as recommended by the hospital and ventilator manufacturer's procedures for ventilator safety. Patient, ventilator and circuit must be monitored on a regular basis as per established hospital standards of care.

Dispose of all materials in accordance with local, state, and federal regulations. Decontaminate and dispose of all potentially biohazardous material.

This single-use product is not designed or validated to be reused. Reuse may cause a risk of cross-contamination, affect the measurement accuracy, system performance, or cause a malfunction as a result of the product being physically damaged due to cleaning, disinfection, re-sterilization, or reuse. A discontinued use of this product on the same patient, with periods of use separated by a period of non-use, is acceptable. The user must ensure that the product is not damaged or contaminated between usages on the same patient.

Risks and Precautionary Measures related to Phthalates: This instruction pertains to the phthalate symbol marked on the device or its packaging. If this device is used for the treatment of children, or treatment of pregnant or nursing women; please note that the following types of procedures may increase the risk of exposure to phthalates: Exchange transfusion in neonates, total parenteral nutrition in neonates, multiple procedures in sick neonates, haemodialysis in peripubertal males, male foetus and male infant of pregnant women, and lactating women; and massive blood infusion into trauma patients. Although these procedures have the potential for increased risk of exposure, conclusive evidence of human health risks has not been established. As a precautionary measure, to reduce the potential for unnecessary exposures to phthalates, the product must be used in accordance with the instructions for use, and practitioners should refrain from using this product beyond the period of time the product is medically necessary or needed.

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Becton Dickinson İth. İhr. Ltd.