AirLife™ Meconium Suction Device (Kurtis MSD) Instructions for Use

REF	Catalogue number	R_X ONLY	U.S. Federal law restricts this device to sales by or on the order of a physician
LOT	Lot Number		Manufacturer
QTY	Quantity		Single Use
	Not made with natural rubber latex	STERILE R	Sterilized using irradiation
	Consult Instructions for Use	PHT DEHP	This device contains Phthalates
	Do not use if package is damaged		Inner Diameter
OD X Length	Outer Diameter X Length		Caution
	Exp. Date.		

EN

Indications For Use

The MSD (Meconium Suction Device) is indicated for the suctioning of meconium secretions from the oropharynx and trachea of newborn infants. It can be easily adapted to allow endotracheal positive pressure ventilation during resuscitation.

Instructions

- Remove MSD (Meconium Suction Device) from package and bend the tracheal tube portion to desired shape.
- Check to ensure that suction adapter is securely attached to the tracheal tube connector, and that suction adapter is in the desired position.
- Connect suction adapter (A) to suction tubing (not included); turn on suction source and set regulator for not more than 100mm Hg vacuum with the suction tubing occluded.
- Perform laryngoscopy according to standard technique.
- If the glottis is obstructed by secretions, suction the oropharynx with the MSD by occluding the venting port with a finger. **Do not attempt to intubate the patient with the venting port (B) occluded.**
- When the glottis is visualized, insert the tracheal tube portion of the MSD into the trachea with the venting port (B) open.
- After the tracheal tube portion is inserted into the trachea, occlude the venting port and withdraw completely from the trachea. The venting port (B) should remain occluded throughout the withdrawal of the tracheal tube portion. The laryngoscope should remain in place, with the glottis visualized, until the decision is made that no further intubations are indicated and that positive pressure ventilation is not indicated.
- If repeated intubations are indicated, re-insert the

AirLife Australia Holdings Pty Ltd
PO Box 97
North Ryde BC, NSW, 1670
Australia

AirLife
2710 Northridge Dr. NW, Suite A
Grand Rapids, MI 49544 USA
www.myAirLife.com

- tracheal tube portion into the trachea with the venting port (B) open. Withdraw as in step 7. The intubation and suction procedure should be repeated as many times as indicated using clinical judgment.
- If it becomes necessary to flush meconium from the MSD, insert tip into sterile water and occlude venting port (B).
 - If positive pressure ventilation is indicated, with the laryngoscope still in place, insert the tracheal tube portion of the MSD into the trachea with the venting port (B) open, as in step 6. If meconium particulate is evident on the tracheal tube portion from prior suctioning, a new tracheal tube must be used for positive pressure ventilation.
 - Remove the suction adapter from the tracheal tube connector, and attach a resuscitation bag to the connector (C).
 - Occluding the venting port (B) with finger or with attached plug, provide positive pressure ventilation as indicated.

Cautions

- To prevent trauma to the patient's airway, do not occlude the venting port during the introduction of the endotracheal tube (ETT) portion of the MSD into the trachea.
- To prevent trauma to the tracheal mucosa, suction should not exceed 100mm Hg with the suction tubing occluded.
- Do not attempt to provide positive pressure ventilation using the MSD if particulate meconium is evident on the sides of the tube from prior suctioning.
- Contents sterile in unopened, undamaged package.


Risks and Precautionary Measures related to Phthalates:

This instruction pertains to the phthalate symbol marked on the device or its packaging. If this device is used for the treatment of children, or treatment of pregnant or nursing women; please note that the following types of procedures may increase the risk of exposure to phthalates: Exchange transfusion in neonates, total parenteral nutrition in neonates, multiple procedures in sick neonates, haemodialysis in peripubertal males, male foetus and male infant of pregnant women, and lactating women; and massive blood infusion into trauma patients. Although these procedures have the potential for increased risk of exposure, conclusive evidence of human health risks has not been established. As a precautionary measure, to reduce the potential for unnecessary exposures to phthalates, the product must be used in accordance with the instructions for use, and practitioners should refrain from using this product beyond the period of time the product is medically necessary or needed.

Dispose of all materials in accordance with local, state, and federal regulations. Decontaminate and dispose of all potentially biohazardous material.

This single-use product is not designed or validated to be reused. Reuse may cause a risk of cross-contamination, affect the measurement accuracy, system performance, or cause a malfunction as a result of the product being physically damaged due to cleaning, disinfection, re-sterilization, or reuse.

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