AirLife® baby safe flow-inflating resuscitation device

REF	Catalogue number	R_X ONLY	U.S. Federal law restricts this device to sales by or on the order of a physician
LOT	Lot Number		Manufacturer
QTY	Quantity		Single Use
	Not made with natural rubber latex		Caution
	Consult Instructions for Use	SIZE	Size
PHT DEHP	This device contains phthalates		

EN English

Product Specifications with APL Valve Closed:

35 cm H ₂ O APL Valve without Bleed Vent:	35 ± 5 cm H ₂ O at 10 LPM
35 cm H ₂ O APL Valve with Bleed Vent:	6 ± 3 cm H ₂ O at 5 LPM
	35 ± 5 cm H ₂ O at 25 LPM
60 cm H ₂ O APL Valve without Bleed Vent:	60 ± 10 cm H ₂ O at 10 LPM
60 cm H ₂ O APL Valve with Bleed Vent:	6 ± 3 cm H ₂ O at 5 LPM
	60 ± 10 cm H ₂ O at 25 LPM

Note: This infant hyperinflation system has been custom assembled to your hospital's specifications and contains either a 35 cm H₂O APL Valve (Light Blue Cap) or a 60cm H₂O APL Valve (Dark Blue Cap). Ensure that all components are securely connected prior to use.

Directions for Use:

- Attach oxygen tubing to oxygen flowmeter. Set flowmeter to 5-8 LPM.
- Pressure monitoring is recommended to ensure the delivery of appropriate airway pressure. Connect pressure gauge to port (utilize plastic sleeve provided with AirLife Disposable Manometer if necessary) or connect pressure monitoring tubing from port on fresh gas elbow to pressure gauge. If AirLife Disposable Manometer comes pre-attached, confirm secure connection. Assume that the pressure port is tightly capped when not in use.
- Occlude patient port and rotate the valve cap clockwise to completely close APL Valve. Squeeze bag and check pressure gauge to verify that the valve opens to relieve circuit pressure at the desired level (35 or 60 cm H₂O).
- Rotate the valve cap counterclockwise to the fully open position and observe that the pressure gauge drops to a maximum of 5 cm H₂O.

- Open patient port and initiate ventilation in accordance with your hospital's procedure. Adjust oxygen flow rate and APL Valve as needed to obtain desired ventilatory pressure and minute volume.

Warning:

- This device is designed for use by trained personnel.
- Do not occlude the APL Valve exhaust outlet. Occlusion may result in the generation of high airway pressures.

If product contains AirLife Disposable Manometer:

Intended Use:

The disposable pressure manometer is used to provide visual indication of a patient's airway pressure during ventilation. It may be attached to the pressure manometer port of the resuscitation bag or hyperinflation bag.

Instructions for Use:

- Locate the pressure sensing port.
- Confirm secure attachment of the disposable manometer directly to port (utilizing plastic sleeve if provided).
- Occlude the patient port on the resuscitator and pressurize the manometer to 60 cmH₂O (maximum) and release. Note that the indicator returns smoothly to zero. If indicator does not return to zero, replace the manometer.
- Connect resuscitator to artificial airway and ventilate patient.
- The "O" ring on the body of the manometer may be set as a point of reference to facilitate the desired pressure reading.
- When repeated use is required on the same patient, check the function and accuracy of the manometer prior to each use.
- When AirLife Disposable Manometer is not in use on resuscitator, the cap on the pressure sensing port must be closed.
- Discard after use.

Caution: Check the manometer prior to use, including repeated use on the same patient. A small amount of gas leaks when the manometer is in use. This is normal and should be evaluated for each patient.

Caution: This device must be used only by persons who have received adequate training in the use of airway pressure monitoring and ventilation devices.

Caution: Test circuit prior to use

Caution: During use, the patient's condition must be monitored. If for any reason the performance of the device is erratic, it must be immediately removed and replaced or the pressure port must be capped.

Manometer accuracy:

- (+/-) 1 cmH₂O from 0-10 cmH₂O (hPa)
- (+/-) 2 cmH₂O from 10-40 cmH₂O (hPa)
- (+/-) 3 cmH₂O from 40 cmH₂O (hPa)

This single-use product is not designed or validated to be reused. Reuse may cause a risk of cross-contamination, affect the measurement accuracy, system performance, or cause a malfunction as a result of the product being physically damaged due to cleaning, disinfection, re-sterilization, or reuse.

A discontinued use of this product on the same patient, with periods of use separated by a period of non-use, is acceptable. The user must ensure that the product is not damaged or contaminated between usages on the same patient.

If outer bag is labeled as containing phthalates (PHT) the following warning applies:

Risks and Precautionary Measures related to Phthalates:

This instruction pertains to the phthalate symbol marked on the device or its packaging. If this device is used for the treatment of children, or treatment of pregnant or nursing women; please note that the following types of procedures may increase the risk of exposure to phthalates: Exchange transfusion in neonates, total parenteral nutrition in neonates, multiple procedures in sick neonates, haemodialysis in peripuberal males, male foetus and male infant of pregnant women, and lactating women; and massive blood infusion into trauma patients. Although these procedures have the potential for increased risk of exposure, conclusive evidence of human health risks has not been established. As a precautionary measure, to reduce the potential for unnecessary exposures to phthalates, the product must be used in accordance with the

instructions for use, and practitioners should refrain from using this product beyond the period of time the product is medically necessary or needed.

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