AirLife™ flexible elbow

REF	Catalogue number	R_X ONLY	U.S. Federal law restricts this device to sales by or on the order of a physician
LOT	Lot Number		Legal Manufacturer
QTY	Quantity		Single Use
	Not made with natural rubber latex		Date of Manufacture
	Consult Instructions for Use		

EN

Intended Use:

The flexible elbow is intended to be used as a flexible connection between an endotracheal or tracheostomy tube and a breathing circuit which is connected to an artificial nose.

Directions for Use:

Connect flexible elbow securely to the endotracheal tube, tracheostomy tube, and/or breathing circuit (as applicable). Compensation for internal dead space of the flexible elbow may be necessary. Pressure check system to ensure leak-free connections. Always utilize appropriate alarms and visually monitor patients on life support equipment. Replace every 24 hours or as needed to prevent accumulation of secretions. If any complications are observed, such as mucous plugging, proper airway care must be instituted immediately.

Do not use a flexible elbow in the following situations:

- On patients who cannot tolerate the additional mechanical dead space of this device.
- On patients who cannot tolerate the additional airway resistance of this device.

Dispose of all materials in accordance with local, state, and federal regulations. Decontaminate and dispose of all potentially biohazardous material.

This single-use product is not designed or validated to be reused. Reuse may cause a risk of cross-contamination, affect the measurement accuracy, system performance, or cause a malfunction as a result of the product being physically damaged due to cleaning, disinfection, re-sterilization, or reuse.

AirLife Australia Holdings Pty Ltd
PO Box 97
North Ryde BC, NSW, 1670
Australia

AirLife
2710 Northridge DR. NW, Suite A
Grand Rapids, MI 49544 USA
www.myAirLife.com

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