AirLife™ filtered small volume nebulizer system

Instructions for Use

REF	Catalogue number	Rx ONLY	U.S. Federal law restricts this device to sales by or on the order of a physician
LOT	Lot Number		Manufacturer
QTY	Quantity		Date of Manufacture
	Not made with natural rubber latex		Expiration Date
	Consult Instructions for Use		Single Use
	Contains Phthalates		Caution

EN English

Intended Use: A device designed to produce an aerosol using baffling to control particle size, and project it into a gas stream.


1. Remove components from polybag.
2. Attach blue "T" piece assembly (A) to blue adaptor of the wye/filter assembly (B). Ensure opening on blue "T" piece is pointing down (C).
3. Add medication by injecting solution into nebulizer at the blue top (D). Tighten top securely.
4. Complete nebulizer assembly by attaching the blue "T" piece into the blue top of nebulizer (D).
5. Attach one end of supply tubing (E) to the prescribed gas source and the other end to the inlet at the base of the nebulizer (F).
6. Initiate flow to 5-7 LPM.

Note: Ensure all one-way valves are pointing in the proper direction for gas flow.

AU REP
AirLife Australia Holdings PTY LTD
PO Box 97
North Ryde BC, NSW, 1670
Australia

AirLife
2710 Northridge Dr. NW, Suite A
Grand Rapids, MI 49544 USA
www.myAirLife.com
Made in Mexico

AirLife®

Risks and Precautionary Measures related to Phthalates:

This instruction pertains to the phthalate symbol marked on the device or its packaging. If this device is used for the treatment of children, or treatment of pregnant or nursing women; please note that the following types of procedures may increase the risk of exposure to phthalates: Exchange transfusion in neonates, total parenteral nutrition in neonates, multiple procedures in sick neonates, haemodialysis in peripuberal males, male foetus and male infant of pregnant women, and lactating women; and massive blood infusion into trauma patients. Although these procedures have the potential for increased risk of exposure, conclusive evidence of human health risks has not been established. As a precautionary measure, to reduce the potential for unnecessary exposures to phthalates, the product must be used in accordance with the instructions for use, and practitioners should refrain from using this product beyond the period of time the product is medically necessary or needed.

Dispose of all materials in accordance with local, state, and federal regulations. Decontaminate and dispose of all potentially biohazardous material.

This single-use product is not designed or validated to be reused. Reuse may cause a risk of cross-contamination, affect the measurement accuracy, system performance, or cause a malfunction as a result of the product being physically damaged due to cleaning, disinfection, re-sterilization, or reuse.

A discontinued use of this product on the same patient, with periods of use separated by a period of non-use, is acceptable. The user must ensure that the product is not damaged or contaminated between usages on the same patient.

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