

GreenLight™ Fiber Optic Laryngoscope Handle

Instructions for Use

	Catalogue number		Manufacturer
	Lot Number		Date of Manufacture
	Quantity		Type BF equipment
	Follow Instructions For Use		Do Not Place In Trash
	Consult instructions for use		Storage Temperature
	Caution		Storage Humidity Limitation
	U.S. Federal law restricts this device to sales by or on the order of a physician		Battery Size
	Atmospheric Pressure Limitation		Use-by Date
	Restriction of Hazardous Substance (RoHS) Indicates the product is RoHS compliant		Conforms To: AAMI/IEC 60601-1 CENELEC EN60601-1 IEC 60601-1 CB Scheme IEC 60601-1-2 IEC 60601-1-6 IEC 62366-1 Certified To: CAN/CSA C22.2 No. 60601-1

EN English

INTENDED USE:

The GreenLight Fiber Optic Laryngoscope Handle is to be used for tracheal tube intubations during anesthesia, intensive care, or emergency medical airway management.


Safety Information

Only trained personnel shall use a laryngoscope for intubation.

BLADES: Compatible with all Green Standard laryngoscope handles.

For best results, use GreenLight single-patient-use blades from AirLife.

Laryngoscope handles should be tested prior to use, after replacing batteries, and after cleaning/disinfecting.

AirLife Australia Holdings PTY LTD
PO Box 97
North Ryde BC, NSW, 1670
Australia

AirLife
2710 Northridge Dr. NW, Suite A
Grand Rapids, MI 49544 USA
www.myAirLife.com
Made in Mexico

AirLife®

Medical Electrical Equipment needs special precautions regarding EMC and needs to be put into service according to the EMC information provided in the Accompanying Documents.

Portable and Mobile RF Communications Equipment can affect Medical Electrical Equipment.

Caution: Do not use this instrument in the presence of intense magnetic fields.
WARNING: No modification of this equipment is allowed.

Notes:

Wrong battery installation can lead to short circuit of the handle or improper performance.

Add Expiration Date 2 years after Date of Manufacture.

Storage Transport and Use Environment

Ambient Temperature range: 0°C–30°C

Relative Humidity range: 30%–95%

Atmospheric Pressure range 54 kPa–106 kPa

Care and Maintenance

Warning: Batteries should be removed from handles being stored for more than 7 days to avoid possible leakage of materials from the batteries.

Warning: Do not mix old and new batteries, do not mix battery brands and do not mix Lithium and Alkaline batteries.

Disposal:

Electronical and Electronic Waste (EEE) - Disposal of electronic devices in compliance with WEEE or dispose per local and state environmental regulations.

Any serious incident that has occurred in relation to the device should be reported to the manufacturer at productquality@myAirLife.com and the competent authority of the Member State in which the user and/or patient is established.

TEST PROCEDURE

1. Connect a Green Standard laryngoscope blade to the handle.
2. Pull blade to operating position.
3. Ensure that the light shines brightly. If the light appears dim, replace batteries immediately.

DISASSEMBLY AND REASSEMBLY

Disassembly:

1. Grasp handle so that the bottom is facing up.
2. Turn battery cap counterclockwise (CCW) to remove.
3. Invert handle slightly and allow the Light/Battery Cartridge to slide into other hand.

Reassembly:

1. Insert Light/Battery Cartridge into main handle and screw cap back on.
2. Perform TEST PROCEDURE.

INSTALLATION AND REPLACEMENT OF BATTERIES

To install batteries:

1. Disassemble handle.
2. Install new batteries with the negative (-) end in first.
3. Reassemble handle.
4. Perform TEST PROCEDURE.

To remove batteries:

1. Disassemble handle.
2. Insert finger or suitable object into the hole in the back of the Light/battery cartridge.
3. Grasp battery and remove.

CLEANING PROCEDURE

PRE-PROCESSING INSTRUCTIONS

1. Ensure all pre-processing instructions are followed prior to cleaning and disinfection or cleaning and sterilization.
2. Initiate cleaning and disinfection or cleaning and sterilization of device as soon as possible after use (recommended within 1 hour).
3. Transport devices via the institution's established transport procedure.
4. Remove excess gross soil as soon as possible after use by rinsing or wiping the device.
5. Light/Battery Cartridge shall be cleaned and disinfected utilizing a minimum of one (1) cleaning method and one (1) disinfection method outlined below.
6. Main Handle/Battery Cap shall be cleaned and disinfected or cleaned and sterilized utilizing a minimum of one (1) cleaning method and one (1) disinfection or sterilization method outlined below.

LIGHT/BATTERY CARTRIDGE CLEANING INSTRUCTIONS

Warning: Light/Battery Cartridge must be removed from Main Handle prior to executing the "Light/Battery Cartridge" cleaning and disinfection process.

Warning: Do not use bleach (sodium hypochlorite), betadine, or peroxide solutions. The instrument may be damaged by these solutions.

Caution: Do not soak or sterilize the Light/Battery Cartridge. Do not allow cleaning solution to enter cartridge.

LIGHT/BATTERY CARTRIDGE MANUAL CLEANING (OPTION 1)

1. Using an unfolded Cavicide® (Isopropanol = 17.2%, Ethylene Glycol Monobutyl Ether (2-Butoxyethanol) = 1-5%, Diisobutylphenoxyethoxyethyl-dimethylbenzylammonium chloride = 0.28%) branded wipe (or equivalent), remove heavy soil as needed.
2. Using an additional fresh wipe that has been unfolded, wipe the device to wet all surfaces.
3. Use additional wipes as needed ensuring the surfaces remain wet for a minimum of 2 minutes.
4. Rinse the device by wiping with disposable, lint-free cloths dampened with reverse osmosis/deionized (RO/DI) water. Repeat two (2) additional times for a total of three (3) rinse cycles, utilizing fresh RO/DI water each time.
5. Allow the device to dry.
6. Visually examine each device for cleanliness.
7. If visible soil remains, repeat cleaning procedure until the device is thoroughly clean.

LIGHT/BATTERY CARTRIDGE DISINFECTION INSTRUCTIONS

Warning: Light/Battery Cartridge must be removed from Main Handle prior to executing the "Light/Battery Cartridge" cleaning and disinfection process.

Warning: Do not use bleach (sodium hypochlorite), betadine, or peroxide solutions. The instrument may be damaged by these solutions.

Caution: Do not soak or sterilize the Light/Battery Cartridge. Do not allow disinfection solution to enter cartridge.

LIGHT/BATTERY CARTRIDGE MANUAL LOW-LEVEL DISINFECTION

1. Using an unfolded Cavicide® (Isopropanol = 17.2%, Ethylene Glycol Monobutyl Ether (2-Butoxyethanol) = 1-5%, Diisobutylphenoxyethoxyethyl-dimethylbenzylammonium chloride = 0.28%) branded wipe (or equivalent), remove heavy soil as needed.
2. Using an additional fresh wipe that has been unfolded, wipe the entire device to wet all surfaces.
3. Use additional wipes as needed ensuring the surfaces remain wet for a minimum of 2 minutes.
4. Rinse the device by wiping with disposable, lint-free cloths dampened with reverse osmosis/deionized (RO/DI) water. Repeat two (2) additional times for a total of three (3) rinse cycles, utilizing fresh RO/DI water each time.
5. Allow the device to dry.

MAIN HANDLE/BATTERY CAP CLEANING INSTRUCTIONS

Warning: Light/Battery Cartridge must be removed prior to cleaning and disinfecting or cleaning and steam sterilizing the Main Handle/Battery Cap.

Warning: Do not use bleach (sodium hypochlorite), betadine, or peroxide solutions. The instrument may be damaged by these solutions.

MAIN HANDLE/BATTERY CAP MANUAL CLEANING (OPTION 1)

1. Using a disposable, lint-free, cloth dampened with Sporidicin® (1.56% phenol and 0.06% sodium phenate based disinfectant) Brand disinfectant (or equivalent), wipe the device to remove visible soil.
2. Ensure the device remains wet with the cleaning solution for a minimum of 10 minutes.
3. Rinse the device by immersing the device in Purified Water for 1 minute. Repeat two (2) additional times for a total of three (3) rinse cycles, utilizing fresh Purified Water each time. During the rinsing process, ensure the device is agitated.
4. Allow the device to dry.
5. Visually examine each device for cleanliness.
6. If visible soil remains, repeat cleaning procedure until the device is thoroughly clean.

MAIN HANDLE/BATTERY CAP MANUAL CLEANING (OPTION 2)

1. Using an unfolded Cavicide® (Isopropanol = 17.2%, Ethylene Glycol Monobutyl Ether (2-Butoxyethanol) = 1-5%, Diisobutylphenoxyethoxyethyl-dimethylbenzylammonium chloride = 0.28%) branded wipe (or equivalent), remove heavy soil as needed.
2. Using an additional fresh wipe that has been unfolded, wipe the device to wet all surfaces.
3. Use additional wipes as needed ensuring the surfaces remain wet for a minimum of 2 minutes.
4. Rinse the device by wiping with disposable, lint-free cloths dampened with reverse osmosis/deionized (RO/DI) water. Repeat two (2) additional times for a total of three (3) rinse cycles, utilizing fresh RO/DI water each time.
5. Allow the device to dry.
6. Visually examine each device for cleanliness.
7. If visible soil remains, repeat cleaning procedure until the device is thoroughly clean.

MAIN HANDLE/BATTERY CAP STEAM STERILIZATION INSTRUCTIONS

Sterilization for United States Market
Prevacuum Steam Sterilization Cycle (U.S. FDA Compliant – Wrapped)

- Conditioning Pulses: 3
- Exposure Temperature: 132°C
- Exposure Time: 4 minutes
- Dry Time: 30 minutes
- Sterilization Configuration: FDA Cleared Sterilization Wrap (2 layer-1 ply, or 1 layer – 2 ply – examples: cellulose, polypropylene, muslin)

Sterilization for United States Market
Prevacuum Steam Sterilization Cycle – Immediate Use Sterilization (U.S. FDA Compliant – Wrapped)

- Conditioning Pulses: 3
- Exposure Temperature: 132°C
- Exposure Time: 4 minutes
- Sterilization Configuration: FDA Cleared Sterilization Wrap (2 layer-1 ply, or 1 layer – 2 ply – examples: cellulose, polypropylene, muslin)

Sterilization for International Market
Prevacuum Steam Sterilization Cycle (Wrapped)

- Conditioning Pulses: 3
- Exposure Temperature: 132°C - 135°C
- Exposure Time: 3 - 18 minutes
- Dry Time: 30 minutes
- Sterilization Configuration: Sterilization Wrap (2 layer-1 ply, or 1 layer – 2 ply – examples: cellulose, polypropylene, muslin)

Sterilization for International Market
Prevacuum Steam Sterilization Cycle – Immediate Use Sterilization (Wrapped)

- Conditioning Pulses: 3
- Exposure Temperature: 132°C - 135°C
- Exposure Time: 3 - 18 minutes
- Sterilization Configuration: Sterilization Wrap (2 layer-1 ply, or 1 layer – 2 ply – examples: cellulose, polypropylene, muslin)

MAIN HANDLE/BATTERY CAP DISINFECTION INSTRUCTIONS

Warning: Light/Battery Cartridge must be removed prior to cleaning and disinfecting or cleaning and steam sterilizing the Main Handle/Battery Cap.

Warning: Do not use bleach (sodium hypochlorite), betadine, or peroxide solutions. The instrument may be damaged by these solutions.

MAIN HANDLE/BATTERY CAP MANUAL LOW LEVEL DISINFECTION (OPTION 1)

1. Using an unfolded Cavicide® (Isopropanol = 17.2%, Ethylene Glycol Monobutyl Ether (2-Butoxyethanol) = 1-5%, Diisobutylphenoxyethoxyethyl-dimethylbenzylammonium chloride = 0.28%) branded wipe (or equivalent), remove heavy soil as needed.
2. Using an additional fresh wipe that has been unfolded, wipe the device to wet all surfaces.
3. Use additional wipes as needed ensuring the surfaces remain wet for a minimum of 2 minutes.
4. Rinse the device by wiping with disposable, lint-free cloths dampened with reverse osmosis/deionized (RO/DI) water. Repeat two (2) additional times for a total of three (3) rinse cycles, utilizing fresh RO/DI water each time.
5. Allow the device to dry.

MAIN HANDLE/BATTERY CAP MANUAL HIGH LEVEL DISINFECTION (OPTION 2)

1. Prepare the disinfectant solution CIDEX®-OPA (0.55% ortho-Phthalaldehyde based disinfectant) according to the manufacturer's directions.
2. Completely immerse the device in the disinfectant solution (0.55% ortho-Phthalaldehyde based disinfectant) and allow device to soak for a minimum of 12 minutes.
3. Ensure no air bubbles are present by agitating the device.
4. Rinse by immersing the device in purified water for 1 minute. Repeat two (2) additional times for a total of three (3) rinse cycles, utilizing fresh purified water each time. During the rinsing process, ensure the device is agitated.
5. Allow the device to dry.

MAIN HANDLE/BATTERY CAP MANUAL LOW LEVEL DISINFECTION (OPTION 3)

1. Using a disposable, lint-free, cloth dampened with Sporidicin® (1.56% phenol and 0.06% sodium phenate based disinfectant) Brand disinfectant (or equivalent), wipe the device to remove visible soil.
2. Ensure the device remains wet with the disinfectant solution for a minimum of 10 minutes.
3. Rinse by immersing the device in Purified Water for 1 minute. Repeat two (2) additional times for a total of three (3) rinse cycles, utilizing fresh Purified Water each time. During the rinsing process, ensure the device is agitated.
4. Allow the device to dry.

POST-PROCESSING INSTRUCTIONS

1. After cleaning and disinfection or cleaning and sterilization, install batteries and perform TEST PROCEDURE.

Technical Specifications

Battery Type:

Standard Handle - Size AA (2x) Alkaline (IEC-LR6) or Lithium Ion (IEC-FR6)
Stubby Handle – Size CR123 Lithium Ion (IEC-CR17345)

Equipment Classification

Type: The laryngoscope system is BF equipment per IEC60601-1
Applied Part: Patient end of laryngoscope blade

Service: Device is not serviceable.

Not made with natural rubber latex.

Accessories:

Device	Accessory	Device	Accessory
4558GSP	4559LED-LED Cartridge	4559GSP	4559LEDS-LED Cartridge
4558GSPKR	4559LEDKR-LED Cartridge	4559GSPKR	4559LEDSKR-LED Cartridge
4000 and 4600 Series GreenLight single-use laryngoscope blades			
4552VS Vital Signs laryngoscope handle sheath			

EMC:

Warning: The use of accessories, transducers and/or cables other than those specified, with the exception of those sold by the manufacturer as replacement parts for internal components, may result in increased emissions or decreased immunity of the equipment or system.

Guidance and Manufacturer's Declaration - Emissions		
The GreenLight is intended for use in the electromagnetic environment specified below. The customer or user of the GreenLight should ensure that it is used in such an environment.		
Emissions Test	Compliance	Electromagnetic Environment – Guidance
RF Emissions CISPR 11	Group 1 Class B	The GreenLight uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment

Warning: The equipment or system should not be used adjacent to or stacked with other equipment and that if adjacent or stacked use is necessary, the equipment or system should be observed to verify normal operation in the configuration in which it will be used.

Guidance and Manufacturer's Declaration – Immunity			
The GreenLight is intended for use in the electromagnetic environment specified below. The customer or user of the GreenLight should ensure that it is used in such an environment.			
Immunity Test	IEC 60601 Test Level	Compliance Level	Electromagnetic Environment – Guidance
ESD IEC 61000-4-2	±8 kV Contact ±15 kV Air	±8 kV Contact ±15 kV Air	Floors should be wood, concrete or ceramic tile. If floors are synthetic, the r/h should be at least 30%
Power Frequency 50/60Hz Magnetic Field IEC 61000-4-8	30 A/m	30 A/m	Power frequency magnetic fields should be that of a typical commercial or hospital environment.

Guidance and Manufacturer's Declaration – Emissions			
The GreenLight is intended for use in the electromagnetic environment specified below. The customer or user of the GreenLight should ensure that it is used in such an environment.			
Immunity Test	IEC 60601 Test Level	Compliance Level	Electromagnetic Environment – Guidance
Radiated RF IEC 61000-4-3	10 V/m 80 MHz – 2.7 GHz 80% Am at 1 kHz	10 V/m 80 MHz – 2.7 GHz 80% Am at 1 kHz	Home Healthcare Facility Environment

Guidance and manufacturer's declaration—electromagnetic immunity to RF wireless communications equipment. The laryngoscope is intended for use in the electromagnetic environment specified below. The customer or the user of the laryngoscope should ensure that it is used in such an environment.						
Test Frequency	Band ¹	Service ¹	Modulation ²	Maximum Power	Distance	Immunity Test Level
MHz	MHz			W	Meters	(V/m)
385	380 - 390	TETRA 400	Pulse modulation ² 18 Hz	1.8	0.3	27
450	430 - 470	GMRS 460, FRS 460	FM ³ ± 5 kHz deviation 1 kHz sine	2	0.3	28
710 745 780	704 - 787	LTE Band 13, 17	Pulse modulation ² 217 Hz	0.2	0.3	9
810 870 930	800 - 960	GSM 800/900, TETRA 800, IDEN 820, CDMA 850, LTE Band 5	Pulse modulation ² 18 Hz	2	0.3	28
1720 1845 1970	1700 - 1900	GSM 1800; CDMA 1900; GSM 1900; DECT; LTE Band 1, 3, 4, 25; UMTS	Pulse modulation ² 217 Hz	MHz	MHz	MHz
2450	2400 - 2570	Bluetooth, WLAN, 802.11 b/g/n, RFID 2450, LTE Band 7	Pulse modulation ² 217 Hz	2	0.3	28
5240 5500 5785	5100 - 5800	WLAN 802.11a/n	Pulse modulation ² 217 Hz	0.2	0.3	9
NOTE If necessary to achieve the IMMUNITY TEST LEVEL, the distance between the transmitting antenna and the ME EQUIPMENT or ME SYSTEM may be reduced to 1 m. The 1 m test distance is permitted by IEC 61000-4-3.						
¹ For some services, only the uplink frequencies are included.						
² The carrier shall be modulated using a 50 % duty cycle square wave signal.						
³ As an alternative to FM modulation, 50 % pulse modulation at 18 Hz may be used because while it does not represent actual modulation, it would be worst case.						

Recommended Separations Distances for the GreenLight

The GreenLight is intended for use in the electromagnetic environment in which radiated disturbances are controlled. The customer or user of the GreenLight can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF Communications Equipment and the GreenLight as recommended below, according to the maximum output power of the communications equipment.				
Rated maximum output power of transmitter (W)	Separation distance according to frequency of transmitter (m)			
	150 kHz to 80 MHz outside ISM bands	150 kHz to 80 MHz in ISM bands	80 MHz to 800 MHz	800 MHz to 2.5 GHz
	$d=1.2\sqrt{P}$	$d=1.2\sqrt{P}$	$d=1.2\sqrt{P}$	$d=2.3\sqrt{P}$
0.01	0.12	0.12	0.12	0.23
0.1	0.38	0.38	0.38	0.73
1	1.20	1.20	1.20	2.30
10	3.79	3.79	3.79	7.27
100	12.00	12.00	12.00	23.00
For transmitters rated at a maximum output power not listed above, the recommended separation distance d in meters (m) can be determined using the equation applicable to the frequency of the transmitter, where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer. NOTE 1 At 80 MHz and 800 MHz, the separation distance for the higher frequency range applies. NOTE 2 The ISM (industrial, scientific and medical) bands between 150 kHz and 80 MHz are 6.765 MHz to 6.795 MHz; 13.553 MHz to 13.567 MHz; 26.957 MHz to 27.283 MHz; and 40.66 MHz to 40.70 MHz. NOTE 3 An additional factor of 10/3 is used in calculating the recommended separation distance for transmitters in the ISM frequency bands between 150 kHz and 80 MHz and in the frequency range 80 MHz to 2.5 GHz to decrease the likelihood that mobile/portable communications equipment could cause interference if it is inadvertently brought into patient areas. NOTE 4 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.				

LIMITED WARRANTY

AirLife will repair or replace, free of charge, any parts of the GreenLight or Vital Signs handles, including the LED light source that fails through causes other than loss, misuse, neglect, damage in shipment, or normal wear for a period of 5 years from the date of manufacture shown on the handle.

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