
Misty Max 10™ Nebulizer

R_X ONLY

Clean, Ready to Use

0°C/
+32°F
+50°C/
+122°F
STORAGE
TEMP

AirLife®

AirLife
2710 Northridge Dr. NW, Suite A
Grand Rapids, MI 49544 USA
www.myAirLife.com

REF

002435

AirLife™**EN Misty Max 10™ Nebulizer** • With 7 foot (2.1m) Crush Resistant Oxygen Tubing, Baffled Tee Adapter (22 mm I.D./ 22 mm O.D.), Mouthpiece and 6 inch (15 cm) Flextube

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**AirLife®**

EN AirLife Misty Max 10 Nebulizer

With 7 foot (2.1m) Crush Resistant Oxygen Tubing, Baffled Tee Adapter (22 mm I.D./ 22 mm O.D.), Mouthpiece and 6 inch (15 cm) Flextube

Indications for Use

The Misty Max 10 nebulizer is used to deliver prescribed medication in aerosol form to the patient's respiratory tract.

Directions for Use

1. Unscrew nebulizer bottle from cap and insert medication into bottle, as prescribed by physician. Reattach bottle to cap. Medication may be added to nebulizer from the top of the nebulizer cap without unscrewing the two parts.
 2. Using a push-twist motion, insert the nebulizer into the base of the tee-adapter.
 3. If using with an incentive spirometer, attach 22 mm I.D. end of tee to 22 mm O.D. outlet of spirometer.
If using without an incentive spirometer, attach the 6" Flextube to the 22 mm O.D. end of the tee.
 4. Attach mouthpiece to other end of tee.
 5. Using a push-twist motion, attach the blue rigid connector on the tubing to the port of the nebulizer bottle. Attach the other end of the tubing to the gas source.
 6. Set gas flow as prescribed by a health care professional and confirm aerosol mist is flowing.
- The nebulizer is suitable for use in lung ventilator breathing systems.
 - If nebulizer is to be used in conjunction with oxygen mixtures in excess of 23%, please refer to the safe handling instructions provided by your oxygen supplier.

Note: The unit was designed to be used with or without a spirometer having a patient outlet of 22 mm O.D. It can be directly connected to a spirometer with the tee adapter or by flextube from the tee to spirometer.

Recommended Flow Rate: 5-10 liters per minute.

8 L/min Flow**	
**Measured utilizing Andersen eight-stage cascade impactor (ACI) with USP inlet, nebulizer filled with 3 mL of 0.083% albuterol sulfate in normal saline.	
MMAD	1.61
GSD	2.18
% < 4.7 µm	85.2%

- Using a solution, suspension, or emulsion different from that recommended by the manufacturer can alter the particle size distribution curve, the mass median aerodynamic diameter (MMAD), aerosol output, and/or aerosol output rate, which can then be different from those disclosed by the manufacturer
- The disclosures listed for nebulizer performance are based upon testing that utilizes adult ventilator patterns and are likely to be different from those stated for the pediatric or infant populations.

Cleaning Instructions

Unscrew the nebulizer cap and bottle. Remove the one-piece jet by pulling and twisting the jet off the jet stem. Wash all components in warm soapy water and rinse well. Air dry or hand dry with a clean, lint free cloth. Reseat the jet by placing the jet over the jet stem and snap into place. Reattach the nebulizer cap and bottle. Device will withstand 50 cleaning cycles.

Caution

Ensure the bottle is reattached to the cap after cleaning

Storage

Recommended normal storage temperature: 0°C/+32°F to +50°C/+122°F.

Contains Phthalates (DEHP)

Risks and Precautionary Measures related to Phthalates: This instruction pertains to the phthalate symbol marked on the device or its packaging. If this device is used for the treatment of children, or treatment of pregnant or nursing women; please note that the following types of procedures may increase the risk of exposure to phthalates: Exchange transfusion in neonates, total parenteral nutrition in neonates, multiple procedures in sick neonates, haemodialysis in peripuberal males, male fetus and male infant of pregnant women, and lactating women; and massive blood infusion into trauma patients. Although these procedures have the potential for increased risk of exposure, conclusive evidence of human health risks has not been established. As a precautionary measure, to reduce the potential for unnecessary exposures to phthalates, the product must be used in accordance with the instructions for use, and practitioners should refrain from using this product beyond the period of time the product is medically necessary or needed. Dispose of all materials in accordance with local, state, and federal regulations. Decontaminate and dispose of all potentially biohazardous material.

Specifications

- Recommended Driving Gas: Oxygen or Compressed Air
- Recommended Driving Gas Flow: 5-10 L/min
- Maximum Recommended Fill Volume: 10 mL
- Maximum A-weighted Sound Pressure Level = 63.0 dB
- Empty Mass of Nebulizer = 0.015 kg

Symbols Glossary

	Catalogue number		U.S. Federal law restricts this device to sales by or on the order of a physician
	Lot Number		Clean, Ready to Use
	Consult instructions for use		Quantity
	Not made with natural rubber latex		Expiration Date
	Manufacturer		Caution
	Contains Phthalates		Do not re-use
	Temperature Limit		

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EN Recommended Nebulizer Set Flow Rates 5 - 10 L/min				
Drug Description	Aerosol Characteristic ¹	5 L/min	8 L/min	10 L/min
Albuterol Sulfate Dose 3 mg/3 mL (0.1% by volume) or 3000 µg of drug substance	MMAD (µm) ²	5.5 ± 0.2	3.5 ± 0.2	2.7 ± 0.2
	GSD ³	2.0 ± 0.1	2.1 ± 0.1	2.2 ± 0.1
	% Particles (<2 µm) ⁴	13 ± 2	24 ± 2	35 ± 4
	% Particles (2 to 5 µm) ⁴	28 ± 2	40 ± 1	40 ± 1
	% Particles (>5 µm) ⁴	60 ± 2	36 ± 3	25 ± 3
	% Respirable Particles (<5 µm) ⁴	40 ± 2	64 ± 3	75 ± 3
	Aerosol Output Rate [µg/min]	98 ± 6	136 ± 9	132 ± 7
	Total Delivered Mass of Aerosol [µg] ⁵	648 ± 17	553 ± 17	530 ± 15
	% Fill Volume Emitted per minute	3.3 ± 0.2	4.5 ± 0.3	4.4 ± 0.2
	% Fill Volume Emitted ⁶	21.6 ± 0.01	18.4 ± 0.01	17.7 ± 0.01
	Residual Volume at Sputter + 1 minute [mL] ⁷	1.3 ± 0.03	1.5 ± 0.05	1.5 ± 0.03
	Driving Gas Pressure [psig]	12 psig	26 psig	36 psig

¹ Performance data using Cascade Impactor (NGI) with an extraction flow of 15 L/min and a representative adult breathing patient for 0.1% Albuterol Sulfate solution measured according to BS ISO 27247:2013. Treatment time was defined as onset of audible sputtering plus one minute. Values represent the Mean and 95% confidence interval of the mean.

² MMAD is Mass Median Aerodynamic Diameter.

³ GSD is the Geometric Standard Deviation.

⁴ Values are the mass of drug substance as a percentage of the Total Delivered Mass of Aerosol.

⁵ The total delivered mass of aerosol is the mass of aerosol emitted by the nebulizing system using a representative adult breathing patient for the given fill volume (i.e. 3mL) that is emitted by the nebulizer for a treatment time of sputter plus one minute

⁶ The percentage of fill volume emitted is the total delivered mass of aerosol expressed as a percentage of the fill volume (i.e. 3mL) that is emitted by the nebulizer for a treatment time of sputter plus one minute.

⁷ Residual Volume is the volume of liquid remaining in the nebulizer system for a treatment time of sputter plus one minute.

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