
AirLife™

EN Misty Max 10™ Nebulizer Medication System with Filter

The Misty Max 10 nebulizer is used to deliver prescribed medication in aerosol form to the patient's respiratory tract.

Contains

Misty Max 10 Nebulizer
 2.1 m (7') Crush Resistant Oxygen Tubing
 Patient "Y"
 Mouthpiece
 Flextube
 Flow Directional Valves
 Nonconductive Respiratory Filter (Bacterial/Viral Retentive)[†]

[†]Bacterial Filtration Efficiency (B.F.E.) 99.9% when tested with aerosols of *S. aureus* and Viral Filtration Efficiency (V.F.E.) 99.7% when tested with aerosols of *E. coli* bacteriophage (virus) to the method Mil-M-36954C specification.

Specifications

Minimum Recommended Fill Volume: 1.5 ml
 Maximum Recommended Fill Volume: 10 ml
 Recommended Driving Gas: Oxygen or Compressed Air
 Minimum Driving Gas Flow: 5 l/min @ 11 PSIG (76 kPa)
 Recommended Driving Gas Flow: 8 l/min @ 26 PSIG (179 kPa)
 Maximum Driving Gas Flow: 10 l/min @ 35 PSIG (241 kPa)
 MMAD Range* 1.0 - 1.4 µm
 Aerosol Output Range* 26 - 1856 µl NaF
 Aerosol Output Rate* 13 - 72 µl NaF min-1
 Sound Pressure Level* 62 dB

*Measured according to the EN13544-1 Standard

Directions for Use

1. Remove unit from package.
2. Attach mouthpiece to patient end of "Y."
3. Ensure that flexible discs on one-way valves are oriented per the diagram and that components are securely attached.
4. Unscrew nebulizer bottle from cap and insert medication into bottle, as prescribed by physician. Reattach bottle to cap. Medication may be added to nebulizer from the top of the nebulizer cap without unscrewing the two parts.
5. Insert nebulizer into tee adapter.
6. Using a push-twist motion, attach the blue rigid connector on the tubing to the port of the nebulizer bottle. Connect the other end of the tubing to gas source.
7. Instruct patient to exhale through the mouthpiece.
8. Set gas flow as recommended by health care professional.
 - Performance information provided per EN13544 may not apply to drugs supplied in suspension or high viscosity form. In such cases please seek information from the drug manufacturer.
 - If nebulizer is to be used in conjunction with oxygen mixtures in excess of 23%, please refer to the safe handling instructions provided by your oxygen supplier.

Warnings

Do not soak, rinse, wash, or sterilize. Soaking, rinsing, washing, steam, cold or gas sterilization of this device may result in either retention of harmful residual gases in the filter, or dangerously increase its breathing resistance. Do not connect directly to an endotracheal tube.

Storage

Recommended normal storage temperature 0°C/+32°F to +50°C/+122°F.

Contains Phthalates (DEHP). The benefit of treatment with medical respiratory support devices outweighs the remote possibility of exposure to phthalates.

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