Vital Signs® JACKSON-REES CIRCUITS

REF	Catalogue number		U.S. Federal law restricts this device to sales by or on the order of a physician
LOT	Lot Number		
QTY	Quantity		Manufacturer
	Not made with natural rubber latex		Single Use
	Consult Instructions for Use		This device contains Phthalates
	Caution		

EN English

Directions for Use:

1. Test circuit for occlusions and confirm flow to patient elbow prior to use.
2. The use of a manometer (pressure gauge) is recommended to measure delivery pressures.
3. If the device includes exhaust valve, adjust valve and oxygen flow rate for desired patient pressures.
4. If device does not include exhaust valve, user MUST supply exhaust.

Warning: Device designed for use by trained personnel. If specific model excludes exhaust valve, use of a pop off safety valve is recommended.

Do not rinse, soak, wash or sterilize.

Risks and Precautionary Measures related to Phthalates:

This instruction pertains to the phthalate symbol marked on the device or its packaging. If this device is used for the treatment of children, or treatment of pregnant or nursing women; please note that the following types of procedures may increase the risk of exposure to phthalates: Exchange transfusion in neonates, total parenteral nutrition in neonates, multiple procedures in sick neonates, haemodialysis in peripubertal males, male foetus and male infant of pregnant women, and lactating women; and massive blood infusion into trauma patients. Although these procedures have the potential for increased risk of exposure, conclusive evidence of human health risks has not been established. As a precautionary measure, to reduce the potential for unnecessary exposures to phthalates, the product must be used in accordance with the instructions for use, and practitioners should refrain from using this product beyond the period of time the product is medically necessary or needed.

AirLife Australia Holdings Pty Ltd
PO Box 97
North Ryde BC, NSW, 1670
Australia

AirLife
2710 Northridge DR. NW, Suite A
Grand Rapids, MI 49544 USA
www.myAirLife.com

Single Use

This single-use product is not designed or validated to be reused. Reuse may cause a risk of crosscontamination, affect the measurement accuracy, system performance, or cause a malfunction as a result of the product being physically damaged due to cleaning, disinfection, re-sterilization, or reuse.

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