Vital Signs™ Fiber Optic Laryngoscope Handle

Instructions for Use

REF

Catalogue number

Manufacturer

LOT

Lot Number

Date of Manufacture

QTY

Quantity

Do Not Place In Trash

Follow Instructions For Use

Storage Temperature

Consult instructions for use

Storage Humidity Limitation

Caution

Battery Size

U.S. Federal law restricts this device to sales by or on the order of a physician

Use-by Date

Atmospheric Pressure Limitation

Type BF equipment

✓ **RoHS**

Restriction of Hazardous Substance (RoHS) Indicates the product is RoHS compliant

Conforms To:
AAMI/IEC 60601-1
CENELEC EN60601-1
IEC 60601-1 CB Scheme
IEC 60601-1-2
IEC 60601-1-6
IEC 62366-1

Certified To:
CAN/CSA C22.2 No. 60601-1

EN

English

INTENDED USE:

The Vital Signs Fiber Optic Laryngoscope Handle is to be used for tracheal tube intubations during anesthesia, intensive care, or emergency medical airway management.

Safety Information

Only trained personnel shall use a laryngoscope for intubation. **BLADES:** For exclusive use with Vital Signs single-patient-use fiber optic laryngoscope blades, manufactured by Vital Signs, Inc. Not compatible with ISO Green standard blade. Laryngoscope handles should be tested prior to use, after replacing batteries, and after cleaning/disinfecting. Medical Electrical Equipment needs special precautions regarding EMC and needs to be put into service according to the EMC information provided in the Accompanying Documents. Portable and Mobile RF Communications Equipment can affect Medical Electrical Equipment.

Caution: Do not use this instrument in the presence of intense magnetic fields.
WARNING: No modification of this equipment is allowed.

AU REP AirLife Australia Holdings Pty Ltd
PO Box 97
North Ryde BC, NSW, 1670
Australia

AirLife
2710 Northridge Dr. NW, Suite A
Grand Rapids, MI 49544 USA
www.myAirLife.com

Notes:

Wrong battery installation can lead to short circuit of the handle or improper performance.

Light/ battery cartridge must be removed from main handle prior to executing the "Light/ Battery Cartridge" Cleaning and disinfection process.

Add Expiration Date 2 years after Date of Manufacture.

Storage Transport and Use Environment

Ambient Temperature range: 0°C–30°C
Relative Humidity range: 30%–95%
Atmospheric Pressure range 54 kPa–106 kPa

Care and Maintenance

Warning: Batteries should be removed from handles being stored for more than 7 days to avoid possible leakage of materials from the batteries
Warning: Do not mix old and new batteries, do not mix battery brands and do not mix Lithium and Alkaline batteries.

Disposal:

Electronical and Electronic Waste (EEE) - Disposal of electronic devices in compliance with WEEE and dispose per local and state environmental regulations.

Any serious incident that has occurred in relation to the device should be reported to the manufacturer at productquality@myAirLife.com and the competent authority of the Member State in which the user and/or patient is established.

TEST PROCEDURE

1. Connect a Vital Signs fiber optic laryngoscope blade to the handle.
2. Pull blade to operating position.
3. Ensure that the light shines brightly. If the light appears dim, replace batteries immediately.

INSTALLATION AND REPLACEMENT OF BATTERIES

1. Grasp handle so that the bottom is facing up.
2. Turn battery cap counterclockwise (CCW) to remove.
3. Invert handle slightly and allow the batteries to slide into other hand.
4. Install new batteries with the positive (+) end in first into the main handle and screw cap back on.
5. Perform TEST PROCEDURE.

CLEANING PROCEDURE

Warning: Current cleaning, disinfection and/or sterilization techniques will not inactivate prions, e.g. new variants of Creutzfeld-Jacob disease (CJD).
Warning: Do not use bleach (sodium hypochlorite), betadine, or peroxide solutions. The instrument may be damaged by these solutions.
Caution: Do not soak or sterilize device. Do not allow solution to enter into the handle. Do not autoclave.

Handles may be wiped with cloth dampened with 70% isopropyl alcohol or Sporicidin® Brand disinfectant.

Upon completion of all cleaning procedures, install batteries and perform TEST PROCEDURE.

Important: After use, dispose of the blade. Vital Signs blades are for single use only. Do not reuse.

Technical Specifications

Battery Type:
Standard Handle - Size AA (2x) Alkaline (IEC-LR6) or Lithium Ion (IEC-FR6)
Stubby Handle – Size CR123 Lithium Ion (IEC-CR17345)
Service: Device is not serviceable.
Not made with natural rubber latex.

Equipment Classification

Type: The laryngoscope system is BF equipment per IEC 60601-1
Applied Part: Patient end of laryngoscope blade

Accessories:

Vital Signs single-use fiber optic laryngoscope blades
Vital Signs laryngoscope handle sheath

EMC:

Warning: The use of accessories, transducers and/or cables other than those specified, with the exception of those sold by the manufacturer as replacement parts for internal components, may result in increased emissions or decreased immunity of the equipment or system.

Guidance and Manufacturer's Declaration - Emissions		
The Vital Signs laryngoscope handle is intended for use in the electromagnetic environment specified below. The customer or user of the Vital Signs laryngoscope handle should ensure that it is used in such an environment.		
Emissions Test	Compliance	Electromagnetic Environment – Guidance
RF Emissions CISPR 11	Group 1 Class B	The Vital Signs laryngoscope handle uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment

Warning: The equipment or system should not be used adjacent to or stacked with other equipment and that if adjacent or stacked use is necessary, the equipment or system should be observed to verify normal operation in the configuration in which it will be used.

Guidance and Manufacturer's Declaration – Immunity			
The Vital Signs laryngoscope handle is intended for use in the electromagnetic environment specified below. The customer or user of the Vital Signs laryngoscope handle should ensure that it is used in such an environment.			
Immunity Test	IEC 60601 Test Level	Compliance Level	Electromagnetic Environment – Guidance
ESD IEC 61000-4-2	±8 kV Contact ±15 kV Air	±8 kV Contact ±15 kV Air	Floors should be wood, concrete or ceramic tile. If floors are synthetic, the r/h should be at least 30%
Power Frequency 50/60Hz Magnetic Field IEC 61000-4-8	30 A/m	30 A/m	Power frequency magnetic fields should be that of a typical commercial or hospital environment.

Guidance and Manufacturer's Declaration – Emissions			
The Vital Signs laryngoscope handle is intended for use in the electromagnetic environment specified below. The customer or user of the Vital Signs laryngoscope handle should ensure that it is used in such an environment.			
Immunity Test	IEC 60601 Test Level	Compliance Level	Electromagnetic Environment – Guidance
Radiated RF IEC 61000-4-3	10 V/m 80 MHz – 2.7 GHz 80% Am at 1 kHz	10 V/m 80 MHz – 2.7 GHz 80% Am at 1 kHz	Home Healthcare Facility Environment

Guidance and manufacturer's declaration—electromagnetic immunity to RF wireless communications equipment. The laryngoscope is intended for use in the electromagnetic environment specified below. The customer or the user of the laryngoscope should ensure that it is used in such an environment.						
Test Frequency	Band ¹	Service ¹	Modulation ²	Maximum Power	Distance	Immunity Test Level
				W	Meters	(V/m)
385	380 - 390	TETRA 400	Pulse modulation ² 18 Hz	1.8	0.3	27
450	430 - 470	GMRS 460, FRS 460	FM ³ ± 5 kHz deviation 1 kHz sine	2	0.3	28
710 745 780	704 - 787	LTE Band 13, 17	Pulse modulation ² 217 Hz	0.2	0.3	9
810 870 930	800 - 960	GSM 800/900, TETRA 800, IDEN 820, CDMA 850, LTE Band 5	Pulse modulation ² 18 Hz	2	0.3	28
1720 1845 1970	1700 - 1900	GSM 1800; CDMA 1900; GSM 1900; DECT; LTE Band 1, 3, 4, 25; UMTS	Pulse modulation ² 217 Hz	MHz	MHz	MHz
2450	2400 - 2570	Bluetooth, WLAN, 802.11 b/g/n, RFID 2450, LTE Band 7	Pulse modulation ² 217 Hz	2	0.3	28
5240 5500 5785	5100 - 5800	WLAN 802.11a/n	Pulse modulation ² 217 Hz	0.2	0.3	9
NOTE If necessary to achieve the IMMUNITY TEST LEVEL, the distance between the transmitting antenna and the ME EQUIPMENT or ME SYSTEM may be reduced to 1 m. The 1 m test distance is permitted by IEC 61000-4-3.						
¹ For some services, only the uplink frequencies are included.						
² The carrier shall be modulated using a 50 % duty cycle square wave signal.						
³ As an alternative to FM modulation, 50 % pulse modulation at 18 Hz may be used because while it does not represent actual modulation, it would be worst case.						

Recommended Separations Distances for the Vital Signs Laryngoscope handle				
The Vital Signs laryngoscope handle is intended for use in the electromagnetic environment in which radiated disturbances are controlled. The customer or user of the Vital Signs laryngoscope handle can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF Communications Equipment and the Vital Signs laryngoscope handle as recommended below, according to the maximum output power of the communications equipment.				
Rated maximum output power of transmitter (W)	Separation distance according to frequency of transmitter (m)			
	150 kHz to 80 MHz outside ISM bands	150 kHz to 80 MHz in ISM bands	80 MHz to 800 MHz	800 MHz to 2.5 GHz
	$d=1.2\sqrt{P}$	$d=1.2\sqrt{P}$	$d=1.2\sqrt{P}$	$d=2.3\sqrt{P}$
0.01	0.12	0.12	0.12	0.23
0.1	0.38	0.38	0.38	0.73
1	1.20	1.20	1.20	2.30
10	3.79	3.79	3.79	7.27

100	12.00	12.00	12.00	23.00
For transmitters rated at a maximum output power not listed above, the recommended separation distance d in meters (m) can be determined using the equation applicable to the frequency of the transmitter, where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.				
NOTE 1 At 80 MHz and 800 MHz, the separation distance for the higher frequency range applies.				
NOTE 2 The ISM (industrial, scientific and medical) bands between 150 kHz and 80 MHz are 6.765 MHz to 6.795 MHz; 13.553 MHz to 13.567 MHz; 26.957 MHz to 27.283 MHz; and 40.66 MHz to 40.70 MHz.				
NOTE 3 An additional factor of 10/3 is used in calculating the recommended separation distance for transmitters in the ISM frequency bands between 150 kHz and 80 MHz and in the frequency range 80 MHz to 2.5 GHz to decrease the likelihood that mobile/portable communications equipment could cause interference if it is inadvertently brought into patient areas.				
NOTE 4 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.				

LIMITED WARRANTY

Vital Signs will repair or replace, free of charge, any parts of the GreenLight or Vital Signs handles, including the LED light source that fails through causes other than loss, misuse, neglect, damage in shipment, or normal wear for a period of 5 years from the date of manufacture shown on the handle.

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