Vital Signs™ flow-inflating tail-end valve transport circuit

REF	Catalogue number		U.S. Federal law restricts this device to sales by or on the order of a physician
LOT	Lot Number		Manufacturer
SIZE	Size		Date of Manufacture
QTY	Quantity		Single Use
	Not made with natural rubber latex		Caution
	Consult Instructions for Use		This device contains Phthalates.

EN English

Warning:

1. Occlusion of the exhaust outlet/valve may result in the generation of high airway pressures.

Instructions:

1. Connect oxygen tubing to the flow meter.
2. Test the system prior to use for leaks and occlusions by pressurizing the circuit to 30 cmH₂O for thirty seconds.
3. Set the fresh gas flow rate at two to three times minute volume to avoid re-breathing of Carbon dioxide.
4. Use of a pressure manometer is recommended to measure airway delivery pressures.
5. If the device includes an exhaust valve, adjust the valve and flow rate to achieve the desired tidal volume and patient airway pressure.
6. If the device does not include an exhaust valve, the user must supply exhaust.

If product contains AirLife™ Disposable Manometer:

Intended Use:

The disposable pressure manometer is used to provide visual indication of a patient's airway pressure during ventilation. It may be attached to the pressure manometer port of the resuscitation bag or hyperinflation bag.

Instructions for Use:

1. Locate the pressure sensing port.
2. Confirm secure attachment of the disposable manometer directly to port (utilizing plastic sleeve if provided).
3. Occlude the patient port on the resuscitator and pressurize the manometer to 60 cmH₂O (maximum) and release. Note that the indicator returns smoothly to zero. If indicator does not return to zero, replace the manometer.
4. Connect resuscitator to artificial airway and ventilate patient.
5. The "O" ring on the body of the manometer may be set as a point of reference to facilitate the desired pressure reading.
6. When repeated use is required on the same patient, check the function and accuracy of the manometer prior to each use.
7. When AirLife Disposable Manometer is not in use on resuscitator, the cap on the pressure sensing port must be closed.
8. Discard after use.

Manometer accuracy:

(+/-) 1 cmH₂O from 0-10 cmH₂O (hPa)
(+/-) 2 cmH₂O from 10-40 cmH₂O (hPa)
(+/-) 3 cmH₂O from 40 cmH₂O (hPa)

Caution: Do not use in the presence of flammable anesthetics. Prior to use, check pressure relief is 6 cmH₂O +/- 1 cmH₂O under static inflow of 5 liters per minute.

Caution: Check the manometer prior to use, including repeated use on the same patient. A small amount of gas leaks when the manometer is in use. This is normal and should be evaluated for each patient.

Caution: This device must be used only by persons who have received adequate training in the use of airway pressure monitoring and ventilation devices.

Caution: During use, the patient's condition must be monitored. If for any reason the performance of the device is erratic, it must be immediately removed and replaced or the pressure port must be capped.

Dispose of all materials in accordance with local, state, and federal regulations. Decontaminate and dispose of all potentially biohazardous material.

This single-use product is not designed or validated to be reused. Reuse may cause a risk of cross-contamination, affect the measurement accuracy, system performance, or cause a malfunction as a result of the product being physically damaged due to cleaning, disinfection, re-sterilization, or reuse.

A discontinued use of this product on the same patient, with periods of use separated by a period of non-use, is acceptable. The user must ensure that the product is not damaged or contaminated between usages on the same patient.

If outer bag is labeled as containing phthalates (PHT) the following warning applies:

Risks and Precautionary Measures related to Phthalates:

This instruction pertains to the phthalate symbol marked on the device or its packaging. If this device is used for the treatment of children, or treatment of pregnant or nursing women; please note that the following types of procedures may increase the risk of exposure to phthalates: Exchange transfusion in neonates, total parenteral nutrition in neonates, multiple procedures in sick neonates, haemodialysis in peripubertal males, male foetus and male infant of pregnant women, and lactating women; and massive blood infusion into trauma patients. Although these procedures have the potential for increased risk of exposure, conclusive evidence of human health risks has not been established. As a precautionary measure, to reduce the potential for unnecessary exposures to phthalates, the product must be used in accordance

with the instructions for use, and practitioners should refrain from using this product beyond the period of time the product is medically necessary or needed.

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