
REF AH165





AirLife™ Infant Heated Wire Circuit

dual-limb, dual-heat, high-flow circuit (>4 L/min)

- ES Circuitos de cables calentados para bebés AirLife**
circuito de gran flujo de dos fuentes de calor y de dos miembros (>4 l/min)
- FR Circuits pédiatriques chauffés par fil AirLife**
Circuit à débit élevé à double branche et à double source de chauffage (>4 l/min)
- PT Circuitos de Fios Aquecidos Infantis AirLife**
circuito de membro duplo, calor duplo, fluxo elevado (>4 L/min)
- DE AirLife beheizte Atemkreislauf für Kleinkinder**
Zweigliedriger, doppelt beheizter Kreis mit hohem Durchfluss (>4 L/min)
- IT Circuiti riscaldati a filo per lattanti AirLife**
Circuito a flusso elevato, due rami, doppio riscaldamento (>4 l/min)
- SV AirLife uppvärmda kretsar för spädbarn**
Krets med dubbla slangar, dubbla värmeslingor, högt flöde (>4 l/min)
- NL AirLife verwarmd buiscircuit voor zuigelingen**
Dubbel-onderdeel, dubbel-verwarmd, hogestroming-circuit (>4 l/min)
- DA AirLife opvarmet wirekredsløb til spædbørn**
Kredsløb med dobbelt slange, dobbelt opvarmning, højt flow (>4 l/min)

The logo for AirLife, featuring a stylized arch above the word "AirLife" with a registered trademark symbol (®).

EN AirLife Infant Heated Wire Circuits

Indications for Use:

The heated breathing circuit is intended to deliver and warm breathing gases before they enter the patient's airway. The heated breathing circuit is used with Fisher and Paykel heated MR850 humidifier. The AirLife Infant Heated Wire Circuit is used with pediatric population, specifically the neonate (birth to 1 month of age) and infant (greater than 1 month to 2 years of age) pediatric subgroups. The product is single use device, non-sterile and used in professional healthcare environments and intra-hospital transport environments under a doctor's supervision and by skilled clinicians. The AirLife Heated Wire Circuit is used for flow rates greater than 4 LPM.

Note: This breathing circuit is a compatible accessory to the Fisher & Paykel MR850 humidifier.

Setup:

1. **INSPIRATORY LIMB (purple):** Connect the inspiratory limb to the humidifier chamber outlet port. Plug appropriate heater wire adapter into mating port on inspiratory limb.
2. **EXPIRATORY LIMB (white, dual-limb only):** Connect the expiratory limb to the exhalation port of the breathing machine. Plug appropriate heater wire adapter into mating port on expiratory limb.
3. **HUMIDIFIER LIMB:** The humidifier limb connects the breathing machine outlet port to the humidifier chamber inlet port.
4. **PRESSURE SENSING LINE:** The pressure sensing line connects the pressure port on the patient connector to the pressure port on the breathing machine. If the line is not required, use the appropriate plug to seal the pressure port on the patient connector.
5. Install temperature probe(s) into the temperature probe port(s), removing any plugs prior to inserting the probe. For application with the MR850, ensure that the wedge on the chamber probe is aligned with, and fully inserted into, the notch on the temperature probe port.
6. **UNHEATED EXTENSION (tubing void of heated wire following the temperature probe—removable):** The unheated extension is designed for use in an incubator setting. The temperature probe in the inspiratory limb should sit outside of the incubator during care. If excessive condensation occurs in the unheated extension it can be removed. The unheated extension is removed when caring for the infant in an open care situation. The unheated extension connects between the inspiratory limb and the patient.

Test:

- Ensure that all connections are tight. Some connections are deliberately assembled only finger tight and must be firmly secured.
- Test circuit prior to use by occluding the patient connection port and pressure test the circuit to ensure that there are no leaks. Also, check that there are no occlusions by allowing the gas to flow and ensure that the gas is emitted from the patient connection port. Block the patient connection port and ensure that gas is flowing through the expiratory limb (dual-limb only).
- Check that the heated wire is evenly distributed along the circuit and not bunched or kinked.

Cautions:

- Do not use this circuit where gas temperature at the outlet of the humidifier exceeds 68 degrees Celsius / 154 degrees Fahrenheit.
- Do not use this circuit with humidifiers or heater wire adapters other than those specified in the applications section of the product label.
- Do not reconfigure circuits into single limb configurations.
- Do not use the heated wire circuit without gas flow. Heated circuits must have a minimum gas flow (please refer to product label) through the tubing at all times. Operating circuits with lower flow rates (minute volume + background / bias flow) may result in reduced patient humidification and increased condensation in the circuit tubing.
- Do not place material on or around the heated wire tubing. Objects such as heavy tapes, towels, or bed linens may over insulate the circuits, impede normal heat convection, and cause damage to the tubing or interruption of gas delivery to the patient.
- Re-use may degrade the performance of the product or contribute to cross contamination.
- Place circuit below patient level to prevent condensate from draining into the patient's airway.

Warnings:

- Avoid contact with patient's skin.
- Do not stretch the tubing.
- Do not soak, rinse, wash, re-use, or sterilize this product.
- Ensure the circuit is gas tight and that there are no leaks.
- This product is intended to be used for a maximum of 30 days
- Only for use with the F&P MR850 humidifier

Patient Population:

Newborn (Neonate)	From birth to 1 month of age
Infant	Greater than 1 month to 2 years of age

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AirLife Australia Holdings Pty Ltd