
REF AH180

AirLife™ Adult Heated Wire Circuit

Dual-limb, Dual-heated circuit

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EN AirLife Adult Heated Wire Circuits

Indications for Use:

The AirLife Adult Heated Wire Circuit is intended to warm breathing gases before they enter a patient's airway. The AirLife Adult Heated Wire Circuit is used with the adult patient population that requires mechanical ventilation, positive pressure breathing or general medical gases. The product is a single use device, non-sterile and used in professional healthcare environments and intra-hospital transport environments under a doctor's supervision and by skilled clinicians. The AirLife Adult Heated Wire Circuit is compatible to the Fisher & Paykel MR850 humidifier.

Setup:

1. **INSPIRATORY LIMB (purple):** Connect the inspiratory limb to the humidifier chamber outlet port. Plug connector into appropriate heated wire adapter.
2. **EXPIRATORY LIMB (white, dual-limb only):** Connect the expiratory limb to the exhalation port of the breathing machine. Plug appropriate heated wire adapter into port.
3. **HUMIDIFIER LIMB:** The humidifier limb connects the breathing machine outlet port to the humidifier chamber inlet port.
4. Install temperature probe(s) into the temperature probe port(s), removing any plugs prior to inserting the probe. For application with the MR850, ensure that the wedge on the chamber probe is aligned with, and fully inserted into, the notch on the temperature probe port.

Test:

- Ensure that all connections are tight. Some connections are deliberately assembled only finger tight and must be firmly secured.
- Test circuit prior to use by occluding the patient connection port and pressure test the circuit to ensure that there are no leaks. Also, check that there are no occlusions by allowing the gas to flow and ensure that the gas is emitted from the patient connection port. Block the patient connection port and ensure that gas is flowing through the expiratory limb (dual-limb only).
- Check that the heater wire is evenly distributed along the circuit and not bunched or kinked.

Cautions:

- Do not use this circuit where gas temperature at the outlet of the humidifier exceeds 68 degrees Celsius / 154 degrees Fahrenheit.
- Do not use this circuit with humidifiers or heated wire adapters other than those specified in the applications section of the product label.
- Do not reconfigure circuits into single limb configurations.
- Do not use the heated wire circuit without gas flow. Heated circuits must have a minimum gas flow (please refer to product label) through the tubing at all times. Operating circuits with lower flow rates (minute volume + background / bias flow) may result in reduced patient humidification and increased condensation in the circuit tubing.
- Do not place material on or around the heated wire tubing. Objects such as heavy tapes, towels, or bed linens may over insulate the circuits, impede normal heat convection, and cause damage to the tubing or interruption of gas delivery to the patient.
- Re-use may degrade the performance of the product or contribute to cross contamination.
- The heated wire circuit has been tested at an altitude of ≤2000 m under environmental conditions of 21-25 degrees Celsius ambient temperature and 20-60 percent relative humidity. Operation outside of these ranges may vary results.
- After use properly dispose of product in accordance with local hospital protocols and all applicable Federal and State regulations.
- The heated wire circuit and heater base are Medical Electrical Equipment that need special precautions regarding EMC and need to be installed and put into service according to the EMC information provided in the accompanying documents.
- Portable and mobile RF communications equipment can affect Medical Electrical Equipment.
- Essential Performance of the heated wire circuit and heater base is to prevent drying out of airway secretions in patients whose upper airways have been bypassed. During therapy measured gas temperatures can fluctuate around set gas temperatures without affecting essential performance. In a single fault condition of exposure to electrical fast transient/bursts up to ±1kV, an audible alarm will sound and power to the heater core will be lost. Measured gas temperatures can temporarily fluctuate but do not affect essential performance.

Warnings:

- Avoid prolonged contact with patient's skin.
- Do not stretch the tubing.
- Do not soak, rinse, wash, re-use, or sterilize this product.
- Ensure the circuit is gas tight and that there are no leaks.
- This product is intended to be used for a maximum of 30 days.

Application

REF No.	Resistance to Flow	Compliance	Rated Flow (EN12342)
AH180 (Dual Limb)	R(I) @ 30 l/min: 0.02 hPa/l/min (cmH ₂ O/l/min) R(E) @ 30 l/min: 0.02 hPa/l/min (cmH ₂ O/l/min)	C @ 60 hPa: 2.2ml/hPa (ml/cmH ₂ O)	(I) 30 l min ⁻¹ , ≤.002 kPa (E) 30 l min ⁻¹ , ≤.002 kPa

REF No.	Humidifier	Wire Adapter	Temperature/ Flow Probe	Flow Rate	Rated Flow	Length	Maximum Operating Pressure	Minimum Inside Diameter
AH180	F&P MR850	900MR806	900MR869	>3 L/min	30 L/min	1.52 m	8 kPa	19.15 mm

- The use of Wire Adapters or Temperature/Flow Probes other than those specified may result in increased emissions or decreased electromagnetic immunity of the heated wire circuit and MR850 heater base.
- The heated wire circuit and heater base should not be used adjacent to or stacked with other equipment with the exception of a mechanical ventilator. If adjacent or stacked use is necessary, check to ensure that no part of the heater base or its cables will be near other AC cables. The heater base should be observed to verify normal operation in the configuration in which it will be used.

Symbols Glossary

	Catalogue number		Manufacturer
	Lot Number		U.S. Federal law restricts this device to sales by or on the order of a physician
	Not made with DEHP		Quantity
	Not made with natural rubber latex		Use-by Date
	Read manual		Date of Manufacture
	For Single use only		Caution

DUAL LIMB CHART

- (A) 22mm CONNECTOR TO PATIENT
- (B) TEMPERATURE PROBE PORT TO PATIENT
- (C) INSPIRATORY LIMB
- (D) WIRE CLIP
- (E) TEMPERATURE PROBE PORT TO CHAMBER
- (F) 22mm CONNECTOR TO CHAMBER
- (G) 22mm CONNECTOR TO VENTILATOR
- (H) DRY LINE
- (I) PATIENT ELBOW
- (J) EXPIRATORY LIMB
- (K) DUAL LIMB HANGAR
- (L) PATIENT WYE

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AirLife Australia Holdings Pty Ltd