REF 001618

AirLife™

EN Universal IPPB Manifold and Nebulizer (Quick-Neb)

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
EN Universal IPPB Manifold and Nebulizer (Quick-Neb)

Rx ONLY Clean, Ready to Use

With 39" (99 cm), 22 mm ID Corrugated Tubing, Two Supply Lines, Mouthpiece and 6" (15 cm), 22 mm ID Flextube

Features of this unit:

- Tapered tubing connectors on body.
- Removable anti-spill plate.
- High output nebulizer.
- Removable gas collection head with 360° swivel action.
- Accepts all tubing from 19 mm to 22 mm ID.
- Polyethylene-EVA tubing is standard on all setups.
- Crack resistant nebulizer bottle with screw threads.
- Drive line adapter for use with single supply line respirators.
- For use with both single and double supply line respirators.

Directions for use:

1. Single Supply Line Setup: Figure (A) represents the proper assembly for manifolds using a single supply line. Attach reducing valve with 18 cm (7") tubing to stem on top of nebulizer bottle. Attach free end of 18 cm (7") tubing to gas collection head stem. Attach a supply line to reducing valve stem. Connect the other end to drive line located on the respirator. Included is an extra drive line adapter for use on those machines that do not include an external adapter. Place the ribbed section into tube and then insert in machine as needed.
2. Double Supply Line Setup: Figure (B) represents the proper assembly for manifolds using two supply lines. Attach nebulizer line to stem on top of nebulizer bottle and the other end to nebulizer drive stem on the respirator. Adjust nebulizer flow rate before beginning treatment.
3. General Information: Manifolds are equipped with a mouthpiece and 15 cm (6") flextube. If monitoring is requested, spirometer attachments will fit directly over exhalation port. Any 22 mm ID connector will fit as well. To use support stem, unscrew bracket on support arm, place ball on stem into support and tighten.
4. Filling of Nebulizer: To fill the nebulizer, unscrew the bottle from the manifold assembly. Fill the nebulizer bottle with medication as prescribed by physician's orders. Replace bottle on the manifold assembly.
5. Removal of Cap and Diaphragm: Figure (C) represents the method for removal of the gas collection cover and diaphragm. Do not lift straight up. Lift at a tilted angle by grasping the exhalation port of the gas collection cover. Remove diaphragm cap for access to diaphragm. If the machine stops or a leak develops in the diaphragm, remove the cover and cap, Figure (D), and reseal the diaphragm. If the unit continues to malfunction, replace the unit. If nebulizer malfunctions, replace the unit. Do not use pins to open jet holes.
6. Assure that all hose connectors are tight.

Caution: All respirator setups should be made in accordance with a physician's orders. Assure that the nebulizer bottle is screwed on tightly or exhalation failure may result.

Contains Phthalates (DEHP). The benefit of treatment with medical respiratory support devices outweighs the remote possibility of exposure to phthalates.

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The logo features a stylized arch above the word "AirLife" in a serif font, followed by a registered trademark symbol (®).
AirLife Australia Holdings PTY LTD