[EN] AirLife™ EdithTrach, disposable

[REF]	557005000	AirLife™ Edith-Trach HME for Spontaneously Breathing Patients, Disposable			
[REF]					
[LOT]	Lot number				
[QTY]	Quantity				
[MFG]	Manufacturer				
[EXP]	Expiry date				


U.S. Federal law restricts this device to sale by or on the order of a physician.

Intended use

EdithTrach is a single use heat and moisture exchanger which transfers heat and humidity from the gas expired by the patient to the inspired gases. EdithTrach is intended for use on patients who are breathing spontaneously through a tracheostomy tube.

Material
Housing: Polypropylene (PP)
HME element: Plastic

Not made with natural rubber latex.

Instructions for use
Read these instructions completely before using the product.
a) Cover lid with suction possibility
b) Oxygen port
c) Patient side

Connect EdithTrach to the tracheostomy tube.
For oxygen therapy connect an oxygen tube to the oxygen port of EdithTrach. EdithTrach shall be changed every 24 hours or more frequently as required.
Shelf life: Max. 3 years. Expiration date printed on package.

Contraindications
EdithTrach shall not be used on patients producing fulminating frothy secretions within their airways and lungs.
EdithTrach shall not be used on patients with very small tidal volumes, for example neonates or very small children.

Precautions
Before use, verify that EdithTrach has no occlusions and that air will flow through it.
During the use of EdithTrach, the patient shall be closely monitored and proper airway care administered if complications arise.
During oxygen therapy, make sure that EdithTrach is positioned in a way that avoids kinking of oxygen tube.

Warning
Single-use products are not designed to be reused. Reuse may cause a risk of cross-contamination, affect the measurement accuracy and/or system performance, or cause a malfunction as a result of the product being physically damaged due to cleaning, disinfection, re-sterilization and/or reuse.

Technical data		
Minimum Vt (ml)	60	
Maximum Vt (ml)	1000	
Weight (g)	6	
Dead space (ml)	16	
Pressure drop, kPa (cm H ₂ O)		
- at 60 l/min	0.02 (0.2)	
Moisture to patient (mg H ₂ O)¹	24.1²	
Moisture loss (mg H ₂ O)¹	13.5²	

[CS]		
Technické údaje		
Minimum Vt (ml)	60	
Maximum Vt (ml)	1000	
Weight (g)	6	
Dead space (ml)	16	
Pressure drop, kPa (cm H ₂ O)		
- at 60 l/min	0.02 (0.2)	
Moisture to patient (mg H ₂ O)¹	24.1²	
Moisture loss (mg H ₂ O)¹	13.5²	
Connector		
- patient side	15 F	

¹ Measured at Vt 500 ml
² Tested with dry gas (not room air).

Ordering information
For further information, please contact your local sales representative or visit www.myAirLife.com.

[FI]		
Technical data		
Minimum Vt (ml)	60	
Maximum Vt (ml)	1000	
Weight (g)	6	
Dead space (ml)	16	
Pressure drop, kPa (cm H ₂ O)		
- at 60 l/min	0.02 (0.2)	
Moisture to patient (mg H ₂ O)¹	24.1²	
Moisture loss (mg H ₂ O)¹	13.5²	
Connector		
- patient side	15 F	

¹ Measured at Vt 500 ml
² Tested with dry gas (not room air).

Ordering information
For further information, please contact your local sales representative or visit www.myAirLife.com.

[FR]		
Technical data		
Minimum Vt (ml)	60	
Maximum Vt (ml)	1000	
Weight (g)	6	
Dead space (ml)	16	
Pressure drop, kPa (cm H ₂ O)		
- at 60 l/min	0.02 (0.2)	
Moisture to patient (mg H ₂ O)¹	24.1²	
Moisture loss (mg H ₂ O)¹	13.5²	
Connector		
- patient side	15 F	

¹ Measured at Vt 500 ml
² Tested with dry gas (not room air).

Ordering information
For further information, please contact your local sales representative or visit www.myAirLife.com.

[BG] AirLife™ EdithTrach, disposable

[REF]	557005000	AirLife™ Edith-Trach HME for Spontaneously Breathing Patients, Disposable			
[REF]					
[LOT]	Lot number				
[QTY]	Quantity				
[MFG]	Manufacturer				
[EXP]	Expiry date				

U.S. Federal law restricts this device to sale by or on the order of a physician.

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Before use, verify that EdithTrach has no occlusions and that air will flow through it.
During the use of EdithTrach, the patient shall be closely monitored and proper airway care administered if complications arise.
During oxygen therapy, make sure that EdithTrach is positioned in a way that avoids kinking of oxygen tube.

Warning
Single-use products are not designed to be reused. Reuse may cause a risk of cross-contamination, affect the measurement accuracy and/or system performance, or cause a malfunction as a result of the product being physically damaged due to cleaning, disinfection, re-sterilization and/or reuse.

[DE]		
Technische Daten		
Minimum Vt (ml)	60	
Maximum Vt (ml)	1000	
Weight (g)	6	
Dead space (ml)	16	
Pressure drop, kPa (cm H ₂ O)		
- at 60 l/min	0.02 (0.2)	
Moisture to patient (mg H ₂ O)¹	24.1²	
Moisture loss (mg H ₂ O)¹	13.5²	
Connector		
- patient side	15 F	

[ES]		
Características técnicas		
Minimum Vt (ml)	60	
Maximum Vt (ml)	1000	
Weight (g)	6	
Dead space (ml)	16	
Pressure drop, kPa (cm H ₂ O)		
- at 60 l/min	0.02 (0.2)	
Moisture to patient (mg H ₂ O)¹	24.1²	
Moisture loss (mg H ₂ O)¹	13.5²	
Connector		
- patient side	15 F	

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