EN	Vital-Signs™ ECG Leadwire Sets 300-Series
 	
REF	Catalogue number
545315-HEL	Vital Signs™ 300-Series ECG Leadwire, 5-lead, Grabber, IEC, 75 cm 30 in
545316-HEL	Vital Signs™ 300-Series ECG Leadwire, 5-lead, Grabber, IEC, 125 cm 41 in & 75 cm 30 in
545317-HEL	Vital Signs™ 300-Series ECG Leadwire, 5-lead, Grabber, AHA, 75 cm 30 in
545318	Vital Signs™ 300-Series ECG Leadwire, 5-lead, Grabber, AHA, 125 cm 41 in & 75 cm 30 in
545327-HEL	Vital Signs™ 300-Series ECG Leadwire, 5-lead, Snap, AHA, 75 cm 30 in
545328	Vital Signs™ 300-Series ECG Leadwire, 5-lead, Snap, AHA, 125 cm 41 in & 75 cm 30 in

REF	Catalogue number	Non-sterile
LOT	Lot number	Medical device
QTY	Quantity	Follow instructions for use
Manufacturer	Manufacturer	Temperature limit
Date of manufacture	Date of manufacture	Not made with natural rubber latex.
MR unsafe	UDI	Unique device identifier
Waste electrical & electronic equipment	 	U.S. Federal law restricts this device to sale by or on the order of a physician.

Intended use
Use leadwire sets and cables™ with Datex-Ohmeda® GE Healthcare monitors or modules with a blue front panel connector either for measuring ECG or impedance supporting that feature. The cables are only compatible with 300 series lead sets.

*REF: 545300-HEL, 545301-HEL, 545302-HEL, 545303-HEL

Instructions for use	EKG Cables and leadwire sets
Defibrillator cables and cables to be used with 300 Series ECG Leadwire Sets with clips or snaps.	
ECG Cable accessories	
- Leadwire Set Combiner keeps all individual lead wires together.	
- Leadwire Set Organizer keeps leads in order.	
- Bedsheet Clips attach the cap cable to the bed sheet.	

These products are indicated for use by qualified medical personnel only. Please refer to the monitor user documentation for detailed information.

NOTE: The 3-lead set can be connected to the 5-lead cable to measure a 12-lead ECG.

Any serious incident that has occurred in relation to the device should be reported to the manufacturer at productquality@myairlife.com and the competent authority of the Member State in which the user and/or patient is established.	
Maintenance	
Before use, check that the product is intact and clean. Manufacturer has successfully tested the enclosed reprocessing instructions for 80 cycles	
Warnings	
- MR unsafe - Do not use this product in MRI environment.	
- Do not let the conductive product come into contact with any other conductive parts used during earth.	
- Ensure proper contact of the electrocuregny unit neutral electrode to the patient to avoid burns and skin irritation. Do not place electrodes too close to the electrocuregny unit neutral electrode or electrocuregny cables.	
- Do not allow humidity or liquid to enter the product.	
- Do not use product if electrical damage before use. Do not use product if you see any visual signs of damage. A damaged product should be discarded immediately. Never try to repair a damaged product, do modifications to product or use a product repaired by others.	
- Products are supplied non-sterile and must be unpacked, cleaned and disinfected prior to first use. Products must be reprocessed prior to each use. Thorough cleaning is a necessary requirement for subsequent effective disinfection.	
- Only use equipment and product-specific validated procedures for cleaning and disinfection.	
- The used equipment and accessories must be maintained and checked regularly according to the manufacturer's instructions.	
- Follow applicable laws, standards and guidelines for reprocessing must be followed.	
- The use of personal protective equipment is the responsibility of the user.	
- Manufacturer is not responsible for damages resulting from improper use.	
- Manufacturer is not responsible for damage to the user's equipment to comply with the validated parameters during every reprocessing.	
- Use any accessory combinations and monitor as mentioned in this IFU for cleaning and disinfection of the product.	
- Do not use the product if the product is damaged or if the product is used to occur as they might cause airway obstruction, irritation, decubitus, falling of other items or people, etc.	

Pre-Processing Instructions	
1. Ensure all pre-processing instructions are followed prior to cleaning and disinfecting.	
2. Thorough cleaning and disinfection of device as soon as possible after use (recommended within 1 hour).	
3. Transport devices using the institution's established transport procedure.	
4. Remove excess gross soil as soon as possible after use by wiping the device.	
5. Device shall be cleaned and disinfected utilizing a minimum of one (1) cleaning method and one (1) disinfection method outlined below.	

Cleaning Instructions	
1. Ensure all pre-processing instructions are followed prior to cleaning.	
2. Prepare a cleaning solution of pH detergent, utilizing tap water containing a temperature range of 27 °C to 44 °C (81 °F to 111 °F), per manufacturer's instructions.	
3. Use a soft, clean cloth saturated with the cleaning solution to manually clean all surfaces of the devices including cracks and crevices for a minimum of 5 minutes. In conjunction with this step, use an appropiately sized brush (e.g. standard soft bristle toothbrush or lumen brush) to disinfected hard-to-reach areas paying particular attention to cracks and crevices ensuring a slightly wet surface is maintained at room temperature for the duration of this disinfection step.	
4. Using a sterile, soft cloth saturated with purified water (PURW) at 27°C – 44°C, manually rinse by wiping the device for 30 seconds. Repeat two (2) additional times for a total of three (3) rinse cycles, utilizing fresh purified water each time.	
5. Drying: Dry the device with a sterile, lint-free cloth.	

Intermediate level disinfection instructions (option 2)	
1. Prepare a 70% IPA disinfectant solution according to the manufacturer's directions.	
2. Using a sterile, gauze saturated with the disinfection solution, manually disinfect all contaminated surfaces of the device for a minimum of 15 minutes.	
3. In conjunction with this step, use an appropriately sized brush (e.g. standard soft bristle toothbrush or lumen brush) to disinfected hard-to-reach areas paying particular attention to cracks and crevices ensuring a slightly wet surface is maintained at room temperature for the duration of this disinfection step.	
4. Using a sterile, gauze saturated with purified water (PURW) at 27°C – 44°C, manually rinse by wiping the device for 30 seconds. Repeat two (2) additional times for a total of three (3) rinse cycles, utilizing fresh purified water each time.	
5. Drying: Dry the device using a sterile, lint-free cloth.	

Disposal	
Reprocess products prior to disposal. Dispose electrical and electronic devices in compliance with local and state environmental regulations. This product does not contain PVC.	
Technical description	
- This product is classified as defibrillation-proof Type CF applied part.	
- The product meets the requirements of the standard, <i>ECG trunk cables and patient leads/wires (ANSI/AAMI EC53)</i> . This product is intended for multi-patient use.	
- Not made with natural rubber latex.	

Ordering information	
For further information, please contact your local sales representative or visit www.mylife.com .	

Intended use	
Use leadwire sets and cables™ with Datex-Ohmeda® GE Healthcare monitors or modules with a blue front panel connector either for measuring ECG or impedance supporting that feature. The cables are only compatible with 300 series lead sets.	
*REF: 545300-HEL, 545301-HEL, 545302-HEL, 545303-HEL	

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These products are indicated for use by qualified medical personnel only. Please refer to the monitor user documentation for detailed information.

NOTE: The 3-lead set can be connected to the 5-lead cable to measure a 12-lead ECG.

Any serious incident that has occurred in relation to the device should be reported to the manufacturer at productquality@myairlife.com and the competent authority of the Member State in which the user and/or patient is established.	
Maintenance	
Before use, check that the product is intact and clean. Manufacturer has successfully tested the enclosed reprocessing instructions for 80 cycles	
Warnings	
- MR unsafe - Do not use this product in MRI environment.	
- Do not let the conductive product come into contact with any other conductive parts used during earth.	
- Ensure proper contact of the electrocuregny unit neutral electrode to the patient to avoid burns and skin irritation. Do not place electrodes too close to the electrocuregny unit neutral electrode or electrocuregny cables.	
- Do not allow humidity or liquid to enter the product.	
- Do not use product if electrical damage before use. Do not use product if you see any visual signs of damage. A damaged product should be discarded immediately. Never try to repair a damaged product, do modifications to product or use a product repaired by others.	
- Products are supplied non-sterile and must be unpacked, cleaned and disinfected prior to first use. Products must be reprocessed prior to each use. Thorough cleaning is a necessary requirement for subsequent effective disinfection.	
- Only use equipment and product-specific validated procedures for cleaning and disinfection.	
- The used equipment and accessories must be maintained and checked regularly according to the manufacturer's instructions.	
- Follow applicable laws, standards and guidelines for reprocessing must be followed.	
- The use of personal protective equipment is the responsibility of the user.	
- Manufacturer is not responsible for damages resulting from improper use.	
- Manufacturer is not responsible for damage to the user's equipment to comply with the validated parameters during every reprocessing.	
- Use any accessory combinations and monitor as mentioned in this IFU for cleaning and disinfection of the product.	
- Do not use the product if the product is damaged or if the product is used to occur as they might cause airway obstruction, irritation, decubitus, falling of other items or people, etc.	

Pre-Processing Instructions	
1. Ensure all pre-processing instructions are followed prior to cleaning and disinfecting.	
2. Thorough cleaning and disinfection of device as soon as possible after use (recommended within 1 hour).	
3. Transport devices using the institution's established transport procedure.	
4. Remove excess gross soil as soon as possible after use by wiping the device.	
5. Device shall be cleaned and disinfected utilizing a minimum of one (1) cleaning method and one (1) disinfection method outlined below.	

Cleaning Instructions	
1. Ensure all pre-processing instructions are followed prior to cleaning.	
2. Prepare a cleaning solution of pH detergent, utilizing tap water containing a temperature range of 27 °C to 44 °C (81 °F to 111 °F), per manufacturer's instructions.	
3. Use a soft, clean cloth saturated with the cleaning solution to manually clean all surfaces of the devices including cracks and crevices for a minimum of 5 minutes. In conjunction with this step, use an appropiately sized brush (e.g. standard soft bristle toothbrush or lumen brush) to disinfected hard-to-reach areas paying particular attention to cracks, crevices and other hard to clean areas.	
4. Using a sterile, soft cloth saturated with purified water (PURW) at 27°C – 44°C, manually rinse by wiping the device for 30 seconds. Repeat two (2) additional times for a total of three (3) total rinse cycles.	
5. Drying: Dry the device with a soft, sterile cloth.	
6. Visually examine each device for cleanliness.	
7. If visible soil remains, repeat cleaning procedure until the device is thoroughly clean.	

Disinfection	
Intermediate level disinfection instructions (option 1)	
1. Prepare a 2% Glutaraldehyde disinfectant solution (CIDEX® or equivalent) according to the manufacturer's directions.	

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