EN	GB	CS	DA	DE	EL	ES	ET
Vital Signs™ Leg Attachment Pad Instructions for Use REF 2464AAO-20 Leg Attachment Pad							
REF Catalogue number	MR unsafe						
MD Medical device	Not made with natural rubber latex						
QTY Quantity							
QTY Quantity	Single use						
Exp. date	Temperature limit						
Follow instructions for use	U.S. Federal law restricts this device to sale by or to the order of a physician.						
LOT Lot number							
Manufacturer	UDI Unique device identifier						
Consult instruction for use							

Intended use
Leg Attachment Pad is intended for use in conjunction with Fetal Spiral Electrode cables (REF 1590AAO and 1591AAO) with the GE Comenetics™ monitor and Fetal Spiral Electrode for Comenetics™ (REF 7000AAO) in the healthcare setting. The Attachment Pad is used as a neutral electrode on expectant mothers for internal fetal heart rate monitoring of their fetus during labor.

This product is intended to be used by a healthcare professional in the hospital setting (labor/ward unit, free-standing birthing centers, emergency department).

Instructions for use
1. Prepare Fetal Spiral Electrode (REF 7000AAO) and Fetal Spiral Electrode Cable (REF 1590AAO or 1591AAO) for fetal heart rate monitoring, as advised.
2. Prepare patient's skin prior to placing Leg Attachment Pad.
3. Take Leg Attachment Pad from its pouch.
4. Connect Fetal Spiral Electrode Cable (1590AAO or 1591AAO) to the Leg Attachment Pad.
5. Remove release liner from the back of the Leg Attachment Pad.
6. Secure the Leg Attachment Pad on the skin of the patient's leg/low abdomen.
7. Start monitoring fetal heart rate.
8. After monitoring, remove the Leg Attachment Pad from the patient's leg/low abdomen.
9. After use, dispose Leg Attachment Pad.

Any serious incident that has occurred in relation to the device should be reported to the manufacturer at productquality@myairlife.com and the competent authority of the Member State in which the user and/or patient is established.

Warnings
- MR unsafe. Do not use this product in MRI environment.
- Not made with natural rubber latex.
- Do not reuse, sterilize or clean the product. Reuse may cause a risk of cross-contamination, affect the measurement accuracy and/or system performance, or cause a malfunction as a result of the product being physically damaged due to cleaning, disinfection, re-sterilization and/or reuse.
- Reuse of this product is against the intended use.
- Deviation from the defined procedures within this instruction is in the responsibility of the user.
- Discard removed release liner appropriately. Release liner may cause oral injury obstruction.

Cautions
- Check on patient initially and occasionally over time and avoid prolonged use to prevent pain, pressure, swelling, hematoma, bruising, etc. from occurring.
- To ensure appropriate and safe functioning of the product and system, check security of attachment to patient initially and occasionally over time.
- Inspect the product for physical damage before use. Do not use product if you see any visual signs of damage. A damaged product should be discarded immediately. Never try to repair a damaged product, do modifications to product or use a product repaired by others.

Storage
Storage Temperature +5.0°C/+41.0°F to +32.0°C/+89.6°F.
Disposal
Dispose per local and state environmental regulations.

Technical description
- Not made with natural rubber latex.
- Do not use more than 24 hours.
- Leg Attachment Pad is part of FEGC applied part, classified as Type BF.

Ordering information
For further information, please contact your local sales representative or visit www.myAirLife.com.

	
AirLife Finland Oy Kortaneentie 2 FI-00510 Helsinki, Finland www.myAirLife.com	

Assembled in the US, with US and foreign components.
Customer Service: productquality@myairlife.com
Australian Sponsor: Austrial Australia Holdings Pty Ltd.
PO BOX 97
NORTH RYDE BC NSW 1670 Australia

100007386 Rev A 2024/10	
-------------------------	--

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