Medisorb™ EX, Disposable

Instructions for Use

REF 2079796-001, Medisorb™ EX, Disposable

REF	Catalogue number		Consult instructions for use
LOT	Lot number		Follow instructions for use
QTY	Quantity		Keep dry
	Manufacturer		Fragile, handle with care
	Date of manufacture		This side up
	Expiry date		Keep away from sunlight
	Temperature limit		Danger

R_X ONLY U.S. Federal law restricts this device to sale by or on the order of a physician.

Intended use

The Medisorb™ EX (later absorber) is a disposable product intended for use with the GE Healthcare Carestation 600 series.

The absorber should only be used with air, oxygen, nitrous oxide, halothane, enflurane, isoflurane, desflurane and sevoflurane.

Warning

- The device is intended to be used under constant attention of qualified professional healthcare personnel.
- The absorber is not an antimicrobiological filter or inhibitor. Based on this information, it is the responsibility of the clinician to take appropriate measures regarding patient cross contamination when required.
- Vital Signs™ recommends that the Vital Signs filter is used in order to minimize the cross contamination and the risk of foreign particles reaching the patient.
- The absorber is a sealed unit, and should not be opened for any reason. No modification is allowed.

Technical description

Materials

Casing: SBC (Styrene-butadiene copolymer)
Grids and Protective seal: PP (Polypropylene)

Soda lime

Soda lime: Medisorb
Granule size: 2–5 mm
Moisture content: 12-19%
12<pH<14
Sodium hydroxide (NaOH, 215-185-5): <4%
Calcium hydroxide (Ca(OH)₂, 215-137-3): >75%

Weight

Total weight: 1500 g
Soda lime weight: 1200 g

Volume

Canister: 1900 ml
Max. soda lime: 1400 ml
Resistance (at 60 l/min)
insp: <0.53 kPa / exp: <0.55 kPa

Instructions for use

1. Read these instructions carefully before using the product.
2. Observe that when handling the absorber channeling of the soda lime might occur. The channeling can reduce the CO₂ absorption capacity of the product resulting in early CO₂ breakthrough. Channeling can be removed by tapping the canister gently on a hand or on a table 3-4 times in the orientation in which it will be installed (see Figure 1).


Figure 1

3. Before using the absorber, remove the seal.
4. Make sure the connecting surfaces are free of dust, the absorber is properly installed and the breathing circuit is checked for leakage and occlusion prior to use.

Connecting to breathing system

1. Slide the absorber canister on the lifter dish (see Figure 2).

2. Pull up the handle to lock the canister to the breathing system (see Figure 3).


Figure 2


Figure 3

Replacing during ventilation

The Carestation 600 series have a CO₂ bypass system allowing replacement of the absorber during ventilation without disturbing the gas flow to the patient.

Warning

- Removal of the absorber while the anesthesia system is delivering gas to a patient allows carbon dioxide to build up in the breathing system.

Removing from breathing system

1. Pull down the handle to unlock the canister from the breathing system.
2. Lift the canister and remove it from the lifter dish.

When to change the absorber

At absorber exhaustion:

A gradual color change from white to violet of the soda lime indicates absorption of carbon dioxide. The color change is only a rough indicator and its intensity may vary from one procedure to another. When the absorber reaches exhaustion of its absorption capacity, there may be some amount of soda lime in the bottom of the absorber that does not change color. This is a natural occurrence.

Operating room personnel must determine the remaining absorptive capacity of the absorber immediately at the end of a surgical procedure.

A typical carbon dioxide absorption capacity of the absorber is 12 hours for inspired CO₂ to reach 0.5% with a minute ventilation of 5 l/min, (a respiratory rate of 10 breaths/min and a tidal volume of 500 ml) a fresh gas flow of 500 ml/min and a carbon dioxide output from the patient of 190 ml/min. Hours of use depend on fresh gas flow and patient carbon dioxide output.

Warnings

- The color change of the soda lime is only a rough indicator. Use carbon dioxide monitoring to determine when to change the absorber.
- If not in use for several hours, soda lime may regain its original color giving a misleading indication of effectiveness.

At absorber desiccation:

Absorber should be change if desiccated.

Warnings

- Excessive use of oxygen flush induces dehydration of the absorber and inhibits its ability to absorb carbon dioxide.
- If the absorber is left mounted overnight, make sure that all gas flows are stopped. Otherwise the absorber may dry out and lose its activity.
- If the absorber dries out completely, it may give off CO (carbon monoxide) when exposed to anesthetic agents.
- For safety, always replace the absorber when its state of desiccation is unknown.

Soda lime

The following safety messages apply to the soda lime:

Warnings

- Hazard and Precautionary statements of the soda lime: H314, Causes severe skin burns and eye damage. P260, Do not breathe dust/fume/gas/mist/vapours/spray. P280, Wear protective gloves/protective clothing/eye protection/face protection.
- P305/P351/P338, IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing.
- P303/361/353, IF ON SKIN (or hair): Remove/Take off immediately all contaminated clothing. Rinse skin with water/shower.
- P310, Immediately call a POISON CENTER or doctor/physician.
- Avoid contact of contents of absorber with eyes, skin or clothing. Avoid creating and breathing airborne dust that may cause irritation.
- Inhaled anesthetic agents in contact with strong bases are reported to be able to degrade to produce breakdown products of unknown toxicity in humans. For further information, contact the manufacturer of the inhaled anesthetic agent in question.
- The absorber must not be used with chloroform or trichloroethylene. The soda lime may react with chloroform/trichloroethylene, producing sodium

formate/dichloroethylene, carbon monoxide and/or phosgene, which are harmful for the patient.

Caution

- For information about the soda lime, refer to the appropriate Safety Data Sheet.

Use of Accessories

The following safety messages apply to the use of accessories:

Warnings

- It is recommended to use water traps in the patient circuit during low flow anesthesia (fresh gas flow less than 1 l/min) and when excessive condensed water accumulates in the limbs. Use water traps in both expiratory and inspiratory limbs to remove water condensation.

- Use only Vital Signs validated supplies and accessories approved by AirLife for proper performance of the anesthesia system. Using products other than those approved may result in improper performance of the anesthesia system.

Disposal

Disposal of a used absorber must be determined by the end user in compliance with federal state, provincial and local regulations. Specific disposal methods vary with the chemical (e.g., anesthetic agent, acid-base vapors, etc.) or biological agents used with the soda lime by the end user.

Warnings

- Keep all absorbers that have been used with highly flammable anesthetics away from heat, sparks and open flames, as residual amounts of these materials will be present and may cause fire hazard.

Caution

- Local authority regulations must be investigated before disposal of the product.

Storage

Unopened absorber package must remain unopened and be stored in a dry and clean environment at a temperature from 0°C/+32°F to +35°C/+95°F. Avoid direct sunlight (UV light). Avoid freezing which may reduce CO₂ absorption performance. Once the absorber package is opened, its contents should be used immediately as the performance will deteriorate if storage is continued. The package must be protected from physical damage and water. Under these conditions, the absorber will retain its carbon dioxide absorption capacity for two years. The expiration date of the absorber is printed on the product.

Ordering information

For further information, please contact your local sales representative or visit myAirLife.com.

AirLife Finland Oy
Kuortaneenkatu 2
FI-00510 Helsinki, Finland
www.myAirLife.com

Made in UK

Australian Sponsor
Airlife Australia Holdings Pty Ltd.
PO BOX 97
NORTH RYDE BC NSW 1670
Australia

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