Medisorb™ Pre-Packed Rounds, Disposable

Instructions for Use

REF 427000100 Medisorb™ Pre-Packed Rounds, Disposable

REF	Catalogue number		Consult instructions for use
LOT	Lot number		Follow instructions for use
QTY	Quantity		Keep dry
	Manufacturer		Fragile, handle with care
	Date of manufacture		This side up
	Expiry date		Keep away from sunlight
	Temperature limit		Danger
R_X ONLY	U.S. Federal law restricts this device to sale by or on the order of a physician.		

Intended use

The Medisorb™ Pre-Packed Round is a disposable carbon dioxide absorbent round product intended for use with anesthesia machines. The Pre-Packed Round should only be used with air, oxygen, nitrous oxide, halothane, enflurane, isoflurane, desflurane and sevoflurane. The device is intended to be used under constant attention of qualified professional healthcare personnel.

Technical description

Dimensions

Weight: 1.1 kg

Materials

Casing: Polystyrene

Volume

Canister 1250 ml

12<pH<14

Sodium hydroxide (NaOH, 215-185-5): <4%

Calcium hydroxide (Ca(OH)₂, 215-137-3): >75%

Absorbent

Granule size: 2.5–5 mm

Moisture content: 12-19%

Instructions for use

Read these instructions carefully before using the product.

Preparing the Pre-Packed Round

- Ensure that the washers on the round holder are correctly located and that the sealing surfaces are free from particles.
- Remove the protective plastic shrink-wrap from the Pre-Packed Round.
- Place Pre-Packed Round in the holder with the arrow pointing upwards.
- Tighten the holder in accordance with the holder manufacturer's instructions.
- Check the breathing circuit for leakage prior to use. Make sure the connecting surfaces are free of dust, the Pre-Packed Round is properly installed, and the breathing circuit is checked for leakage and occlusion prior to use.

Using the Pre-Packed Round

Warning

It is contraindicated to remove the Pre-Packed Round from the machine during an anesthetic procedure. Removal of the Pre-Packed Round while the anesthesia system is delivering gas to a patient may cause leak and/or may allow carbon dioxide to build up in the breathing system.

When to change Pre-Packed Round

The absorption of carbon dioxide is indicated by a gradual color change of the soda lime. The intensity of the color change may vary from one procedure to another. The color change of the soda lime is only a rough indicator. Use carbon dioxide monitoring to determine when to change the Pre-Packed Round. Color change is reversible and the exhausted material will change back to its original color. Saturated Pre-Packed Round should therefore be disposed of immediately after use.

Operating room personnel must determine the remaining absorptive capacity of the Pre-Packed Round immediately at the end of the surgical procedure. Excessive use of oxygen flush induces dehydration of the soda lime and inhibits its ability to absorb carbon dioxide.

There may be some amount of absorbent in the bottom of the absorber that does not change color; this is a natural occurrence.

It is recommended to use water traps in the patient circuit during low flow anesthesia (fresh gas flow less than 1 l/min) and when excessive condensed water accumulates in the limbs. Use water traps in both expiratory and inspiratory limbs to remove water condensation.

Maintenance

Disposal

Saturated Pre-Packed Round should be disposed of immediately after use. Keep all Pre-Packed Rounds that have been used with highly flammable anesthetics away from heat, sparks and open flames, as residual amounts of these materials will be present and may cause a fire hazard. Disposal of a used Pre-Packed Round must be determined by the end user in compliance with federal state, provincial and local regulations. Specific disposal methods vary with the chemical (e.g., anesthetic agent, acid-base vapors, etc.) or biological agents used with the absorbent by the end user.

Storage

The Pre-Packed Round package must be closed and stored in a dry and clean environment at a temperature from 0°C/+32°F to +35°C/+95°F. Avoid direct sunlight (UV light). Avoid freezing which may reduce CO₂ absorption performance. Once the Pre-Packed Round protective plastic shrink-wrap has been removed contents should be used immediately as the performance will deteriorate if storage is continued. The package must be protected from physical damage and water. Under these conditions, the Pre-Packed Round will retain its carbon dioxide absorption capacity for two years. The expiration date of the Pre-Packed Round is printed on the product.

Warnings

- Hazard and Precautionary statements of the soda lime: H314, Causes severe skin burns and eye damage. P260, Do not breathe dust/fume/gas/mist/vapours/spray. P280, Wear protective gloves/protective clothing/eye protection/face protection.
- P305/P351/P338, IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing.
- P303/P361/P353, IF ON SKIN (or hair): Remove/Take off immediately all contaminated clothing. Rinse skin with water/shower.
- P310, Immediately call a POISON CENTER or doctor/physician.
- Constant attention by a qualified professional is needed whenever this product is used.
- The Pre-Packed Round is a sealed unit and should not be opened for any reason. No modification is allowed.
- Avoid contact of contents of Pre-Packed Round with eyes, skin or clothing. Avoid creating and breathing airborne dust that may cause irritation. In the event of skin or eye contact, immediately rinse the affected area with water and seek medical assistance.
- The color change of the soda lime is only a rough indicator. Use carbon dioxide monitoring to determine when to change the Pre-Packed Round. Color change is reversible and the exhausted material will change back to its original color. Saturated Pre-Packed Round should therefore be disposed of immediately after use.
- Excessive use of oxygen flush induces dehydration of the soda lime and inhibits its ability to absorb carbon dioxide.
- If the Pre-Packed Round is left mounted overnight, make sure that all gas flows are stopped. Otherwise the absorber may dry out and lose its activity. If the Pre-Packed Round dries out completely, it may give off CO (carbon monoxide) when exposed to anesthetic agents. Use carbon dioxide monitoring to determine when to change the Pre-Packed Round.
- For safety, always replace the absorber when its state of desiccation is unknown.
- The Pre-Packed Round must not be used with chloroform or trichloroethylene. The soda lime may react with chloroform/trichloroethylene slightly producing sodium formate/dichloroethylene, carbon monoxide and/or phosgene, which are harmful for the patient.
- Inhaled anesthetic agents in contact with strong bases are reported to be able to degrade to produce breakdown products of unknown toxicity in humans. For further information contact the manufacturer of the inhaled anesthetic agent in question.
- The Pre-Packed Round with soda lime is not an antimicrobiological filter or inhibitor. Based on this information, it is the responsibility of the clinician to take appropriate measures regarding patient cross contamination when required.
- Vital Signs™ recommends that the Vital Signs filter is used in order to minimize the cross contamination and the risk of foreign particles reaching the patient.
- It is recommended to use water traps in the patient circuit during low flow anesthesia (fresh gas flow less than 1 l/min) and when excessive condensed water accumulates in the limbs. Use water traps in both expiratory and inspiratory limbs to remove water condensation.
- Use only Vital Signs validated supplies and accessories approved by AirLife for proper performance of the anesthesia system. Using products other than those approved may result in improper performance of the anesthesia system.
- Keep all Pre-Packed Rounds that have been used with highly flammable anesthetics away from heat, sparks and open flames, as residual amounts of these materials will be present and may cause fire hazard.

Cautions

For information about the soda lime, refer to the appropriate Safety Data Sheet.

Local authority regulations must be investigated before disposal of the product.

Ordering information

For further information, please contact your local sales representative or visit myAirLife.com.

AirLife Finland Oy
Kuortaneenkatu 2
FI-00510 Helsinki, Finland
www.myAirLife.com

Made in UK

Australian Sponsor
Airlife Australia Holdings Pty Ltd.
PO BOX 97
NORTH RYDE BC NSW 1670
Australia

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