EN

Medisorb™ Loose Fill Twinpack 4.5 kg (9.9 lb.), Disposable, Carbon Dioxide Absorbent

Instructions for Use

REF 8570043 Medisorb™ Loose Fill Twinpack 4.5 kg (9.9 lb.), Disposable, Carbon Dioxide Absorbent

REF	Catalogue number		Consult instructions for use
LOT	Lot number		Follow instructions for use
QTY	Quantity		Keep dry
	Manufacturer		Fragile, handle with care
	Date of manufacture		This side up
	Expiry date		Keep away from sunlight
	Temperature limit		Danger

R_x ONLY U.S. Federal law restricts this device to sale by or on the order of a physician.

Intended use

The Medisorb™ Loose Fill Twinpack (later Medisorb carbon dioxide absorbent) is a carbon dioxide absorbent intended for use with anesthesia systems.

The Medisorb Carbon dioxide absorbent should only be used with air, oxygen, nitrous oxide, halothane, enflurane, isoflurane, desflurane and sevoflurane.

The product is intended to be used under constant attention of qualified professional healthcare personnel.

Technical description

Absorbent	Material
Granule size: 2.5–5 mm	Canister: Polyethylene
Moisture content: 12-19%	
12<pH<14	
Sodium hydroxide (NaOH, 215-185-5): <4%	
Calcium hydroxide (Ca(OH) ₂ , 215-137-3): >75%	

Instructions for use

Read these instructions carefully before using the product.

When to change the Medisorb carbon dioxide absorbent soda lime

The absorption of carbon dioxide is indicated by a gradual color change of the soda lime. Discard the product as soon as the color change has occurred. Intensity of the color change may vary from one procedure to another. The color change of the soda lime is only a rough indicator. Use carbon dioxide monitoring to determine when to change the Medisorb carbon dioxide absorbent. The color change is reversible and the exhausted material will change back to its original color. Saturated Medisorb carbon dioxide absorbent should therefore be disposed of immediately after use.

Operating room personnel must determine the remaining absorptive capacity of the Medisorb carbon dioxide absorbent immediately at the end of a surgical procedure. Excessive use of oxygen flush induces dehydration of the Medisorb carbon dioxide absorbent and inhibits its ability to absorb carbon dioxide.

There may be some amount of absorbent that does not change color in the bottom of the absorber; this is a natural occurrence.

It is recommended to use water traps in the patient circuit during low flow anesthesia (fresh gas flow less than 1l/min) and when excessive condensed water accumulates in the limbs. Use water traps in both expiratory and inspiratory limbs to remove water condensation.

Maintenance

Disposal

Saturated Medisorb carbon dioxide absorbent should be disposed immediately after use. Keep all Medisorb carbon dioxide absorbent that has been used with highly flammable anesthetics away from heat, sparks and open flames, as residual amounts of these materials will be present and may cause a fire hazard. Disposal of the used Medisorb carbon dioxide absorbent must be determined by the end user in compliance with federal state, provincial and local regulations. Specific disposal methods vary with the chemical (e.g., anesthetic agent, acid-base vapors, etc.) or biological agents used with the absorbent by the end user.

Storage

The unopened Medisorb carbon dioxide absorbent package must be closed and stored in a dry and clean environment at a temperature from 0°C/+32°F to +35°C/+95°F. Avoid direct sunlight (UV light). Avoid freezing which may reduce CO₂ absorption performance. Once a Medisorb carbon dioxide absorbent package is opened, its contents should be used immediately as the performance will deteriorate if storage is continued. The package shall be protected from physical damage and water. Under these conditions the Medisorb carbon dioxide absorbent in unopened Twin Pack containers will retain its carbon dioxide absorption capacity for five years from date of manufacture. The expiration date of the Medisorb carbon dioxide absorbent is printed on the product package.

Warnings

- Hazard and Precautionary statements of the soda lime: H314, Causes severe skin burns and eye damage. P260, Do not breathe dust/fume/gas/mist/vapours/spray. P280, Wear protective gloves/protective clothing/eye protection/face protection.
- P305/P351/P338, IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing.
- P303/361/353, IF ON SKIN (or hair): Remove/Take off immediately all contaminated clothing. Rinse skin with water/shower.
- P310, Immediately call a POISON CENTER or doctor/physician.
- Constant attention by a qualified professional is needed whenever this product is used.
- Avoid contact of the Medisorb carbon dioxide absorbent with eyes, skin or clothing. Avoid creating and breathing airborne dust that may cause irritation. In the event of skin or eye contact, immediately rinse the affected area with water and seek medical assistance.
- Use protective gloves and protective glasses when filling containers with Medisorb carbon dioxide absorbent. Do not use the very last fraction of soda lime from the Medisorb carbon dioxide absorbent package, since it contains more dust.
- The color change of the soda lime is only a rough indicator. Use carbon dioxide monitoring to determine when to change the Medisorb carbon dioxide absorbent. The color change is reversible and the exhausted material will change back to its original color. Saturated Medisorb carbon dioxide absorbent should therefore be disposed immediately after use.
- Excessive use of oxygen flush induces dehydration of Medisorb carbon dioxide absorbent and inhibits its ability to absorb carbon dioxide.
- If an absorber with Medisorb carbon dioxide absorbent is left mounted overnight, make sure that all gas flows are stopped. Otherwise the absorber may dry out and lose its activity. If the Medisorb carbon dioxide absorbent dries out completely, it may give off CO (carbon monoxide) when exposed to anesthetic agents. Use carbon dioxide monitoring to determine when to change the Medisorb carbon dioxide absorbent.
- For safety, always replace the absorber when its state of desiccation is unknown.
- The Medisorb carbon dioxide absorbent must not be used with chloroform or trichloroethylene. Medisorb carbon dioxide absorbent may react with chloroform/trichloroethylene slightly producing sodium formate/dichloroethylene, carbon monoxide and/or phosgene, which are harmful for the patient.
- Inhaled anesthetic agents in contact with strong bases are reported to be able to degrade to produce breakdown products of unknown toxicity in humans. For further information contact the manufacturer of the inhaled anesthetic agent in question.
- Medisorb carbon dioxide absorbent is not an antimicrobiological filter or inhibitor. Based on this information, it is the responsibility of the clinician to take appropriate measures regarding patient cross contamination when required.
- Vital Signs™ recommends that the Vital Signs filter is used in order to minimize the cross contamination and the risk of foreign particles reaching the patient.
- It is recommended to use water traps in the patient circuit during low flow anesthesia (fresh gas flow less than 1l/min) and when excessive condensed water accumulates in the limbs. Use water traps in both expiratory and inspiratory limbs to remove water condensation.
- Use only Vital Signs validated supplies and accessories approved by AirLife for proper performance of the anesthesia system. Using products other than those approved may result in improper performance of the anesthesia system.
- Keep all Medisorb carbon dioxide absorbent that has been used with highly flammable anesthetics away from heat, sparks and open flames, as residual amounts of these materials will be present and may cause fire hazard.

Cautions

For information about the soda lime, refer to the appropriate Material Safety Data Sheet.

Local authority regulations must be investigated before disposal of the product.

Ordering information

For further information, please contact your local sales representative or visit myAirLife.com.

AirLife Finland Oy
Kuortaneenkatu 2
FI-00510 Helsinki, Finland
www.myAirLife.com

0537

Made in UK

Australian Sponsor
Airlife Australia Holdings Pty Ltd.
PO BOX 97
NORTH RYDE BC NSW 1670
Australia

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