EN

English

Salter Labs™ Medium Concentration Oxygen Mask

DEVICE DESCRIPTION:
The device is non-sterile, disposable and for single-patient multi-use. An oxygen mask is a device placed over a patient's nose or mouth to administer oxygen.

INTENDED PURPOSE:
A Medium Concentration Oxygen Mask is a device placed over a patient's nose or mouth to administer oxygen.

INDICATIONS FOR USE:
For patients requiring supplemental oxygen. To treat hypoxia, relieve shortness of breath and pre-oxygenate for surgical procedures in spontaneously breathing patients.

ENVIRONMENT:
Hospitals, sub-acute, medical clinics, pre-hospital, home, surgical centers, skilled nursing facilities.

INTENDED USER:
Trained medical professionals.

PATIENT TARGET GROUP:
Infant, pediatric, adult.

EXPECTED CLINICAL BENEFITS:
The oxygen mask delivers up to 15 LPM of oxygen to reverse or prevent hypoxemia, and provide relief of shortness of breath in oral and nasal breathers.

CONTRAINDICATIONS:
None Known.

- ⚠️ WARNINGS:**
- If oxygen is in use, do not use near flame or heat source.
 - Position tubing to avoid strangulation.
 - Device contains components that may present choking hazard.
 - Patient may become hypoxic if oxygen flow is interrupted.
 - Use a flow rate ≥ 5 LPM to prevent CO2 accumulation and re-breathing.
 - Keep excess tubing loosely coiled and out of the way to prevent kinking and tripping hazard.
 - To reduce the risk of misconnections and patient injury, always trace tubing from gas source to the medical device before connecting.

- ⚠️ CAUTIONS:**
- Consult physician if patient develops infection, skin irritation or material sensitization.
 - Do not place anything on supply tubing that may obstruct flow.
 - This is a single-patient multi-use product, reuse on more than one patient may cause a risk of cross-contamination.

RESIDUAL RISKS:
Refer to Warnings and Cautions.

NECESSARY USER QUALIFICATIONS:
The end-user should receive training on the proper use of the device and be familiar with the indications for use, warnings, and cautions as stated in the IFU.

Rx Only: Federal law (USA) restricts this device to sale by or on the order of a licensed healthcare practitioner.

INSTRUCTIONS:
Follow instructions to avoid risk of infection and contamination.

- 1) Connect tubing to oxygen source.
- 2) Place mask over nose and mouth, adjust elastic strap and nose bar for a snug fit.
- 3) Check for leaks at connections and for proper placement prior to use.
- 4) Add humidification if nasal dryness or irritation occurs.

This product is disposable and is not intended for prolonged use greater than 14 days.

CLEANING/REUSE INSTRUCTIONS:
A discontinued use of this product on the same patient, with periods of use separated by a period of non-use, is acceptable. The user must ensure that the product is not damaged or contaminated between usages on the same patient.

SAFE DISPOSAL:
Decontaminate and dispose of all potentially biohazardous material. Dispose of device in accordance with local, state or national regulations.

INCIDENT REPORTING:
Notice to the user that any serious incident that has occurred in relation to the device should be reported to the manufacturer at productquality@myairlife.com, and the competent authority of the Member State in which the user and/or patient is established. If the patient's physical symptoms worsen or experience a sudden change in their condition (e.g., increased shortness of breath, fever, dizziness) call the physician.

PERFORMANCE CHARACTERISTICS:

- For use with Flow Rate 3 - 15 LPM
- For use with Bonded O2 Tubing adaptor

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