








Material Specification:7DayCatheter.INIFU

Material Specification	7DayCatheter.INIFU	Style	Flat with 4 x 1/4” rounded corners (RC), machine trim; Aqueous coating (STARKOTE AQ-705) on front side of Insert
Revision	1	Size (mm) (OD) (W)	N/A
Description	IFU - BALLARD 7DAY Closed Suction System	Size (mm) (OD) (D)	N/A
Item Number	7DayCatheter.INIFU	Size (Inches) (OD) (D)	N/A
Previous item number	N/A	Size (mm) (OD) (L)	N/A
SAP Number	7130020	Size (Inches) (OD) (W)	N/A
Published Date	DCR-17047	Size (Inches) (L)	4.875
Item Description	IFU - BALLARD 7DAY Closed Suction System	Other Requirements	N/A
Material Type	IFU / MANUAL / BOOKLET	Primary Vendor	Litho Flexo Graphics Inc / Salt Lake City, UT
OS Size (Inches) (OD) (W)		Alternate Vendor	N/A
OS Size (Inches) (OD) (D)	N/A	Size (Inches) (ID) (W)	N/A
OS Size (Inches) (OD) (L)	N/A	Size (Inches) (ID) (D)	N/A
Printing Details	N/A	Flute	N/A
Roll Core Size	N/A	Size (Inches) (ID) (L)	N/A
Size (mm) (L)	123.825	Size (mm) (W)	323.85
Target Min. Compression	N/A	Size (mm) (D)	N/A
Size (Inches) (W)	12.75	Size (mm) (ID) (D)	N/A
Manufacturer's Joint	N/A	Size (Inches) (OD) (L)	N/A
Size (Inches) (D)	N/A	Size (mm) (ID) (W)	N/A
ECT	N/A	Master Text	(EN) English - United States
Material Description	100 # Bookgloss White Paper	Translation Language(s)	N/A
Printing Colors	Process Black	Size (mm) (ID) (L)	N/A
Roll Diameter	N/A	PDF Viewable Files	DayCatheter.INIFU, IFU - BALLARD 7DAY Closed Suction System.pdf
Bar Code	N/A	Latest Revision	Initial release per DCR-17047
Mill Part Number	N/A	UPC	N/A
Label Shell Available	N/A	Percentage (%) of Post-Consumer Recycled Content	N/A
Dieline Available	N/A	Material Contains the Following Synthetics	N/A
Primary Vendor Location	Salt Lake City, UT	Percentage (%) of Pre-Consumer Recycled Content	N/A
Primary Vendor Contact	Paulette Camden	Is Material Compostable?	N/A
Dieline	N/A	Is Material Biodegradable?	N/A
Adhesive	N/A	Is Material Laminated?	N/A
Board Grade	N/A	Is Material Coated with Plastic?	N/A
Alternate Vendor Contact	N/A	Is Material Labelled with Recycling Information?	N/A
Primary Vendor Name	Litho Flexo Graphics Inc	Is Material Recyclable?	N/A
Alternate Vendor Location	N/A	Market Clearance Number	N/A
Primary Vendor Phone	801-994-1130	Weight (g)	N/A
Alternate Vendor Phone	N/A	Thickness Unit	N/A
Revision History	New per DCR-17047	Thickness (Pts)	N/A
Drawing Number	N/A	Special Considerations	N/A
UKCA	No	Master Specification	N/A
EUR24	No		

BALLARD™ 7DAY
CLOSED SUCTION SYSTEM



Interval
volume

Diameter

Length

Single
use only

STERILE R

Do not use if
package is damaged

Do not
resterilize

Not made with
natural rubber latex

Product is NOT
made with DEHP

R_{only}

Caution

MR MRI Safe

Endotracheal
Length

Tracheostomy
Length

50°C

-20°C

7DayCatheter.INIFU Rev 1 2025/08

SunMed, LLC,
a Subsidiary of SunMed Holding LLC,
2710 Northridge Dr. NW, Suite A
Grand Rapids, MI 49544 USA
www.myAirLife.com
Made in Mexico

These instructions apply to the following families of
BALLARD products:

- Double Swivel Elbow
 - T-Piece
- These instructions also apply to the following
configurations of the above families of products:
- Endotracheal Length
 - Tracheostomy Length

Device Description:
A Closed Suction Catheter is a product that enables
maintaining the patency of a patient's artificial airway
without breaking the ventilation circuit during suctioning.
It is a protected suction catheter inside a sterile plastic
sleeve.

Intended Use:
The BALLARD Closed Suction Catheter is a product that
enables maintaining the patency of a patient's artificial
airway without breaking the ventilation circuit during
suctioning.

Indications for Use:
The BALLARD Closed Suction Catheter is intended to
suction a patient's artificial airway without breaking the
ventilation circuit to allow for a patient's airway. The target
population includes adult patients with artificial airways.

Environment:
• Hospital
• Home

MR Safety:
• MRI Safe

Patient Target Group:
• Adult

Intended User:
The BALLARD 7Day Closed Suction Catheter devices
are intended to be used by trained medical professionals,
or trained caregivers.

Contraindications:
None known.

- Warnings:**
1. Cap on BALLARD T-Piece prevents continuous flow
therapy. Remove cap before starting continuous flow
therapy. Failure to remove cap prior to continuous
flow therapy may result in serious injury or death.
 2. Do not trim or cut the endotracheal tube (not
supplied) while the BALLARD catheter is attached,
otherwise the BALLARD catheter may also be cut
and that portion of the catheter may be aspirated into

the lower respiratory tract of the patient and may
cause death or serious injury.

3. Do not reuse, reprocess, or resterilize this medical
device. Reuse, reprocessing, or resterilization may
1) adversely affect the known biocompatibility
characteristics of the device, 2) compromise the
structural integrity of the device, 3) lead to the device
not performing as intended, or 4) create a risk of
contamination and cause the transmission of
infectious diseases resulting in patient injury, illness,
or death.

Cautions:

1. Inspect the BALLARD 7Day catheter package before
opening. Do not use product if packaging has been
compromised. Non-sterile contents may cause
infection.
2. Excess fluid in heat and moisture exchanger (HME)
may increase gas flow resistance. When introducing
fluid into T-Piece, ensure that fluid does not enter
HME.
3. Always place the thumb valve in the locked position
when not in use to prevent inadvertent activation.
4. Single patient use only.
5. The BALLARD 7Day catheter is intended to be used
for 7 days before changing. Change more frequently
if catheter becomes heavily soiled during use.
6. Remove red wedge plug from BALLARD 7Day
catheter prior to use.
7. Inspect sodium chloride vial prior to opening. Do not
use product if vial has been compromised.
Compromised contents may cause infection.
8. Do not use 54 cm (21.3 inch) catheters on
tracheostomy patients. Mucosal damage may result.
9. Select the appropriate size BALLARD 7Day catheter.
Most experts suggest that the catheter selected
should occupy no more than one half of the internal
diameter of the artificial airway.
10. Do not leave the catheter within the airway. Always
pull back until the black stripe is visible within the
sleeve. Any catheter left extended into the airway will
cause increased airway resistance.
11. Use appropriate regulated vacuum levels. Most
experts suggest -80 to -120 mmHg (-10.7 to -15.9
kPa).

12. Use appropriate suction technique. Most experts
suggest that the entire suction procedure should last
no longer than 10 to 15 seconds and that actual
duration of negative pressure should be no longer
than 5 to 8 seconds per episode.
13. Always use caution and good clinical judgement no
matter what ventilator mode is in use. If the clinician
notes any signs of suction intolerance such as
oxygen desaturation, negative ventilator system
pressures, patient stress or excessive discomfort,
adjustments to the ventilator settings may need to be
made. These adjustments (please refer to the
ventilator's instructions for use) may include
manipulation of the inspiratory trigger sensitivity,
inspiratory volume or flowrate, and selection of a
different ventilator mode; or may require the use of
an alternate suction technique. Failure to follow the
above precautions may increase the risk of positive
and negative barotrauma.

Necessary User Qualifications:
Federal law (USA) restricts this device to sale by or on
the order of a licensed healthcare practitioner.

Instructions for Use:

- Setup:**
1. Select appropriate size BALLARD 7Day catheter.
 2. Attach thumb control valve to suction tubing.
 3. Depress and hold thumb valve and simultaneously
adjust vacuum regulator to desired level.
 4. Release thumb control valve and attach BALLARD
7Day catheter between patient and the ventilator
circuit.

- Suggested Suction Procedure:**
1. Stabilize BALLARD 7Day catheter and endotracheal
(ET) adapter with one hand then push the catheter
into the endotracheal tube with the thumb and
forefinger of the opposite hand (Fig. 1).
 2. Advance BALLARD 7Day catheter to desired depth.
 3. Depress and hold thumb control valve, then gently
withdraw catheter. Stop withdrawal when black
marking ring is visible inside sleeve (Fig. 2).
 4. Release thumb control valve.
 5. Repeat steps 1 through 4 above as necessary.

Catheter Irrigation Instructions:

1. Be sure the black marking ring is visible in the sleeve

- (Fig. 2). Open the cap on the irrigation port.
2. Introduce fluid slowly into the port, simultaneously
depress the thumb control valve (Fig. 3).
3. Continue to irrigate until the BALLARD 7Day
catheter is clear (Fig. 4).
4. Close the cap on the port.
5. Lift and turn the thumb control valve 180 degrees to
lock in position (Fig. 5).
6. Place the BALLARD 7Day catheter and suction
tubing alongside the breathing circuit.

**Control Depth Suction: (Endotracheal length
catheters printed with numerical markings only):**

1. Align any printed depth number on the catheter with
the similar number printed on the endotracheal tube,
or
2. Observe the printed number on the endotracheal
tube closest to the endotracheal tube adapter. Add 8
cm to this number. Advance catheter until the sum
(depth plus 8) appears in the window directly across
from irrigation port connector.
3. Catheter tip will be within 1 cm of the end of the
artificial airway.

Thumb Control Valve Operation:

1. The thumb control valve can be locked to prevent
inadvertent or accidental suction. To lock, lift white
part of thumb control valve and rotate 180 degrees.
To unlock, repeat this action (Fig. 5).

Day Sticker Usage:

1. Apply the appropriate day sticker to the thumb
control valve.

**These instructions also apply to the following
configurations of the above families of products:**

• Tracheostomy Length

Tracheostomy Patients:

1. Use tracheostomy 30 cm (12 inch) BALLARD 7Day
catheter for patients with tracheostomy artificial
airway only, per hospital policy. If 30 cm BALLARD
7Day catheter is used on endotracheal artificial
airway, ineffective suction may result.
2. Do not use 54 cm (21.3 inch) catheter on
tracheostomy patients. Mucosal damage may result.

4.875"

Aug 19, 2025

Preprinted Insert-Instructions for use

12.750" (W) x 4.875" (L)

7DayCatheter.INIFU

Prints 1/c Process Black