SKIN TEMPERATURE SENSOR

	MR UNSAFE
	MEDICAL DEVICE
	STERILIZED WITH ETHYLENE OXIDE
	SINGLE STERILE BARRIER SYSTEM
	DO NOT REUSE
	DO NOT RESTERILIZE
	DO NOT USE IF PACKAGING IS DAMAGED
	NOT MADE WITH NATURAL RUBBER LATEX
	FEDERAL U.S.A. LAW RESTRICTS THIS DEVICE TO SALE OR USE BY OR ON THE ORDER OF A PHYSICIAN OR PROPERLY LICENSED PRACTITIONER.

IMPORTANT
 Please read all warnings and instructions before use. Correct application is essential for proper product function and to reduce the risk of injury.

INTENDED USE
The AirLife Skin Temperature Sensor is to be used for routine monitoring of a patient's skin surface temperature.

INTENDED USERS
The AirLife Skin Temperature Sensor is intended to be used only by a licensed practitioner or trained healthcare professionals.

PATIENT TARGET GROUP(S)
The intended patient population is pediatric and adults.

PRODUCT DESCRIPTION
The AirLife Skin Temperature Sensor consists of a temperature-sensing wire set with an adhesive sensor cover that facilitates application to the patient's skin surface. The temperature monitoring sensor located at the distal end of the wire is affixed to the center of the disk-shaped foam pad, which has a reflective backing and an adhesive layer to hold the sensor in place.

The temperature monitoring sensor is located at the tip of the sensor. The blue connector at the proximal end interfaces with manufacturer-branded cables for connection to a monitoring unit to display temperature. The temperature monitoring sensor is available in YSI 400 or YSI 700 series configurations for connection to a compatible monitoring unit. Whether the sensor is a 400 or 700 series configuration is indicated on the device labeling. Consult the user manual for the patient monitor to ensure compatibility. The sensors are designed to interface with manufacturer-branded cables for connection with YSI 400 or 700 series compatible monitors, including the following patient monitors and equivalent models: Mindray® Passport, Philips® IntelliVue®, Siemens/Dräger® Infinity®, and GE® Datex-Ohmeda® brands.

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INDICATIONS
The AirLife Skin Temperature Sensor is indicated to monitor the patient's skin surface temperature.

CONTRAINDICATIONS
None known.

ADVERSE REACTIONS
While adverse reactions are rare, the following have been associated with use of the skin temperature sensor during use or immediately after use:

- electrical burns,
- skin irritation,
- inflammation, or
- allergic reaction to materials.

WARNINGS
 The AirLife Skin Temperature Sensor should only be used on healthy skin.

- **DO NOT** apply the adhesive cover on this sensor to a skin surface which has been damaged by abrasion, laceration, incision, burning, irritation, chafing or inflammation.
- **DO NOT** use this skin sensor on neonatal patients.
- **DO NOT** use sensor without sensor cover. Use without sensor cover can cause inaccurate readings, and could potentially lead to patient injury.
- Carefully inspect the packaging and device for damage or defect. **DO NOT** use if the packaging is open or the device is damaged.

CAUTIONS

- Device is intended for use by a licensed practitioner.
- Take care to ensure the cable and connector do not get wet and proper techniques are used during electrosurgical procedures to reduce radio frequency interference, current and potential burns to the patient.
- An adequate electrosurgical dispersive ground electrode close to the active surgical site should be properly connected.
- **DO NOT** intertwine the cables, especially the monitor cables, with the electrosurgical unit's cables.
- The operation of the patient temperature monitor may be temporarily affected during electrosurgical activations. Unusual temperature readings should be checked.
- If temperature does not read properly, check all connections. If problem persists, discard and use a new sensor.
- **DO NOT** use an expired device.
- In case of adverse skin reaction, discontinue use.
- The skin sensor is intended to be used with manufacturer-branded interface cables and sensor cover. Use with an incompatible cable may affect performance.

SAFETY IN MRI NOT EVALUATED
The AirLife Skin Temperature Sensor's safety in MRI has not been evaluated, hence, the sensor is considered MR Unsafe.

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DIRECTIONS FOR USE

1. Prepare the skin for adhesive pad by removing body oils and hair.
2. Verify the compatibility of the sensor, interface cable, and patient monitor before use.
3. Remove the AirLife Skin Temperature Sensor from its package.
4. Peel off release liner and apply sensor with sensor cover to patient's skin.
5. Connect the temperature sensor to the appropriate reusable interface cable from AirLife. Connect the cable to the temperature monitor. Secure the interface cable using the drape clip.
6. Follow the directions for use of the patient temperature monitor instrument.
7. When temperature monitoring is completed, carefully remove adhesive cover from patient's skin.
8. Once temperature monitoring is no longer needed, disconnect the sensor (at the connector) and discard the sensor. To disconnect, grasp both connectors firmly and pull. **DO NOT** pull on cable or wire.

OPERATING SPECIFICATIONS

- Operating Temperature: 25° C to 45° C
- Rated Output Range: 35° C to 42° C
- Rated Extended Output Range: 25° C to 35° C and 42° C to 45° C
- The reference body site is the skin surface temperature
- The measuring site(s) are intact skin
- Accuracy in Rated Output Range: +/- 0.2° C
- Accuracy in Rated Extended Output Range: +/- 0.3° C
- Time response:
 - The temperature sensor has a direct mode of operation
 - The maximum heating transient time is 12 seconds
 - The maximum cooling transient time is 10 seconds

DISPOSAL
Dispose of any used device according to local, state, and federal laws and regulations. For safe disposal of devices, follow your facility's protocol.

	KEEP DRY
	KEEP AWAY FROM SUNLIGHT

In addition to the competent authority in the country where the patient resides, serious incidents must be reported to DeRoyal Industries, Inc.

WARRANTY
DeRoyal® products are warranted for one hundred twenty (120) days from the date of shipment from DeRoyal as to product quality and workmanship. **DEROYALS WRITTEN WARRANTIES ARE GIVEN IN LIEU OF ANY IMPLIED WARRANTIES, INCLUDING WARRANTIES OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE.**

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DeRoyal uses the symbols and meanings from standards ISO 15523-1. These symbols are placed next to the text explaining their meaning in this instructions for use (IFU), a complete symbols glossary is available online at deroyal.com/symbols or by contacting customer service.

DISTRIBUTED BY:
AirLife
2710 Northridge Dr. NW, Suite A
Grand Rapids, MI 49544 USA
www.myAirLife.com

MANUFACTURED BY:
DeRoyal Industries, Inc.
200 DeBusk Lane
Powell, TN 37849 USA
888.938.7828 or (001) 865.938.7828
www.deroyal.com

CE 2460

MedEnvoy Global B.V.
Prinses Margrietplantsoen 33 – Suite 123
2595 AM the Hague
The Netherlands

UK Responsible Person
MedEnvoy UK Limited
85, Great Portland Street, First Floor
London, W1W 7LT
United Kingdom

EMERGO AUSTRALIA
Level 20 Tower II
Darling Park
201 Sussex Street
Sydney, NSW 2000
Australia

PART#74-14902G-VY | REVISED 2/25

©2025 DeRoyal Industries, Inc. All Rights Reserved.
DeRoyal and the DeRoyal logo are registered trademarks of DeRoyal Industries, Inc.
AirLife is a trademark of SunMed Group Holdings, LLC.
Mindray is a registered trademark of Shenzhen Mindray, Bio-Medical Electronics Co. Ltd.
Siemens, Draeger, and Infinity are registered trademarks of Siemens Trademark GmbH Co.
Philips and IntelliVue are registered trademarks of Koninklijke Philips N.V.
GE and Datex-Ohmeda are trademarks or registered trademarks of General Electric Company.

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