
DeRoyal uses the symbols and meanings from standards ISO 15823-1. These symbols are placed next to the text explaining their meaning in this instructions for use (IFU), a complete symbols glossary is available online at deroyal.com/symbols or by contacting customer service.

AirLife®

DISTRIBUTED BY:
AirLife
2710 Northridge Dr. NW, Suite A
Grand Rapids, MI 49544 USA
www.myAirLife.com

MANUFACTURED BY:
DeRoyal Industries, Inc.
200 DeBusk Lane
Powell, TN 37849 USA
888.938.7828 or (001) 865.938.7828
www.deroyal.com

MedEnvoy Global B.V.
Prinses Margrietplantsoen 33 – Suite 123
2595 AM the Hague
The Netherlands

UK Responsible Person
MedEnvoy UK Limited
85, Great Portland Street, First Floor
London, W1W 7LT
United Kingdom

EMERGO AUSTRALIA
Level 20 Tower II
Darling Park
201 Sussex Street
Sydney, NSW 2000
Australia

PART#74-14904G-VY 1 REVISED 2/2025

©2025 DeRoyal Industries, Inc. All Rights Reserved.
DeRoyal and the DeRoyal logo are registered trademarks of DeRoyal Industries, Inc.
AirLife is a trademark of AirLifeSunMed Group Holdings, LLC.

TYMPANIC TEMPERATURE PROBE

	MR UNSAFE
	MEDICAL DEVICE
	NON-STERILE
	DO NOT REUSE
	DO NOT USE IF PACKAGE IS DAMAGED
	NOT MADE WITH NATURAL RUBBER LATEX
RX ONLY	FEDERAL U.S. LAW RESTRICTS THIS DEVICE TO SALE OR USE BY OR ON THE ORDER OF A PHYSICIAN OR PROPERLY LICENSED PRACTITIONER.

IMPORTANT
Please read all warnings and instructions before use. Correct application is essential for proper product function and to reduce the risk of injury.

INTENDED USE
The AirLife Tympanic Temperature Probe is to be used for routine monitoring of the core body temperature.

INTENDED USERS
The AirLife Tympanic Temperature Probe is intended to be used only by a licensed practitioner or trained healthcare professionals.

PATIENT TARGET GROUP(S)
The intended patient population is for pediatric and adult.

PRODUCT DESCRIPTION
The AirLife Tympanic Temperature Probe consists of a temperature-sensing wire set to which an ear piece is affixed at the distal end for application to the aural canal. The ear piece may be comprised of a foam ear plug and cotton ball or just a cotton ball.

The connector at the proximal end interfaces with manufacturer-approved cables for connection to a monitoring unit to display temperature. The probes are designed to interface with manufacturer-branded cables for connection with YSI 400 series compatible monitors. Consult the user manual for the patient monitor to ensure compatibility.

INDICATIONS
The AirLife Tympanic Temperature Probe is indicated to be used for the monitoring of core body temperature via the aural canal.

CONTRAINDICATIONS
The use of this device is contraindicated for ears in which pathological conditions exist, such as infections, inflammations, tumors, perforated eardrums, or polyps. It is contraindicated in neonates and in patients whose ear is too small or with anatomical shape and size not suitable for insertion of the temperature probe.

WARNINGS
• This is a single use device. It is not intended for disinfection and/or subsequent re-use. Re-use could result in microbial contamination and potential patient injury.
• This device has not been evaluated for reprocessing or sterilization. Reprocessing and/or sterilization may damage the device, rendering

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it unusable and/or may lead to device failure.
• Never use force to insert the probe into the ear to overcome resistance.
• **DO NOT** insert probe to the point where the probe tip touches the tympanic membrane.
• Carefully inspect the packaging and device for damage or defect.
• **DO NOT** use if the packaging is open or the device is damaged.

ADVERSE REACTIONS
While adverse reactions are rare, the following are potential adverse reactions associated with the use of temperature probes during insertion or while the device is in use:
• bleeding from the tympanic membrane or ear canal
• ear ache
• laceration of the ear canal
• otitis externa
• perforation of the tympanic membrane
• decreased hearing
• allergic reaction to materials

CAUTIONS
• Perform an otoscopic examination to verify that the ear canal is free of obstruction and pathology.
• In an awake patient, the probe should be partially or completely withdrawn when discomfort is reported.
• The patient's head should not lie on the side where the probe is inserted.
• Take care to ensure that the cable and connector do not get wet and proper techniques are used during electro-surgical procedures to reduce radio frequency interference current and potential burn to the patient.
• An adequate electro-surgical dispersive ground electrode close to the active surgical site should be properly connected.
• **DO NOT** intertwine the cables, especially the monitor cables with the electro-surgical unit's cables.
• The operation of the patient temperature monitor may be temporarily affected during electro-surgical activation. Unusual temperature readings should be checked.
• The AirLife Tympanic Temperature Probe is intended to be used with manufacturer-approved interface cables. Use with an incompatible cable may affect performance.

SAFETY IN MRI NOT EVALUATED
The AirLife Tympanic Temperature Probe's safety in MRI has not been evaluated, hence, the probe is considered MR Unsafe.

DIRECTIONS FOR USE
1. Perform an otoscopic examination of the ear where temperature is to be monitored. Remove any excess cerumen (ear wax).
2. Verify the compatibility of the probe, interface cable, and patient monitor before use.
3. Remove the AirLife Tympanic Temperature Probe from its package.
4. Gently pull on the pinna (outer ear) to "straighten" the ear canal.
5. Gently insert the probe (cotton end) into the ear canal avoiding contact with the tympanic membrane. The ear plug's flexible structure serves as an anchor and an insulator for the temperature sensor.
6. **DO NOT OVER-INSERT.** Probe tip does not have to touch the tympanic membrane for this temperature probe to work properly. Stop to withdraw or reposition when resistance is encountered.
7. Connect the temperature probe to the appropriate reusable cable. Connect the cable to the patient monitor. Secure the cable using the drape clip.
8. The sensor lead wire may be taped to the patient's head for

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Before using this device, visit our website for the latest revision.

convenience.
9. Follow the directions for use of the patient temperature monitor instrument.
10. When temperature monitoring is completed, carefully remove the probe.
11. Once temperature monitoring is no longer needed, disconnect the probe (at the connector) and discard the probe. To disconnect, grasp both connectors firmly and pull. **DO NOT** pull on cable or wire.
12. Perform a second otoscopic examination of the ear canal and tympanic membrane according to hospital protocol.

OPERATING SPECIFICATIONS
• Operating Temperature: 25° C to 45° C
• Rated Output Range: 35° C to 42° C
• Rated Extended Output Range: 25° C to 35° C and 42° C to 45° C
• The reference body site is the core body temperature
• The measuring site(s) is the tympanic membrane
• Accuracy in Rated Output Range: +/- 0.2° C
• Accuracy in Rated Extended Output Range: +/- 0.3° C
• Time response:
• The temperature probe has a direct mode of operation
• The maximum heating transient time is 15 seconds
• The maximum cooling transient time is 10 seconds

DISPOSAL
Dispose of any used device according to local, state, and federal laws and regulations. For safe disposal of devices, follow your facility's protocol.

	KEEP DRY
	KEEP AWAY FROM SUNLIGHT

In addition to the competent authority in the country where the patient resides, serious incidents must be reported to DeRoyal Industries, Inc.

WARRANTY
DeRoyal products are warranted for one hundred twenty (120) days from the date of shipment from DeRoyal as to product quality and workmanship. **DEROTAL'S WRITTEN WARRANTIES ARE GIVEN IN LIEU OF ANY IMPLIED WARRANTIES, INCLUDING WARRANTIES OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE.**

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