Device Description: MYOTouch is a hand-held reusable device designed to be used for surgical correction of neonatal anorectal malformations. MYOTouch is intended for prescription use. The device is a 9V battery-powered, one-channel muscle stimulator. The device is supplied with a Sterile Stylus, which connects to the control unit by cable and plug.

Intended Use: Used to identify striated muscle in the treatment of high and low anorectal malformations.

Indications for Use: The MYOTouch Muscle Stimulator is indicated for use as an adjunct in the treatment of high and low anorectal malformations by helping with the identification of striated muscles to be used in anal reconstructions. The MYOTouch device is intended for use by trained personnel who are competent to apply the appropriate stimuli. The surgeon uses a bipolar stimulating probe to electrically stimulate the muscle which will twitch in response to the electrical stimulus provided by the MYOTouch electrical stimulator. The MYOTouch device is used for surgical correction of neonatal anorectal malformations. It is up to the pediatric surgeon to identify the muscle tissue used for anorectal reconstruction. The MYOTouch device is a Muscle Stimulator, and the device is not intended for use in a life support situation. Improper operation of the MYOTouch device can result in personal injury to a patient or to a device operator.

WARNINGS AND CAUTIONS

- CAUTION: Federal Law (USA) restricts this device to sale or on order by a physician or appropriate licensed practitioner.
- CAUTION should be used for patients with suspected or diagnosed heart problems.
- CAUTION should be used for patients with suspected or diagnosed epilepsy.
- CAUTION: Avoid using the MYOTouch device in close proximity to RFID readers or avoid using it concurrently with an RFID reader, as there are emission frequencies that may adversely affect the device.
- CAUTION: Strong electromagnetic fields (electro-surgery/ microwave cookers/ mobile phones) may affect the correct operation of this unit. If it appears to behave unusually, move it away from these devices.
- CAUTION: DO NOT use the device if a battery has leaked acid into the battery compartment. The leaked battery acid could impair the internal circuits for the device. Please contact AirLife to inquire about a device replacement.
- DO NOT: Use the device for nerve localization for anesthetic regional block.
- DO NOT: Use the device on patients with neuromuscular or skin diseases.
- DO NOT: Use the device in the presence of a flammable anesthetic mixture with air, oxygen, or nitrous oxide.
- DO NOT: Modify any components of this equipment.
- DO NOT: Use applied across or through the head, directly on the eyes, covering the mouth, on the front of the neck, (especially the carotid sinus), or from electrodes placed on the chest and the upper back or crossing over the heart.
- DO NOT: Use transcranially. Doing so may influence brain function.
- DO NOT: Stimulate over, or in proximity to, cancerous lesions.
- The MYOTouch device may be hazardous to patients with implanted electrical medical devices.
- Stimulus current must be increased gradually until the supramaximal stimulus is achieved. Applying currents greater than necessary for supramaximal stimulation may increase the risk of skin burns.
- Prior to using the MYOTouch device the patient's skin should be cleaned and completely dried. The affected area for use should be free of excessive hair, scar tissue, or any other lesions.
- The MYOTouch controller should not contact liquids.
- DO NOT: Attach or remove the probe from the controller until the power has been turned OFF.
- The MYOTouch controller must be used with 9V alkaline battery only.
- The operator of the MYOTouch device cannot touch the actual battery and the patient simultaneously.
- The battery should be removed from the MYOTouch controller for long term storage.
- Prior to each use, check the MYOTouch device for proper condition and functionality.
- WARNING: There may be residual muscle paralytics in the patient which may block stimulation of skeletal muscles with the MYOTouch device. Anesthesia check for the presence of paralytics is required prior to the use of the device. Misleading results from device use may result in unnecessary additional surgical dissection causing patient injury.
- WARNING: Simultaneous connection to high frequency surgical equipment may result in burns and damage to the stimulator.
- WARNING: The stylus tips will discharge an electrical shock to any conductive material that they contact.
- WARNING: To minimize the risk of thermal burns due to tissue heating, do not keep the probe in one location for any length of time (for example, level 18 for more than 20 seconds and level 20 for more than 1 second could result in tissue burn).
- WARNING: The bipolar probe of device has current densities exceeding 2 mA/cm2 may require special attention of the operator.

Intended Patient Population

The intended population for the MYOTouch device is as listed below:

- Neonate/Newborn (birth through 28 days)
- Infant (from 29 days to 2 years of age)
- Child (from 2 years to 12 years of age)
- Adolescent (from 12 years to 18 years of age)

Intended User

The intended user is a specialist clinician trained in neurophysiological studies or a pediatric surgeon trained in the specialist reconstruction surgeries for which the MYOTouch is indicated for.

Environment of Use

Clinics, hospitals. The MYOTouch Muscle Stimulator is intended for static tabletop use in a clinical environment, typically a neurosurgery operating theatre.

Contraindications

- Do not use on a patient with a cardiac pacemaker, implanted defibrillator, or other implanted metallic or electronic device. Such use could cause electric shock, burns, electrical interference, or death.

Adverse Reactions

Skin irritation and burns beneath the electrodes have been reported with the use of powered muscle stimulators.

Clinical Benefit

Aids in the identification of striated muscles to be used in anal reconstructions.

INSTALLATION AND INITIAL CHECK

Before installation and use of the MYOTouch device, ensure that this manual has been read and understood.

Before use with a patient, the MYOTouch device should be tested in its operational location using a simulated situation.

Caution: To assure patient safety, on-site training is required.

Check if the 9V battery is inserted into the battery compartment correctly.

Without any button operation, ensure that the plug of stimulator stylus MS 200 leadwire is inserted into the socket of MYOTouch device correctly.

Press the ON/OFF button; check if the BATTERY indicator is illuminated.

Operational Tests

These tests are to be performed only by qualified personnel.


- Switch on the device by pressing the ON/OFF button. Increase the intensity setting to any level higher than 0 by pressing the “+” button on side.
- Open circuit symbol checking: The stimulator stylus tips do not contact any item. Press “SINGLE PULSE” or “CONTINUOUS PULSE ON” button, to check if the PULSE indicator is illuminated in yellow. Also, the triangle warning symbol “add symbol” is flashing on the LCD display.

Output function checking

The following tests are performed by qualified personnel:

- The following load resistors need to be prepared before the test: 100R, 1K, 1.2K each (the power no less than 1/2W)
- The 100R, 1 K load in series and connected to the stylus, the oscilloscope connected to both ends of the 100R resistor (See Circuit below), Turn on the device, with the output current set to 18 steps / 100mA, enable the output by pressing CONTINUOUS PULSE ON button. Check the peak voltage across the 100R resistor is 10V ± 20% (i.e. within 8V~12V), as expected from Ohm’s Law (10V=100mA x 100Ohms).
- Switch the 200µs/500µs slide switch in the battery compartment and observe whether the pulse width of the oscilloscope waveform can switch 200µS/500µS.
- Replace the 1K load resistor with 1.2K, press CONTINUOUS PULSE ON, and a triangle symbol is flashing on the LCD.

Note: Use the test stylus tips to perform the “output function checking”

MS100 CONTROLLER OPERATING INSTRUCTIONS

Switching the unit On and Off

- Press and hold the button identified as  for three seconds to turn the controller On or Off.

NOTES:

- The stimulation amplitude will always start at zero when the controller power is cycled on or off.
- The pulse mode will always reset to SINGLE PULSE when the power is cycled on the device.

Selecting the MYOTouch mode of operation:

- The MYOTouch Muscle Stimulator has two modes that can be used to identify striated muscles. Either the Single Pulse setting, or the Continuous Pulse setting can be used for muscle stimulation. The current output can also be adjusted on either of the pulse settings.

SINGLE PULSE – This pulse setting will deliver an instantaneous pulse of electrical current when the button is pushed. An indicator

TECHNICAL DATA

Model		MYOTouch
Type of Device	BF 	
Membrane Touch Switches	ON/OFF SINGLE PULSE CONTINUOUS PULSE ON CONTINUOUS PULSE OFF +/- (Output intensity adjustment) 200µs/500µs slide switch	
Pulse Characteristics Pulse Width: Frequency: Pulse Type:	200/500µs 50Hz Mono-phasic rectangular	
Output Current	Adjustable 0-112mA at load less than or equal to 1KΩ	
Max. Amplitude	56V@500Ω, 126V@2K, 126V@10K	
Display Pulse LED:	The PULSE LED will illuminate when is generated.	
Battery LED:	Steady <i>green</i> when the device is turned ON	
Battery Power	One 9V alkaline battery Battery symbol on LCD screen when the battery voltage is low.	
Power Supply	Internally powered ME equipment	
Protection Grade	IP20	
Operation Mode	Single Pulse or Continuous Pulse	
Size	113 mm x 61 mm x 25 mm	
Weight (without Battery)	105g	
+ / – (Switch)	Controls the output current	
Software Version	R6_00	

LIGHT-EMITTING DIODE (LED) INDICATORS

Feature Icon	Feature Description
LED color/ label	LED description
Yellow/ PULSE	The PULSE LED will illuminate when pulse is generated.
Green/ BATTERY	The BATTERY LED will illuminate when the device is turned on. Battery symbol on LCD screen when the battery voltage is low.


LIQUID CRYSTAL DISPLAY (LCD) ICONS

Icon Illustrated on the display	Icon Description
	The Battery icon illustrates the battery power level.
SINGLE PULSE	The Single Pulse icon illustrates the device is in the Single Pulse mode.
CONTINUOUS PULSE	The Continuous Pulse icon illustrates the device is in the Continuous Pulse mode.
	The Impedance feature will illuminate on the LCD when the device detects a high resistance.
Setting: X ("X" can be 0–20)	The LCD will display the device output setting.
Output: X mA ("X" can be 0–112)	The LCD will display the device output current level.

DEVICE ADJUSTABLE FEATURES

Feature Icon	Feature Description
	Power Switch – this button is used to switch the device ON and OFF. Press and hold the button for 3 seconds to cycle the power on the device.
	Single Pulse Button – this button is used to deliver a momentary pulse of stimuli current.
	Continuous Pulse On Button – this button allows the device to deliver a continuous pulse of stimuli current.
	Continuous Pulse Off Button – this button stops the device from continuously delivering stimuli current. NOTE: <i>The Continuous Pulse mode must be turned off to switch the device to the Single Pulse mode.</i>

DEVICE ADJUSTABLE FEATURES

Feature Icon	Feature Description
	200µs/500µs slide switch – This switch is located within the battery compartment. You may move the slide to the desired Pulse width (200µs or 500µs)

DEVICE NON-ADJUSTABLE FEATURES

Feature Icon	Feature Description
	The Impedance feature will illuminate on the LCD when the device detects a high resistance (see detail Alarm load below).
Step	Alarm load (KΩ) (Tolerance: ±10%)
1	45
2	16
3	9.1
4	6.8
5	5
6	4
7	3.5
8	2.9
9	2.5
10	2.3
11	2.1
12	2
13	1.8
14	1.6
15	1.5
16	1.4
17	1.3
18	1.2
19	1.1
20	1

of “SINGLE PULSE” will illuminate on the controller display only when the Single Pulse button is pressed.

CONTINUOUS PULSE ON – This pulse setting will continuously deliver electrical current when the “CONTINUOUS PULSE ON” button is pressed on the controller. An indicator of “CONTINUOUS PULSE” will illuminate on the controller display when the Continuous Pulse setting is activated. The device will stay in the Continuous Pulse mode until the Continuous Pulse setting is deactivated by pressing the “CONTINUOUS PULSE OFF” button on the device controller.

CONTINUOUS PULSE OFF – This button on the controller will deactivate the Continuous Pulse setting so the continuous electrical output will stop. The Continuous Pulse setting will need to be de-activated for switching the device from Continuous Pulse to Single Pulse.

200µs/500µs slide switch – This switch is located within the battery compartment. You may move the slide to the desired Pulse width (200µs or 500µs).

Adjusting the stimulation amplitude:

- The MYOTouch Muscle Stimulator can adjust the current output (stimulation amplitude). The device will display the output for each setting identified on the LCD display.

Output Setting	Output (stimulation amplitude)
1	6mA
2	11mA
3	17mA
4	22mA
5	28mA
6	34mA
7	39mA
8	45mA
9	50mA
10	56mA
11	62mA
12	67mA
13	73mA
14	78mA
15	84mA
16	90mA
17	95mA
18	100mA
19	106mA
20	112mA

- The device output setting will reset to 0 when the device power is turned ON or OFF.
- The device output setting is adjusted with the Output Intensity button on the side of the device. The Output Intensity is identified by “+” (plus) and “-” (minus) symbols printed on the switch.
- The plus (+) side of the switch will increase the output setting. The minus (-) side of the switch will decrease the output setting.
- If the output setting is set to zero (0) then the device will not deliver any electrical output current.
- Please be sure to review the “Warning/Caution Statements” prior to operating the device.

ACCESSORIES

The accessories for the MYOTouch Muscle Stimulator are as follows:

MS200 - Sterile Stylus specifications:

- AirLife Stylus part number = MS200
- Approximate Stylus handle diameter = 14 mm
- Approximate Stylus handle length = 130 mm
- Approximate wire length = 1400 mm
- Stylus tip diameter = 1.4 mm + 0.05mm
- Approximate Stylus tip length = 25 mm
 - Length exposed from the stylus handle.

Replacement battery specifications:

- 9 Volt Alkaline battery

WARNING:

- Please do not use accessories and/or detachable parts not specified or authorized by the manufacturer.
- The use of accessories and cables other than those specified by AirLife could result in increased electromagnetic emissions or decreased electromagnetic immunity of this equipment and result in improper operation.

STERILE STYLUS ATTACHMENT INSTRUCTIONS:

- Remove the stylus packaging tray from the carton.
- Peel the lid from the stylus packaging tray.
- Remove the stylus from the packaging tray.
- Plug the stylus into the MYOTouch Muscle Stimulator device.

- The stylus plug is bi-polar so the stylus will operate as intended in either orientation.
- Please align the stylus plug receptacles with the mating pins in the MYOTouch device.

OPERATING INSTRUCTIONS FOR THE MYOTOUCH DEVICE WITH STYLUS ATTACHMENT (MS100 WITH MS200)

- Turn on the power for the MYOTouch muscle stimulator.
- Attach the probe to the MYOTouch controller if it is not already attached.
- Increase the intensity output setting to 1.
- The MYOTouch controller will default to the Single Pulse when turned on. Please select the desired mode for operation (Single Pulse or Continuous Pulse).

CAUTION: The stylus tips will continuously discharge electrical energy when the device is in the Continuous Pulse mode.

- Increase the stimuli output setting to the desire level. If unsure of the needed output, then start at 1 and continue to increase the output to achieve the desired stimuli results.
- When the operator is ready to proceed with stimulation it is usual for the area of interest to be swabbed with saline in order to reduce the electrical resistance between the bipolar probe tips and the tissue of interest.
- To stimulate the muscle touch both tips of the stylus on the muscle. The duration of the tip to muscle contact time only needs to be 1 to 3 seconds.

WARNING: The stylus tips will discharge an electrical shock to any conductive material that they contact.

WARNING: To minimize the risk of thermal burns due to tissue heating, do not keep the probe in one location for any length of time (for example, level 18 for more than 20 seconds and level 20 for more than 1 second could result in tissue burn).

- Continue 1-3 seconds contacts with the muscle until the muscles are identified as desired.
- When completed with the device please turn the power off and discard the used stylus. The stylus is approved for single patient use only.

MAINTENANCE

Replacing the 9-volt battery is the only maintenance item for the MYOTouch device. The battery icon will flash when the battery needs to be replaced.

- The battery can be replaced by opening the battery cover on the back of the MYOTouch device.
- Pull the black ribbon to assist with the removal of the battery.
- Remove the used battery and discard the battery per local regulations.
- Place the black ribbon perpendicular to the long side of the battery. The ribbon should be placed so it can be utilized to assist with the next battery removal process.
- Install the replacement battery per the illustration in the battery compartment.

6. Please be sure the battery is pressed firmly into the battery compartment.
7. Install the battery cover and be sure the black ribbon is tucked into the battery compartment before completely installing the battery cover.
8. The trained professional or pediatric surgeon is the intended operator.
9. Do not service the device while in use.
10. Maintenance and repairs should only be carried out by an authorized agency. The manufacturer will not be held responsible for repairs by unauthorized personnel.
11. The user must not attempt any repairs to the device or accessories.
12. Opening of the equipment by unauthorized agencies is not allowed, and will terminate the right to warranty.
13. If necessary, we will provide circuit diagrams, component part lists or other information that will assist authorized service personnel to repair the device.

WARNING: Don't use accessories and detachable parts not specified or authorized by manufacturer. Otherwise, it may cause damage to the unit or danger to the user or patients. Further information on purchasing accessories can be obtained at info@myAirLife.com.

NOTE: Do not use the device if a battery has leaked acid into the battery compartment. The leaked battery acid could impair the internal circuits for the device. Please contact AirLife to inquire about a device replacement.

CLEANING

 If the stimulator requires cleaning, it should be unplugged from the mains power supply and wiped with a damp (not wet) cloth.

CAUTION: Excessive liquids can damage the device. Please make sure the cleaning cloth is only damp.

 Frykman-Kimble Stimulator Stylus (Part Number: MS200) has been EO sterilized. The Stimulator Stylus is a sterile single patient use device, it does not require cleaning.

CAUTION:

- If the Stimulator Stylus packaging is damaged, please do not use it.
- Do not reuse the Stimulator Stylus (Model: MS200).

 During use the stylus tips can be wiped off with a dry non-conductive lint free wipe.

CAUTION: The stylus can discharge electrical energy if it contacts a conductive material with the device turned on.

SERVICE LIFE

 Please refer to the label information for the expiration date of the product.

 The expected shelf life of Frykman-Kimble Stimulator Stylus (Model: MS200) is 2 years. The expiry date will be marked on the packaging of Stimulator Stylus.

ENVIRONMENTAL CONDITIONS FOR USE AND STORAGE

The recommended environmental conditions for use and storage are as listed below:

Operating Temperature	+5°C to +40°C
Operating Humidity	15% RH to 90% RH
Storage Temperature Range	-10°C to +60°C
Storage Humidity Range	15% RH to 90% RH
Operating & Storage Atmospheric Pressure	700hPa to 1,060hPa

TROUBLESHOOTING

If your machine is not working properly please check the following:

Please contact the manufacturer if the troubleshooting suggestions do not resolve the problem.

Problem	Possible Causes	Solution
No output or a blank display	Low Battery	Replace Battery
	Battery inserted incorrectly	Check if battery inserted correctly
	Damaged battery contact terminal in battery compartment	Contact supplier
Battery light flashes	Low battery	Replace battery
No Muscle stimulation	Inadequate connection between the stylus and the MYOTouch controller	Verify the stylus plug is fully inserted into the receptable on the MYOTouch controller
	Stimuli setting is too low	Increase the stimuli output setting

SAFE DISPOSAL

- Dispose of the device per local guidelines and regulations for electrical devices
- As the Stimulator Stylus (MS200) is a single use item, it should be disposed of in accordance with local clinical disposal procedures.

EMC PRECAUTIONS

Wireless communication equipment such as mobile phones, cordless telephones and, walkie-talkies can affect the device. Wireless communication devices should be kept at a distance of at least 3.3 meters from the MYOTouch device.

(NOTE: As indicated in Table 6 of IEC60601-1-2 for ME EQUIPMENT, a typical cell phone with a maximum output power of 2 W yields d = 3.3 m at an IMMUNITY LEVEL of 3 V/m).

The essential performances are listed below:

1. When the output current is set to 100mA (output setting = 18) at resistance = 1200 Ohm (tolerance ±10%), press the SINGLE PULSE or CONTINUOUS PULSE ON button, the triangle symbol should be shown on the LCD screen, otherwise, it will not be acceptable;
2. When the output current is set to 112mA (output setting = 20) at resistance = 500 Ohm, press the SINGLE PULSE or CONTINUOUS PULSE ON button, the output voltage is within 56V±20%. It is not acceptable if the output voltage is not in the range, or triangle symbol appears.

WARNING: Use of this equipment adjacent to or stacked with other equipment should be avoided because it could result in improper operation. If such use is necessary, this equipment and the other equipment should be observed to verify that they are operating normally.

WARNING: Use of accessories, transducers, and cables other than those specified or provided by the manufacturer of this equipment could result in increased electromagnetic emissions or decreased electromagnetic immunity of this equipment and result in improper operation.

WARNING: Portable RF communications equipment (including peripherals such as antenna cables and external antennas) should be used no closer than 30 cm (12 inches) to any part of the ME equipment, including cables specified by the manufacturer. Otherwise, degradation of the performance of this equipment could result.

Table 1	
Declaration - electromagnetic emission	
Emissions test	Compliance
RF emissions CISPR 11	Group 1
RF emissions CISPR 11	Class B
Harmonic emission IEC 611000-3-2	Not applicable
Voltage fluctuations / flicker emissions IEC 61000-3-3	Not applicable

Table 2		
Declaration - electromagnetic immunity		
Immunity test	IEC 60601 test level	Compliance Level
Electrostatic discharge (ESD) IEC 61000-4-2	±8 kV contact ±2 kV, ±4 kV, ±8 kV, ±15 kV air	±8 kV contact ±2 kV, ±4 kV, ±8 kV, ±15 kV air
Electrical fast transient/ burst IEC 61000-4-4	±2kV for power supply lines ± 1 kV for input/output lines	Not applicable
Surge IEC 61000-4-5	±0.5kV, ± 1 kV line(s) to lines ±0.5kV, ± 1kV, ± 2 kV line(s) to earth	Not applicable
Voltage dips, short interruptions and voltage variations on power supply input lines IEC 61000-4-11	0 % UT; 0.5 cycle At 0°, 45°, 90°, 135°, 180°, 225°, 270° and 315° 0% UT; 1 cycle and 70% UT; 25/30 cycles Single phase: at 0° 0% UT; 250/300 cycles	Not applicable
Power frequency (50/60 Hz) magnetic field IEC 61000-4-8	30 A/m	30 A/m
NOTE: UT is the a.c. mains voltage prior to application of the test level.		

Table 3		
Declaration - electromagnetic immunity		
Immunity test	IEC 60601 test level	Compliance level
Conducted RF IEC 61000-4-6	3 V 0.15 MHz to 80 MHz 6V in ISM bands between 0.15 MHz and 80 MHz	Not applicable
Radiated RF IEC 61000-4-3	10V/m 80 MHz to 2.7 GHz	10V/m

Table 4					
Declaration - IMMUNITY to proximity fields from RF wireless communications equipment					
Immunity test	IEC 60601 test level				Compliance level
	Test Frequency	Modulation	Maximum power	Immunity level	
Radiated RF IEC 61000-4-3	385 MHz	**Pulse modulation: 18 Hz	1.8 W	27 V/m	27 V/m
	450 MHz	*FM+ 5 Hz deviation: 1 kHz sine	2 W	28 V/m	28 V/m
	710 MHz 745 MHz 780 MHz	**Pulse Modulation: 217 Hz	0.2 W	9 V/m	9 V/m
	810 MHz 870 MHz 930 MHz	**Pulse Modulation: 18 Hz	2 W	28 V/m	28 V/m
	1720 MHz 1845 MHz 1970 MHz	**Pulse Modulation: 217 Hz	2 W	28 V/m	28 V/m
	2450 MHz	**Pulse Modulation: 217 Hz	2 W	28 V/m	28 V/m
	5240 MHz 5500 MHz 5785 MHz	**Pulse Modulation: 217 Hz	0.2 W	9 V/m	9 V/m
Note* - As an alternative to FM modulation, 50% pulse modulation at 18 Hz may be used because while it does not represent actual modulation, it would be the worst case. Note** - The carrier shall be modulated using a 50% duty cycle square wave signal.					

INCIDENT REPORTING

Notice to the user and/or patient that any serious incident that has occurred in relation to the device should be reported to the manufacturer and the competent authority of the Member State in which the user and/or patient is established.


EasyMed Instruments Co., LTD.
3F-6F, Block A, No.4, Fengxin Road
Fengxiang Industrial District
Daliang, Shunde, Foshan, 528300 Guangdong
China

Distributed Exclusively by:
AirLife
2710 Northridge Dr. NW, Suite A
Grand Rapids, MI 49544 USA
www.myAirLife.com
Made in China

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