Valved Aerosol Mask

EN Valved Aerosol Mask

DEVICE DESCRIPTION: The device is non-sterile, disposable, and for single-patient use. The aerosol valved mask includes one-way valves placed over the exhalation hole to aid in medication delivery.

INTENDED PURPOSE: Aerosol masks are intended to be used for aerosolized medication delivery applications.

INDICATIONS: Aerosol mask is used in conjunction with a nebulizer to deliver aerosolized medication to the mouth and nose in spontaneously breathing patients.

ENVIRONMENT: Hospitals, sub-acute, medical clinics, pre-hospital, home, surgical centers, skilled nursing facilities.

PATIENT TARGET GROUP: Pediatric, adult.

EXPECTED CLINICAL BENEFITS: When used to deliver bronchodilators, clinical benefit is to relieve bronchospasms, decrease the frequency of exacerbations and improve lung function.

CONTRAINDICATIONS: None Known.

WARNINGS:

- If oxygen is in use, do not use near flame or heat source.
- Device contains components that may present choking hazard.
- Use a flow rate ≥ 5 LPM to prevent CO2 accumulation and re-breathing.
- Contains Octamethylcyclotetrasiloxane (D4) above 0.1% w/w.

CAUTIONS:

- Clinician or caregiver should avoid inhaling aerosolized medication.
- Consult physician if patient develops infection, skin irritation, or material sensitization.

RESIDUAL RISKS: Refer to Warnings and Cautions.

NECESSARY USER QUALIFICATIONS: This device is restricted to sale by or on the order of a physician. There are no additional requirements for special facilities, special training, or particular qualifications for the use of this device.

INSTRUCTIONS: Follow instructions for proper care of device to prevent risk of contamination and injury.

- Place mask over nose and mouth.
 - Adjust elastic strap and nose bar for snug fit.
 - Insert connector into nebulizer or connect to aerosol tubing.
 - Set prescribed flow rate.
 - Verify aerosol mist is visible.
- This product is disposable and is not intended for prolonged use greater than seven days.

SAFE DISPOSAL: Dispose of device in accordance with local, state or national regulations.

INCIDENT REPORTING: Contact your country’s Competent Authority and Manufacturer to report any serious incident.

PERFORMANCE CHARACTERISTICS:

- For use with Flow Rate 5 - 8 LPM
- For use with 22 mm connector

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