**IFU SPECIFICATIONS ARE NOT INTENDED TO PRINT.
THE SPECIFICATIONS BOX IS FOR REFERENCE ONLY.
PLEASE REMOVE BEFORE PRINTING ARTWORK.**

IFU SPECIFICATION BOX

AirLife FG:	MW100-0-120	Customer FG#:	N/A
RLS:	RLS-20011	RIS:	RIS-105000
DCR #:	DCR-17960		
Document Name:	105000, IFU - Salter Mini Whiskers MW100		
DocumentType: (Select by filling box with black)	☒ Paper Insert / IFU ☐ Printed Pouch ☐ Printed Box ☐ Product Imprint ☐ Label ☐ Other:		
Dimensions:	Flat:	11" (W) x 17" (H) ±0.5"	
	Folded:	2.75" (W) x 5.67" (H)	
Material Notes:	White paper		
Binding: (Select by filling box with black)	☐ Stapled:	(if applicable, note location of staple here)	☐ Folded ☐ Other:
Ink Color: (Fill box with ink color and note PMS or CMYK. Delete any extra boxes.)	☒ Black		
Printer Information	Printer (Supplier or Site):	Team Tec	
Document Control Number:	105000	Revision Number:	H
Graphic Designer:	Grace Contreras	Date:	2025/11

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Securement Device

EN English

DEVICE DESCRIPTION:

Mini Whiskers Securement Devices use a gentle hydrocolloid adhesive and hook and loop strips to secure nasal cannulas or feeding tubes on premature infants. Hydrocolloid provides gentle adhesion for 24-36 hours. Mini Whiskers are non-sterile, single-use, disposable devices.

INTENDED PURPOSE:

Mini Whiskers securement device aids the placement and securement of nasal cannulas and feeding tubes in infants less than 2000g. Not for use with endotracheal tubes.

ENVIRONMENT FOR USE:

Acute care environments, e.g. hospital, NICU.

MR SAFETY:

MR Safe

INTENDED USERS:

Trained medical professionals.

EXPECTED CLINICAL BENEFITS:

The device suspends the nasal cannula or feeding tube away from the infant's skin, which

minimizes friction and pressure preventing pressure sores and tissue damage.

By securely holding the tubing in place, the device reduces the risk of accidental dislodgement. This ensures continuous and uninterrupted therapy, such as oxygen or nutrition delivery.

CONTRAINDICATIONS:

None known

WARNINGS:

- Do not use to secure endotracheal tubes.
- Device contains components that may present choking hazard.
- Not for use over broken skin or sores.

CAUTIONS:

- Do not use if known sensitivity to hydrocolloids.
- If skin irritation is noted, discontinue use.
- This product is disposable and is not intended for prolonged or extended use.

RESIDUAL RISKS:

Residual risks are covered in the warning and caution statements.

NECESSARY USER QUALIFICATIONS:

The end-user should receive training on the proper use of the device and be familiar with the indications for use, warnings, and cautions as stated in the IFU.

Rx Only: Federal law (USA) restricts this device to sale by or on the order of a licensed healthcare practitioner.

INSTRUCTIONS:

- Prepare the skin. Clean infant's skin with water and dry completely.
- Remove adhesive loop strip from the protective liner. Wrap the loop strip wrap around the cannula or feeding tube in a candy cane style or flag style.
- Remove Mini Whiskers Securement Device from the protective liner and align the hydrocolloid adhesive side onto the skin. Gently press down on device to adhere to the skin.
- Secure the tubing. Attached the loop-wrapped tubing to the hook strip on the Mini Whiskers Securement Device. Ensure the hook and loop surfaces engage and hold tubing in place.
- To reposition the tubing, gently pull the wrapped tubing away from the hook material and realign as needed.
- Change every 36 hours, or sooner if the device starts to lift away from skin or if skin irritation is observed.

Note: Use your hands to gently warm the hydrocolloid strip before applying to improve the adhesion.


SAFE DISPOSAL:

Dispose of all materials in accordance with local, state, and federal regulations.

INCIDENT REPORTING:

Notice to the user that any serious incident that has occurred in relation to the device should be reported to the manufacturer at productquality@myairlife.com, and the competent authority of the Member State in which the user and/or patient is established.

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MT Promedt Consulting Ltd.
MT Promedt Consulting GmbH
MedEnvoy UK Limited