












SYSTÈME DE CODE COULEUR HANDTEVY POUR AIR-Q®SP3G

ARTICLE	TAILLE	PCI (LIVRES)	PCI (KG)	HANDTEVY	OETT MAXIMUM	OUVERTURE DE LA BOUCHE	LONGUEUR DE LA SONDE	SONDE OG MAXIMUM	CONNECTEUR AIR Q3
2114-60005	0	< 4 livres	< 2 kg	—	3,0 mm	5 mm	6.86 cm	5 Fr	
2114-60055	0,5	4 à 9 livres	2 à 4 kg	2 kg	4,0 mm	8 mm	8.05 cm	6 Fr	
2114-60055	0,5	4 à 9 livres	2 à 4 kg	3 à 5 kg	4,0 mm	8 mm	8.05 cm	6 Fr	
2114-60105	1	9 à 15 livres	4 à 7 kg	6 kg	4,5 mm	11 mm	9.29 cm	8 Fr	
2114-60155	1,5	15 à 37 livres	7 à 17 kg	7 à 17 kg	5,0 mm	14 mm	10.9 cm	8 Fr	
2114-60155	1,5	15 à 37 livres	7 à 17 kg	7 à 17 kg	5,0 mm	14 mm	10.9 cm	8 Fr	
2114-60155	1,5	15 à 37 livres	7 à 17 kg	7 à 17 kg	5,0 mm	14 mm	10.9 cm	8 Fr	
2114-60155	1,5	15 à 37 livres	7 à 17 kg	7 à 17 kg	5,0 mm	14 mm	10.9 cm	8 Fr	
2114-60205	2	37 à 66 livres	17 à 30 kg	17 à 30 kg	5,5 mm	17 mm	11.61 cm	10 Fr	
2114-60205	2	37 à 66 livres	17 à 30 kg	17 à 30 kg	5,5 mm	17 mm	11.61 cm	10 Fr	
2114-60205	2	37 à 66 livres	17 à 30 kg	17 à 30 kg	5,5 mm	17 mm	11.61 cm	10 Fr	
2114-60305	3	66 à 132 livres	30 à 60 kg	30 à 60 kg	7,0 mm	20 mm	14.25 cm	14 Fr	