

REF 0556

UNITS - 1
LENGTH-7' (2.1 m)

RX Only



CE **E0482**



SUGGESTED PROCEDURE:

1. These products are intended for use in delivering oxygen therapy.
2. All products are to be used in accordance with the instructions for use.
3. Remove the product from the package.
4. Attach the device to the oxygen source and make appropriate connection.
5. Ensure that all connections are secure, then start the gas flow.
6. confirm appropriate function prior to use.
7. check for gas flow throughout the device.
8. Discard after single patient use.



CAUTION:

Use of sterilization or cleaning solutions may leave this device nonfunctional.

Patents Pending

Westmed 

5580 S. Nogales Highway
Tucson, Arizona 85706 USA
Phone: 800-975-7987
Fax: 520-294-6061
www.westmedinc.com

EC **REP**

MT Promed Consulting GmbH
Altenhofstr. 80
66386 St. Ingbert, Germany
Product of Mexico
Assembled in Mexico
LABEL PN 77600, Rev. 05



LOT

Comfort Soft Plus®
ADULT CANNULA
CÁNULA PARA ADULTOS
CÁNULA DE ADULTO
CANULE POUR ADULTE
CANNULA ADULTI
KANÜLEFÜR ERWACHSENE