

REF 6159

**Unit - 1
Rx Only**



SUGGESTED PROCEDURE:

1. Filters may be fitted at patient end or machine end of circuit.

2. Change the filter as required or as per protocol.

CAUTION: For single patient use only. Use of sterilization or cleaning solutions may leave this device nonfunctional. Store in a cool dry place.

Do not use between the patient and a nebulizer or humidifier.

Westmed 

5580 S. Nogales Highway
Tucson, Arizona 85706 USA
Phone: 800-975-7987
Fax: 520-294-6061
www.westmedinc.com



MT Promedt Consulting GmbH
Altenhofstr. 80
66386 St. Ingbert, Germany
Assembled in MEXICO
LABEL PN 78085, Rev. 05



LOT

HEPA FILTER w/CO2 Port

22mm OD/15mm ID x 22mm ID

99.977 % @ 30 LPM Filter Efficiency

99.999 % BFE, 99.999 % VFE

55 ml Dead Space

1.59 cmH₂O @ 30 LPM Flow Resistance

3.55 cmH₂O @ 60 LPM Flow Resistance