

**REF** 0360

UNITS - 1

**RX Only**



CE 0482

 **CAUTION:**

Use of sterilization or cleaning solutions may leave this device nonfunctional.

**SUGGESTED USE:**

1. These products are intended for use in delivering oxygen therapy.
2. All products are to be used in accordance with the applicable prescription.
- 3) Attach the Mask to the supply tube, and then connect the supply tube to the gas source.
- 4) Insure that all connections are secure, and then start gas flow.
- 5) Check for appropriate function prior to use.
- 6) Discard after single patient use.

**Westmed** 

5580 S. Nogales Highway  
Tucson, Arizona 85706 USA  
Phone: 800-975-7987  
Fax: 520-294-6061  
[www.westmedinc.com](http://www.westmedinc.com)



MT Promedt Consulting GmbH  
Altenhofstr. 80  
66386 St. Ingbert, Germany  
Assembled in Mexico  
LABEL PN 78745, Rev. 01

**LOT**

TRACHEOSTOMY MASKS  
MASCARILLAS DE TRAQUEOSTOMÍA  
MÁSCARAS DE TRAQUEOSTOMIA  
MASQUES DE TRACHÉOTOMIE  
MASCHERE TRACHEOSTOMICHE  
TRACHEOTOMIEMASKEN  
ADULT, ADULTOS, ADULTO,  
ADULTE, ADULTI, ERWACHSENE