


**REF** 7-8113-10-62-68

**UNITS - 1**  
**RX Only**



**DISPOSABLE ANESTHESIA  
BREATHING CIRCUITS**  
Anesthesia, Adult, 96"

 **CAUTION:** Use of sterilization or cleaning solutions may leave this device nonfunctional.

**SUGGESTED USE:**

1. Attach the circuit pursuant to standard operating protocol.
2. Prior to patient connection: Check that all connections are secure. Fully extend the tubing (if using expandable circuit). Pressurize to confirm there is no leakage or occlusion of the circuit.

**Westmed** 

5580 S. Nogales Highway  
Tucson, Arizona 85706 USA  
Phone: 800-975-7987  
Fax: 520-294-6061  
www.westmedinc.com



MT Promed Consulting GmbH  
Altenhofstr. 80  
66386 St. Ingbert, Germany  
Patent Pending

Assembled in Mexico  
LABEL PN 79082, Rev. 01



**LOT**