

**REF 9147**

# DISPOSABLE RESPIRATORY CIRCUITS

72" Adult Single Limb Circuit w/Exhalation Valve

**UNITS - 1**

**RX Only**



Not Made  
With Natural  
Rubber Latex



Not Made  
With DEHP



**CE 0482**

**⚠ CAUTION:** Use of sterilization or cleaning solutions may leave this device nonfunctional.

## SUGGESTED USE:

1. Attach the circuit pursuant to standard operating protocol.
2. Prior to patient connection:
  - a) Check that all connections are secure.
  - b) Fully extend the tubing (if using expandable circuit).
  - c) Pressurize to confirm there is no leakage or occlusion of the circuit.



MT Promedt Consulting GmbH  
Altenhofstr. 80  
66386 St. Ingbert, Germany

**Westmed** 

5580 S. Nogales Highway  
Tucson, Arizona 85706 USA

Phone: 800-975-7987

Fax: 520-294-6061

www.westmedinc.com

Assembled in Mexico

LABEL PN 79268, Rev. 01



**LOT**