

**REF** 100619

**UNITS - 1**  
**RX Only**



# HEART<sup>®</sup> DIRECT CONNECT CONTINUOUS NEBULIZER

Adult Kit

**! CAUTION:** Use of sterilization or cleaning solutions may leave this device nonfunctional.

- Contains One Kit Per Package
- Disposable Non-Sterile
- Clean Ready for Use

This Kit Contains:

1 each 100610 HEART<sup>®</sup> Direct Connect

Continuous Nebulizer

1 each 0250 Aerosol Mask

1 each 8' x 22 mm Flex Hose (cuffed every 6")

Aerosol performance is based on using a 50 psi pressure source.

**Westmed**

5580 S. Nogales Highway  
Tucson, Arizona 85706 USA  
Phone: 800-975-7987

Fax: 520-294-6061  
[www.westmedinc.com](http://www.westmedinc.com)

Made in U.S.A.

LABEL PN 79509, Rev. 01



MT Promedt Consulting GmbH  
Altenhofstr. 80  
66386 St. Ingbert, Germany



**LOT**