

**REF** 101021

**IV-HEART® CONTINUOUS NEBULIZER**

Adult Kit

**UNITS - 1**  
**RX Only**



Not Made  
With Natural  
Rubber Latex



Not Made  
With DEHP



**CE** 0482

- Contains One Kit Per Package
- Disposable Non-Sterile
- Clean Ready for Use

This Kit Contains:

- 1 each 101020 IV-HEART® Continuous Nebulizer
  - 1 each 0250 Adult Aerosol Mask
  - 1 each 7' Oxygen Tubing w/Threaded Nut
  - 1 each 8' x 22 mm Flex Hose (cuffed every 6")
- Aerosol performance is based on using a 50 psi pressure source.



MT Promedt Consulting GmbH  
Altenhofstr. 80  
66386 St. Ingbert, Germany

**Westmed**

5580 S. Nogales Highway  
Tucson, Arizona 85706 USA  
Phone: 800-975-7987  
Fax: 520-294-6061  
www.westmedinc.com  
Made in U.S.A.  
LABEL PN 79547, Rev. 01



**LOT**