















Nr etapu	Skuteczna odseteczka (µm)	NebuTech Nebulizator do ostrej opieki		
		Srednia	Próbka 1	Próbka 3
8	0.00	0.0%	0.0%	0.0%
7	0.68	12.5%	11.1%	15.3%
6	0.90	17.4%	15.3%	21.5%
5	1.22	38.0%	32.4%	39.0%
4	2.87	44.5%	41.9%	47.2%
3	4.85	63.6%	61.7%	62.0%
2	7.50	76.2%	72.8%	75.1%
1	8.27	80.7%	76.7%	85.3%



Salter Labs
 30 Spur Drive
 El Paso, TX 79906 USA
 Made in Mexico
 www.salterlabs.com



MT Promedt Consulting GmbH
 Altenhofstrasse 80
 66386 St. Ingbert
 Germany