

SAFET™

T-PIECE RESUSCITATOR



Deliver consistent PIP and PEEP with every breath.

The SafeT™ T-Piece Resuscitator is the safe solution for neonatal resuscitation and CPAP. Delivering continuous oxygen flow alongside pre-set PIP and PEEP pressures enables clinicians to safely provide the necessary ventilatory support during resuscitation of infants. Designed specifically for use on the delicate lungs of premature babies and newborns under 10 kg/22 lbs.

Safely ventilate the fragile lungs of neonates.

Provide consistent ventilatory pressure.

SafeT T-Piece Evaluations Available

Request a sample or product demonstration for your NICU + labor/delivery team today.



**Pre-set
PIP and PEEP
settings prior
to use.**



Integrated, color-coded pressure manometer for effective monitoring of ventilatory pressures.



PEEP dial allows clinician to set PEEP; ventilation is provided by occluding PEEP dial hole intermittently at the desired respiratory rate allowing for inspiration and exhalation to occur.



Use optional tethered override cap on inspiratory pressure controller when pressures greater than 40 cm H₂O are needed.

SIDE-BY-SIDE

Current NRP guidelines recommend the use of either self-inflating, flow-inflating or T-piece resuscitators for manual ventilation of neonates. This side-by-side comparison highlights the additional features of a T-piece resuscitator compared to a BVM and hyperinflation system.



NEONATAL RESUSCITATOR FEATURES + FUNCTIONALITY

	SAFET™	BVM	HYPERINFLATION
Continuous O2 Flow	●		●
Pre-Set + Deliver Consistent PIP	●		
Pre-Set + Deliver Consistent PEEP	●	●	●
Can Provide CPAP	●		●
Pre-set Maximum Pressure + Adjustability During Resus.	●		
Integrated Pressure Manometer	●	●	●
Pressure Relief Valve	●	●	●
Hand Compression Required		●	●
Portable	●	●	●

CONFIGURATIONS



SafeT™ w/ Pop Off Override

SAFET™ T-PIECE RESUSCITATOR

ITEM	DESCRIPTION	MASK SIZE	PK
ST710	SafeT T-Piece Resuscitator w/ Pop Off	X-Small Neonate Mask	10
ST711	SafeT T-Piece Resuscitator w/ Pop Off	Small Neonate Mask	10
ST712	SafeT T-Piece Resuscitator w/ Pop Off	Small Infant Mask	10
ST713	SafeT T-Piece Resuscitator w/ Pop Off	Infant Mask	10
ST714	SafeT T-Piece Resuscitator w/ Pop Off	Neonate Round Mask	10
ST715	SafeT T-Piece Resuscitator w/ Pop Off	Infant Round Mask	10
ST717	SafeT T-Piece Resuscitator w/ Pop Off	No Mask	10

