


SUNMED®

ENDOTRACHEAL TUBES



SunMed® Endotracheal Tubes

SunMed® endotracheal tubes are designed for either oral or nasal use with a smooth Murphy Eye as an alternative ventilation port. Made of latex free, flexible PVC, that are available cuffed, uncuffed, with or without stylet. Cuffed tubes feature a high volume, low pressure barrel cuff. A radiopaque strip is embedded to assess exact location of tube. All have a standard 15 mm ISO connector. These airways are individually packaged.



SunMed® Endotracheal Tubes

Family of Products



ENDOTRACHEAL TUBE - CUFFED

ITEM	I.D.	O.D.	LENGTH	SIZE	PK
1-7333-25	2.5 mm	3.5 mm	150 mm	10 Fr	10
1-7333-30	3.0 mm	4.2 mm	165 mm	12 Fr	10
1-7333-35	3.5 mm	4.9 mm	188 mm	14 Fr	10
1-7333-40	4.0 mm	5.5 mm	210 mm	16 Fr	10
1-7333-45	4.5 mm	6.2 mm	225 mm	18 Fr	10
1-7333-50	5.0 mm	6.9 mm	245 mm	20 Fr	10
1-7333-55	5.5 mm	7.5 mm	275 mm	22 Fr	10
1-7333-60	6.0 mm	8.2 mm	285 mm	24 Fr	10
1-7333-65	6.5 mm	8.8 mm	295 mm	26 Fr	10
1-7333-70	7.0 mm	9.6 mm	310 mm	28 Fr	10
1-7333-75	7.5 mm	10.2 mm	315 mm	30 Fr	10
1-7333-80	8.0 mm	10.9 mm	330 mm	32 Fr	10
1-7333-85	8.5 mm	11.5 mm	330 mm	34 Fr	10
1-7333-90	9.0 mm	12.1 mm	330 mm	36 Fr	10
1-7333-95	9.5 mm	12.7 mm	330 mm	38 Fr	10
1-7333-10	10.0 mm	13.6 mm	330 mm	40 Fr	10




ENDOTRACHEAL TUBE - CUFFED WITH STYLETTE

ITEM	I.D.	O.D.	LENGTH	SIZE	PK
1-7343-30	3.0 mm	4.2 mm	165 mm	12 Fr	10
1-7343-35	3.5 mm	4.9 mm	188 mm	14 Fr	10
1-7343-40	4.0 mm	5.5 mm	210 mm	16 Fr	10
1-7343-45	4.5 mm	6.2 mm	225 mm	18 Fr	10
1-7343-50	5.0 mm	6.9 mm	245 mm	20 Fr	10
1-7343-55	5.5 mm	7.5 mm	275 mm	22 Fr	10
1-7343-60	6.0 mm	8.2 mm	285 mm	24 Fr	10
1-7343-65	6.5 mm	8.8 mm	295 mm	26 Fr	10
1-7343-70	7.0 mm	9.6 mm	310 mm	28 Fr	10
1-7343-75	7.5 mm	10.2 mm	315 mm	30 Fr	10
1-7343-80	8.0 mm	10.9 mm	330 mm	32 Fr	10
1-7343-85	8.5 mm	11.5 mm	330 mm	34 Fr	10
1-7343-90	9.0 mm	12.1 mm	330 mm	36 Fr	10
1-7343-95	9.5 mm	12.7 mm	330 mm	38 Fr	10
1-7343-10	10.0 mm	13.6 mm	330 mm	40 Fr	10



ENDOTRACHEAL TUBE - UNCUFFED

ITEM	I.D.	O.D.	LENGTH	SIZE	PK
1-7330-20	2.0 mm	2.9 mm	145 mm	8 Fr	10
1-7330-25	2.5 mm	3.5 mm	150 mm	10 Fr	10
1-7330-30	3.0 mm	4.2 mm	165 mm	12 Fr	10
1-7330-35	3.5 mm	4.9 mm	185 mm	14 Fr	10
1-7330-40	4.0 mm	5.5 mm	210 mm	16 Fr	10
1-7330-45	4.5 mm	6.2 mm	225 mm	18 Fr	10
1-7330-50	5.0 mm	6.9 mm	245 mm	20 Fr	10
1-7330-55	5.5 mm	7.5 mm	275 mm	22 Fr	10
1-7330-60	6.0 mm	8.2 mm	285 mm	24 Fr	10
1-7330-65	6.5 mm	8.8 mm	295 mm	26 Fr	10
1-7330-70	7.0 mm	9.6 mm	310 mm	28 Fr	10



ENDOTRACHEAL TUBE - REINFORCED, CUFFED

ITEM	I.D.	O.D.	LENGTH	SIZE	PK
1-7363-35	3.5 mm	5.8 mm	195 mm	14 Fr	10
1-7363-40	4.0 mm	6.3 mm	215 mm	16 Fr	10
1-7363-45	4.5 mm	6.8 mm	240 mm	18 Fr	10
1-7363-50	5.0 mm	7.3 mm	256 mm	20 Fr	10
1-7363-55	5.5 mm	7.8 mm	286 mm	22 Fr	10
1-7363-60	6.0 mm	8.4 mm	295 mm	24 Fr	10
1-7363-65	6.5 mm	9.0 mm	301 mm	26 Fr	10
1-7363-70	7.0 mm	9.6 mm	315 mm	28 Fr	10
1-7363-75	7.5 mm	10.2 mm	325 mm	30 Fr	10
1-7363-80	8.0 mm	10.8 mm	337 mm	32 Fr	10
1-7363-85	8.5 mm	11.3 mm	337 mm	34 Fr	10
1-7363-90	9.0 mm	11.0 mm	337 mm	36 Fr	10
1-7363-95	9.5 mm	12.4 mm	337 mm	38 Fr	10
1-7363-10	10.0 mm	12.9 mm	337 mm	40 Fr	10




ENDOTRACHEAL TUBE - PREFORMED, CUFFED, NASAL

ITEM	I.D.	O.D.	LENGTH	SIZE	PK
1-7323-40	4.0 mm	5.5 mm	310 mm	16 Fr	10
1-7323-45	4.5 mm	6.2 mm	310 mm	18 Fr	10
1-7323-50	5.0 mm	6.9 mm	330 mm	20 Fr	10
1-7323-55	5.5 mm	7.5 mm	330 mm	22 Fr	10
1-7323-60	6.0 mm	8.2 mm	360 mm	24 Fr	10
1-7323-65	6.5 mm	8.8 mm	360 mm	26 Fr	10
1-7323-70	7.0 mm	9.6 mm	380 mm	28 Fr	10
1-7323-75	7.5 mm	10.2 mm	380 mm	30 Fr	10
1-7323-80	8.0 mm	10.9 mm	390 mm	32 Fr	10



ENDOTRACHEAL TUBE - PREFORMED, UNCUFFED, NASAL

ITEM	I.D.	O.D.	LENGTH	SIZE	PK
1-7322-30	3.0 mm	4.2 mm	160 mm	12 Fr	10
1-7322-35	3.5 mm	4.9 mm	180 mm	14 Fr	10
1-7322-40	4.0 mm	5.5 mm	200 mm	16 Fr	10
1-7322-45	4.5 mm	6.2 mm	220 mm	18 Fr	10
1-7322-50	5.0 mm	6.9 mm	240 mm	20 Fr	10
1-7322-55	5.5 mm	7.5 mm	270 mm	22 Fr	10



ENDOTRACHEAL TUBE - PREFORMED, CUFFED, ORAL

ITEM	I.D.	O.D.	LENGTH	SIZE	PK
1-7321-40	4.0 mm	5.5 mm	200 mm	16 Fr	10
1-7321-45	4.5 mm	6.2 mm	220 mm	18 Fr	10
1-7321-50	5.0 mm	6.9 mm	240 mm	20 Fr	10
1-7321-55	5.5 mm	7.5 mm	270 mm	22 Fr	10
1-7321-60	6.0 mm	8.2 mm	280 mm	24 Fr	10
1-7321-65	6.5 mm	8.8 mm	290 mm	26 Fr	10
1-7231-70	7.0 mm	9.6 mm	300 mm	28 Fr	10
1-7321-75	7.5 mm	10.2 mm	310 mm	30 Fr	10
1-7321-80	8.0 mm	10.9 mm	320 mm	32 Fr	10
1-7321-85	8.5 mm	11.5 mm	320 mm	34 Fr	10
1-7321-90	9.0 mm	12.1 mm	320 mm	36 Fr	10



ENDOTRACHEAL TUBE - PREFORMED, UNCUFFED, ORAL

ITEM	I.D.	O.D.	LENGTH	SIZE	PK
1-7320-30	3.0 mm	4.2 mm	160 mm	12 Fr	10
1-7320-35	3.5 mm	4.9 mm	180 mm	14 Fr	10
1-7320-40	4.0 mm	5.5 mm	200 mm	16 Fr	10
1-7320-45	4.5 mm	6.2 mm	220 mm	18 Fr	10
1-7320-50	5.0 mm	6.9 mm	240 mm	20 Fr	10
1-7320-55	5.5 mm	7.5 mm	270 mm	22 Fr	10