

# INFU-SURG® CLARITY

LATEX FREE | SINGLE PATIENT USE | DISPOSABLE | DEHP FREE | MR CONDITIONAL | CE MARK

The Infu-Surg® Clarity Pressure Infusion Bag is an innovative, user-friendly device for intravenous and intra-arterial infusion treatment. Designed with a clear panel, the patient label and barcode on the solution bag can be easily scanned. The dedicated fluid bag hook enables efficient loading and unloading of the system without having to remove the infusion bag from the IV pole.



Clear panel enables full visibility of solution bag to monitor fluid level and easily scan patient barcode



360° view of the color-coded gauge facilitates accurate pressure monitoring while the safety pressure relief valve prevents over-inflation



Three-way stopcock ensures precise control of pressure while the large oval bulb allows for quick bladder inflation, all with ONE hand



Solution bag hook enables easy bag changing without having to remove the Infu-Surg from the IV pole; keeps fluid bag in proper position while in use

## INFU-SURG® CLARITY PRESSURE INFUSION BAG

- Transparent panel for clear visibility of solution bag
- Dedicated hook enables quick change of solution bag while providing proper fluid bag positioning in pocket
- Pressure relief valve prevents over-inflation (375 mm Hg)
- Color-coded gauge provides accurate pressure monitoring (0-300 mm Hg)
- Single-handed inflation and deflation with oval-shaped bulb makes it easy to use and requires minimal training
- Welded panel and hook construction are more durable than sewn designs to withstand repeated use
- Suitable for use with external inflation sources
- Three-way stopcock ensures precise control of pressure
- Incredibly reliable – 100% tested

ITEM	DESCRIPTION	VOLUME	PK
4050	Infu-Surg Clarity	500 mL	25
4100	Infu-Surg Clarity	1000 mL	25

LIT-10016 Rev.2