



GREEN SYSTEM FIBER OPTIC LARYNGOSCOPE

LATEX FREE | CLINICALLY CLEAN | SINGLE PATIENT USE | DISPOSABLE

IntuBrite® disposable fiber optic blades and handles are built to ISO Green Standard System specification. The innovative greenline compatible dual LED lights combine white and ultraviolet wavelengths to illuminate airway structures. Its now possible to go disposable and improve the intubation experience.



Fiber Optic Light



Miller Blades



Handle Options



Macintosh Blades

INTUBRITE® GREEN SYSTEM FIBER OPTIC LARYNGOSCOPE

- Green System compatible
- Blended Ultraviolet + White LED lights
- UV light provides unique bright detail
- Improved clarity and definition without glare
- Lightweight, disposable aircraft aluminum handle
- Broad range of combinations

ITEM	DESCRIPTION	SIZE	PK
6020.C	Disposable Fiber Optic Blade	Mac 0	20/box
6021.C	Disposable Fiber Optic Blade	Mac 1	20/box
6022.C	Disposable Fiber Optic Blade	Mac 2	20/box
6023.C	Disposable Fiber Optic Blade	Mac 3	20/box
6029.C	Disposable Fiber Optic Blade	Mac 3 Plus	20/box
6024.C	Disposable Fiber Optic Blade	Mac 4	20/box
6040.C	Disposable Fiber Optic Blade	Miller 000	20/box
6039.C	Disposable Fiber Optic Blade	Miller 00	20/box
6030.C	Disposable Fiber Optic Blade	Miller 0	20/box
6031.C	Disposable Fiber Optic Blade	Miller 1	20/box
6035.C	Disposable Fiber Optic Blade	Miller 1.5	20/box
6032.C	Disposable Fiber Optic Blade	Miller 2	20/box
6033.C	Disposable Fiber Optic Blade	Miller 3	20/box
6034.C	Disposable Fiber Optic Blade	Miller 4	20/box
6015.C	Disposable Fiber Optic Handle	Medium	20/box
6016.C	Disposable Fiber Optic Handle	Pediatric	20/box
6017.C	Disposable Fiber Optic Handle	Stubby	20/box
1058	Storage Bag	-	1/box

LIT-10032 Rev 1