

# DUAL-LIMB ANESTHESIA CIRCUIT

LATEX FREE | SINGLE PATIENT USE | DISPOSABLE

Westmed™ brand Dual-Limb Circuits offer virtually unlimited configurations and customization.

9599 Dual-Limb Anesthesia Circuit - Adult with 96" expandable tubing, parallel wye, elbow with gas sampling port, Mini HMEF, 2 L breathing bag, gas sampling line with swivel elbow, and large adult mask

Elbow with and without gas sampling port

HMEs and HMEFs filters with luer gas sampling ports available

Unique ergonomic mask fits most facial configurations for improved sealing capabilities

Standard and swivel gas sampling lines

Space saving expandable tubing available from 72" to 144". Corrugated tubing also available.

Breathing bag options in 1, 2, and 3 liter

8267 Dual-Limb Anesthesia Circuit - Neonatal with 60" corrugated tubing, parallel wye, elbow with gas sampling port, two B/V filters, and gas sampling line with swivel elbow

Bacterial viral, and HEPA filters with luer gas sampling ports available

## DUAL-LIMB ANESTHESIA CIRCUITS

- Standard and swivel gas sampling ports, lines and filters
- High-quality tubing with exceptional flexibility and drapability
- ISO connections for secure, leak resistant connections
- Available in a wide variety of lengths and configurations to meet your custom needs

# DUAL-LIMB ANESTHESIA CIRCUIT

## ADULT STANDARD DUAL-LIMB ANESTHESIA CIRCUIT

ITEM	DESCRIPTION	WYE	ELBOW	FILTER	HME/HMEF	BAG	GAS SAMPLING LINE	MASK
9596	72" Expandable Circuit	Parallel	Gas Sampling	(2) Bacterial/Viral	–	2 L	10' x .050" ID w/Swivel Elbow	Large Adult 7700 Series
9598	72" Expandable Circuit	Parallel	Gas Sampling	–	Mini HMEF	2 L	10' x .050" ID w/Swivel Elbow	Large Adult 7700 Series
7-8113-10-62-68	96" Expandable Circuit	Parallel	Gas Sampling	Bacterial/Viral	–	3 L	–	–
8661	96" Expandable Circuit	Parallel	Gas Sampling	(2) Bacterial/Viral	–	3 L	10' x .050" ID M/M	–
9129	96" Expandable Circuit	Parallel	Gas Sampling	–	Mini HMEF	2 L EZ Gripper	10' x .050" ID M/M	–
9342	96" Expandable Circuit	Parallel w/Port	Gas Sampling	Bacterial/Viral	–	2 L	–	–
9590	96" Expandable Circuit	Parallel	Gas Sampling	Bacterial/Viral	–	3 L	–	–
9597	96" Expandable Circuit	Parallel	Gas Sampling	(2) Bacterial/Viral	–	2 L	10' x .050" ID w/Swivel Elbow	Large Adult 7700 Series
9599	96" Expandable Circuit	Parallel	Gas Sampling	–	Mini HMEF	2 L	10' x .050" ID w/Swivel Elbow	Large Adult 7700 Series
9063	120" Expandable Circuit	Parallel	Gas Sampling	(2) Bacterial/Viral	–	2 L	10' x .050" ID M/M	–
9203	144" Expandable Circuit	Parallel	Gas Sampling	–	–	–	–	–

## PEDIATRIC STANDARD DUAL-LIMB ANESTHESIA CIRCUIT

ITEM	DESCRIPTION	WYE	ELBOW	FILTER	HME/HMEF	BAG	GAS SAMPLING LINE	MASK
8569	96" Expandable Circuit	Parallel	Gas Sampling	Bacterial/Viral	–	1L	–	–
9196	108" Expandable Circuit	Parallel	Gas Sampling	(2) Bacterial/Viral	–	–	10' x .050" ID M/M	–
9361	120" Expandable Circuit	Parallel	Gas Sampling	Bacterial/Viral	–	1L	10' x .060" ID M/M	–
9204	144" Expandable Circuit	Parallel	Gas Sampling	–	–	–	–	–

## NEONATAL STANDARD DUAL-LIMB ANESTHESIA CIRCUIT

ITEM	DESCRIPTION	WYE	ELBOW	FILTER	HME/HMEF	BAG	GAS SAMPLING LINE	MASK
8267	60" Anesthesia Circuit	Parallel	Gas Sampling	(2) Bacterial/Viral	–	1/2 L	10' x .050" ID w/Swivel Elbow	–


CUSTOMIZE YOUR CIRCUIT

Choose from the following components to create a custom circuit:



BAGS

Three standard bag size options available (1, 2, and 3 L) and a 2 liter EZ gripper bag that is easy to grip with a gloved hand




TUBING

Flexible and drapable breathing tubes available in corrugated or space-saving expandable options




FILTERS

HMEs, HMEFs, bacterial viral, and HEPA filters with luer gas sampling ports available. Sampling ports located at the top of the housing, just as if a common gas sampling elbow were being used. ISO-standard connections for secure, leak-resistant attachment to breathing systems.



FLEX EXTENDERS

5.2" Flex Tubing with and without capped gas sampling port



FACE MASKS

Unique ergonomic mask shape that fits most facial configurations for improved sealing capabilities




GAS SAMPLING LINES

Standard and swivel lines available



WYES AND CONNECTORS

Elbows, parallel and swivel wyes available with and without port



YANKAUER

Sterile Yankauer with bulb tip and non-sterile Yankauer with regular tip available