

MEDISORB®

MEDICAL SODA LIME CO₂ ABSORBENT

LATEX FREE | DISPOSABLE | DEHP FREE

Medisorb® is a CO₂ absorbent designed for use with anesthesia systems in hospitals or surgery centers. It effectively removes CO₂ from breathing gases during anesthesia procedures.



M1173310
Compatible with
GE Aespire™, GE Avance™,
GE Aisys™ and anesthesia systems



Medisorb EF medical
soda lime CO₂ absorbent



427000100
Compatible with GE
Aestiva anesthesia systems



2079796-001
Compatible with GE Carestation™ 600 series

FEATURES

- Low-alkaline formulation available. It's specifically formulated to help reduce concern over anesthetic agent degradation.
- Color change indicates absorption progress.
- D-profile granule shape prevents clumping for better gas flow distribution.
- Compatible with GE Healthcare anesthesia machines.
- Safe to use with desflurane.
- Controlled moisture content prevents desiccation.
- Conforms with APSF recommendations for anesthetic agent degradation.

EF ABSORBENT

ITEM	DESCRIPTION	PK
M1173311	Multi-Absorber EF, Disposable, 0.8 kg	6
2079797-001	Medisorb EF EX, Disposable, 1 kg	6

ABSORBENT

ITEM	DESCRIPTION	PK
M1173310	Multi-Absorber, Disposable, 0.8 kg	6
8570043	Twin Pack, 4.5 kg Canister	2
2079796-001	Medisorb EX, Disposable, 1.2 kg	6
427000100	Medisorb, Disposable, 1.1 kg	12