


GREENLINE® FIBER OPTIC GERMAN PROFILE BLADE


LATEX FREE | REUSABLE

Fiber optic German profile laryngoscope blades are built with an integral, extra large, fiber optic bundle providing a more cool, white light. This integral design enables superior cleaning by eliminating cracks and crevices while preventing cross contamination. These blades are made of 303/304 surgical stainless steel, not chrome plated.



MACINTOSH GERMAN PROFILE

ITEM	DESCRIPTION	SIZE	A	B	PK
5-5272-01	Infant	1	93 mm	15 mm	1
5-5272-02	Child	2	100 mm	100 mm	1
5-5272-13	Strong Curve Medium Adult	3	135 mm	135 mm	1
5-5272-03	Medium Adult	4	135 mm	135 mm	1
5-5272-04	Large Adult	5	155 mm	135 mm	1



MILLER GERMAN PROFILE

ITEM	DESCRIPTION	SIZE	A	B	PK
5-5273-20	Premature	00	65 mm	11 mm	1
5-5273-00	Newborn	0	77 mm	12 mm	1
5-5273-01	Infant	1	110 mm	12 mm	1
5-5273-15	Small Child	1.5	130 mm	14 mm	1
5-5273-02	Child	2	152 mm	14 mm	1
5-5273-03	Medium Adult	3	192 mm	15 mm	1
5-5273-04	Large Adult	4	207 mm	16 mm	1
5-5243-02	Large Hand Adult	3	130 mm	22 mm	1

