

# SUNONE™ SYSTEM FIBER OPTIC GREENLINE®/D™ HANDLE

LATEX FREE | STERILE | SINGLE PATIENT USE | DISPOSABLE

GreenLine®/D™ disposable LED laryngoscope handles are ideal for use with GreenLine® Fiber Optic laryngoscope blades. This single patient use handle helps address secondary infection concerns and costs associated with reprocessing.



## GREENLINE®/D™ DISPOSABLE LED HANDLE

- One complete, fully-contained disposable unit; product recycling available
- Ideal for use with GreenLine®/D™ disposable blades
- Reduce the risk of cross-contamination with single patient use handles
- LED illumination provides exceptional, clear viewing
- Complies with ISO 7376 Standard
- Textured surface ensures secure grip
- Mercury-free, pre-installed batteries

ITEM	DESCRIPTION	SIZE	STERILE	PK
5-0236-69	Ribbed Grip	Large	Clean Pkg	20
5-0236-70	Smooth Grip	Penlite	Clean Pkg	20
5-5338-09	Textured Grip	Medium	Yes	20
5-5338-10	Textured Grip	Penlite	Yes	20
5-5338-11	Textured Grip	Stubby	Yes	20
5-5338-12	Textured Grip	Standard	Yes	20