

CONVENTIONAL LARYNGOSCOPE HANDLES

LATEX FREE | REUSABLE

Chrome plated and stainless steel waterproof handles are available from AirLife™ for use with conventional laryngoscope blades. Conventional Howland Lock is also available.



CHROME PLATED BRASS HANDLES

- Chrome plated brass
- Batteries are not included

STAINLESS STEEL WATERPROOF HANDLES

- Stainless steel waterproof handles may be immersed in disinfectant without removing batteries
- Batteries are not included

ITEM	DESCRIPTION	SIZE	BATTERIES	PK
5-0237-01	Chrome Plated Brass	Penlite	2 "AA"	1
5-0237-02	Chrome Plated Brass	Stubby	2 "AA"	1
5-0237-03	Chrome Plated Brass	Medium	2 "C"	1
5-0237-05	Chrome Plated Brass	Large	2 "D"	1
5-0237-08	Chrome Plated Brass	Mini	Lithium 123	1

ITEM	DESCRIPTION	SIZE	BATTERIES	PK
5-0237-91	Waterproof Handle	Penlite	2 "AA"	1
5-0237-92	Waterproof Handle	Stubby	2 "AA"	1
5-0237-93	Waterproof Handle	Medium	2" C"	1
5-0237-96	Waterproof Handle	Micro	Lithium 123	1