


DOUBLE LUMEN ENDOBRONCHIAL TUBE

STERILE | LATEX FREE | SINGLE PATIENT USE | DISPOSABLE

The Salter Labs® Double Lumen Endobronchial Tube helps with easier and faster endobronchial intubation. Conveniently packaged assembled, the clampless system uses valves to easily isolate one lung ventilation or one lung isolation during surgery. Available in left- and right-sided versions in multiple sizes.



Low-friction suction catheters offer the ability to adjust suction pressure



Clampless Y-connector offers scope and suction catheter access in either open or closed positions

DOUBLE LUMEN ENDOBRONCHIAL TUBE

- Isolates the left or right lung for surgery
- 90° valve eliminates the need for awkward clamps
- Y-connector is preassembled for ease and efficiency
- Scope and suction catheter access in either open or closed positions
- Low friction, adjustable pressure suction catheter has a textured outer wall for easier advancement
- Durable PVC tracheal cuff resists tearing during insertion
- Low profile, smooth tracheal opening can reduce trauma during insertion
- Distal tip angle facilitates insertion into bronchus

DOUBLE LUMEN (LEFT)

ITEM	SIZE	PK
H-PSTDL-28-1	28 Fr	1
H-PSTDL-32-1	32 Fr	1
H-PSTDL-35-1	35 Fr	1
H-PSTDL-37-1	37 Fr	1
H-PSTDL-39-1	39 Fr	1
H-PSTDL-41-1	41 Fr	1

DOUBLE LUMEN (RIGHT)

ITEM	SIZE	PK
H-PSTDR-28-1	28 Fr	1
H-PSTDR-32-1	32 Fr	1
H-PSTDR-35-1	35 Fr	1
H-PSTDR-37-1	37 Fr	1
H-PSTDR-39-1	39 Fr	1
H-PSTDR-41-1	41 Fr	1