


# DOUBLE LUMEN ENDOBRONCHIAL TUBE

LATEX FREE | STERILE | SINGLE USE | DISPOSABLE

The Robertshaw double lumen endobronchial tube is available in left and right designs. This tube helps provide an airtight seal of the trachea and one bronchus. Packaged sterile with dual connectors, suction ports, and suction catheters fully assembled.



## DOUBLE LUMEN ENDOBRONCHIAL TUBE

- Color coded clear (tracheal) and blue (bronchial) for easy differentiation of appropriate lumen, pilot balloon and cuff
- Soft PVC
- High volume, low pressure cuff
- Radiopaque strip
- Graduated depth marks
- Suction catheters included
- Bronchial stylette included
- 15 mm Y-connector included with each tube

## DOUBLE LUMEN (LEFT)

| ITEM      | DESCRIPTION                   | SIZE  | PK |
|-----------|-------------------------------|-------|----|
| 1-7370-35 | Endobronchial Robertshaw Left | 35 Fr | 1  |
| 1-7370-37 | Endobronchial Robertshaw Left | 37 Fr | 1  |
| 1-7370-39 | Endobronchial Robertshaw Left | 39 Fr | 1  |
| 1-7370-41 | Endobronchial Robertshaw Left | 41 Fr | 1  |

## DOUBLE LUMEN (RIGHT)

| ITEM      | DESCRIPTION                    | SIZE  | PK |
|-----------|--------------------------------|-------|----|
| 1-7371-35 | Endobronchial Robertshaw Right | 35 Fr | 1  |
| 1-7371-37 | Endobronchial Robertshaw Right | 37 Fr | 1  |
| 1-7371-39 | Endobronchial Robertshaw Right | 39 Fr | 1  |
| 1-7371-41 | Endobronchial Robertshaw Right | 41 Fr | 1  |

LIT-20040-C1 Rev. 1