

SUNONE™ SYSTEM: GREENLINE®/D™ FIBER OPTIC ALL-METAL BLADES

LATEX FREE | STERILE | SINGLE PATIENT USE | DISPOSABLE

GREENLINE®/D™ ALL-METAL BLADES

- Reduce the risk of cross-contamination with single patient use blades
- Fully disposable; product recycling available
- Five year sterility



5-5336-20 5-5336-00 5-5336-01 5-5336-15 5-5336-02 5-5336-03 5-5336-04



5-5335-04 5-5334-35 5-5334-03 5-5334-02 5-5334-01 5-5334-00

ALL-METAL MILLER

ITEM	DESCRIPTION	LENGTH	HEIGHT	SIZE
5-5336-20	Premature	65 mm	10 mm	00
5-5336-00	Neonate	79 mm	10 mm	0
5-5336-01	Infant	105 mm	11 mm	1
5-5336-15	Small Child	130 mm	13 mm	1.5
5-5336-02	Child	157 mm	13 mm	2
5-5336-03	Medium Adult	197 mm	13 mm	3
5-5336-04	Large Adult	208 mm	16 mm	4

ALL-METAL MACINTOSH


ITEM	DESCRIPTION	LENGTH	HEIGHT	SIZE
5-5335-00	Neonate	80 mm	15 mm	0
5-5335-01	Infant	89 mm	15 mm	1
5-5335-02	Child	110 mm	15 mm	2
5-5335-03	Medium Adult	131 mm	23 mm	3
5-5335-35	Ext. Med. Adult	145 mm	23 mm	3.5
5-5335-04	Large Adult	160 mm	23 mm	4

WIS-HIPPLE

ITEM	DESCRIPTION	LENGTH	HEIGHT	SIZE
5-5334-15	Small Child	115 mm	13 mm	1.5

PHILLIPS

ITEM	DESCRIPTION	LENGTH	HEIGHT	SIZE
5-5339-01	Infant	113 mm	12 mm	0
5-5339-02	Child	166 mm	18 mm	1



LIT-20043-C1 Rev.1