


LATEX FREE ROBERTAZZI AIRWAY

LATEX FREE | STERILE | SINGLE USE | DISPOSABLE

The latex free Robertazzi Nasopharyngeal Airway performs like latex. The soft silicone is flexible, yet maintains its shape. The unique packaging features extra large size notation for immediate visual reference and accurate airway selection. Available in nine sizes to accommodate a variety of airway types.



LATEX FREE ROBERTAZZI AIRWAY

- Flexible silicone performs like latex airways
- Soft material prevents damage to nasal passage
- Blunt tip
- Resists collapse or kinking
- Sets are ideal for crash carts

INDIVIDUAL

ITEM	SIZE	ID	LENGTH	PK
1-5075-20	20 Fr	5.0 mm	114 mm	10
1-5075-22	22 Fr	5.5 mm	114 mm	10
1-5075-24	24 Fr	6.0 mm	120 mm	10
1-5075-26	26 Fr	6.5 mm	121 mm	10
1-5075-28	28 Fr	7.0 mm	132 mm	10
1-5075-30	30 Fr	7.5 mm	139 mm	10
1-5075-32	32 Fr	8.0 mm	146 mm	10
1-5075-34	34 Fr	8.5 mm	157 mm	10
1-5075-36	36 Fr	9.0 mm	170 mm	10

SETS

ITEM	SIZES INCLUDED	PK
1-5075-04	20, 24, 28, 32 Fr and sterile lubricating jelly	1 Set
1-5075-06	22, 24, 26, 28, 30, 32 Fr and sterile lubricating jelly	1 Set
1-5075-09	20, 22, 24, 26, 28, 30, 32, 34, 36 Fr and sterile lubricating jelly	1 Set
1-5075-05	28, 30, 32, 34, 36 Fr and sterile lubricating jelly	1 Set

(1 of each size per set/pack)

ADJUSTABLE FLANGE AIRWAY

LATEX FREE | SINGLE USE | DISPOSABLE



LATEX FREE ADJUSTABLE ROBERTAZZI AIRWAY

- Flexible neoprene performs like latex airways
- Soft material prevents damage to nasal passage
- Soft, adjustable flange
- Blunt tip

ITEM	SIZE	ID	LENGTH	PK
1-5072-12	12 Fr	3.0 mm	95 mm	10
1-5072-14	14 Fr	3.5 mm	95 mm	10
1-5072-16	16 Fr	4.0 mm	95 mm	10
1-5072-18	18 Fr	4.5 mm	95 mm	10
1-5072-20	20 Fr	5.0 mm	115 mm	10
1-5072-22	22 Fr	5.5 mm	120 mm	10
1-5072-24	24 Fr	6.0 mm	120 mm	10
1-5072-26	26 Fr	6.5 mm	130 mm	10
1-5072-28	28 Fr	7.0 mm	135 mm	10
1-5072-30	30 Fr	7.5 mm	135 mm	10
1-5072-32	32 Fr	8.0 mm	145 mm	10
1-5072-34	34 Fr	8.5 mm	160 mm	10
1-5072-36	36 Fr	9.0 mm	175 mm	10